

# InteleOrchestrator

## 4.5 | User Guide



intelerad

[Intelerad.com](https://www.intelerad.com)

## COPYRIGHT

© 2024 Intelerad Medical Systems Incorporated. All Rights Reserved.

No portion of the contents of this publication may be reproduced or transmitted in any form or by any means without the express written permission of Intelerad Medical Systems Incorporated.

Except as expressly provided otherwise in writing, the information provided in this document is provided AS IS, without any condition or warranty, whether written, oral, express, implied, legal, or statutory. No condition or warranty of merchantability or fitness for a particular purpose applies to anything provided by Intelerad. Without limiting the foregoing, Intelerad nor any of its suppliers warrants this documentation is or will be accurate, error free or meets or will meet user requirements.

## TRADEMARKS

Intelerad®, Intelerad Medical Systems, and the Intelerad logo, IntelePACS®, InteleViewer™, Enterprise Viewer™, InteleOne®, InteleArchive™, InteleShare VNA™, InteleShare™, InteleShare Research™, InteleOrchestrator™, InteleHeart™, InteleScreen™, and InSight PACS™ are either registered trademarks or trademarks of Intelerad Medical Systems Incorporated or its affiliates.

## THIRD-PARTY TRADEMARKS

Adobe, Acrobat, and Reader are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries. Firefox is a registered trademark of Mozilla Foundation in the United States and other countries. Google Chrome browser is a registered trademark of Google Inc. in the United States and other countries. Internet Explorer is a registered trademark of Microsoft Corporation in the United States and other countries. Intel, Pentium, Pentium II Xeon, and Pentium III Xeon are trademarks or registered trademarks of Intel Corporation or its subsidiaries in the United States and other countries. Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. NVIDIA is a trademark or registered trademark of NVIDIA Corporation in the United States and/or other countries. Oracle and Java are registered trademarks of Oracle Corporation and or its affiliates. PowerScribe 360 is a trademark or registered trademark of Nuance Communications Inc. or its affiliates in the United States and/or other countries. Sun, Sun Microsystems, and Java are trademarks or registered trademarks of Sun Microsystems, Inc. in the United States and other countries. Safari, Mac, and OS X are trademarks of Apple Inc., registered in the United States and other countries.

All other brand names, product names, or trademarks belong to their respective holders.



Intelerad Medical Systems Incorporated  
800, boul. De Maisonneuve East, 14th floor  
Montreal (Quebec)  
H2L 4L8 Canada

Title: InteleOrchestrator User Guide  
Software version: 4.5 (Multiple Patient Multiple Use)  
Date: 2024-12-24  
Part number: CLSWEN4.5.0UG-O Issue 005

# TABLE OF CONTENTS

<b>Getting Started .....</b>	<b>9</b>
About This Help .....	10
Document Conventions .....	10
Document Revisions .....	11
Obtaining Printed Documentation .....	11
Comments and Questions .....	11
Contacting Intelerad Technical Support .....	12
About InteleOrchestrator .....	13
Logging into InteleOrchestrator .....	13
Understanding the InteleOrchestrator Interface .....	19
About the Worklist Column .....	21
Understanding the Patient View .....	25
Selecting a Shift or Work Hours .....	26
About Extender .....	31
Installing InteleOrchestrator's Extender Application .....	31
Additional Extender Actions .....	35

Setting Your Status to Away .....	36
Viewing Your Current Productivity .....	39
Entering RVU and WU Credits Manually .....	40
Managing Your User Names and Passwords .....	41
<b>Application Settings and User Preferences .....</b>	<b>44</b>
Viewing Your User Settings .....	45
Setting History Indicators .....	47
Restoring Your User Settings Layout .....	49
Adding Your Contact Information .....	52
Profile Settings .....	54
Creating a Profile .....	54
Setting Your Display Mode .....	59
Setting Your Default Page .....	60
Filtering Your Reading Queue .....	62
Including All Shifts On Your Worklist .....	64
Exam Double Click .....	66
Start Viewer .....	67
Start Dictation .....	69
Enabling Combinable or Consecutive Reports .....	71
Selecting Sounds .....	75
Setting Search Result Page Size .....	77
Worklist Max Show Count .....	79
Open Peer Review in Window .....	81
Changing the InteleOrchestrator Theme .....	83
Setting the Default Image Viewer and Dictation System .....	85
<b>Worklists and Exams .....</b>	<b>90</b>
Managing My Reading Queue .....	91

Customizing My Reading Queue .....	91
To reorder My Reading Queue worklists: .....	92
Manually Creating the My Reading Queue .....	93
To manually create your My Reading Queue: .....	93
To edit My Reading Queue after it is created manually: .....	94
Opening Exams .....	95
Managing Exam Assignments .....	96
Using Auto-Next .....	99
Using the Auto-Next Optimizer .....	106
Creating Your Own Worklist .....	107
Customizing the Work Area with Worklists .....	109
Searching in InteleOrchestrator .....	112
Quick Search .....	113
Searching in Exam History .....	113
Searching In A Worklist .....	115
Advanced Search .....	116
Using Folders .....	118
<b>Collaboration and Communication .....</b>	<b>126</b>
Collaborating With Other Radiologists .....	127
Sending Notifications and Alerts .....	127
Adding Notes to Exams .....	131
Assigning a Note to the Operations/Support Staff .....	140
<b>InteleOrchestrator Messaging .....</b>	<b>141</b>
About InteleOrchestrator Messaging .....	142
Opening Messaging .....	142
Understanding the Messaging Page .....	143
Opening the Messaging App Window .....	144

Changing your Messaging Status .....	145
To change your Messaging status to Busy: .....	146
Setting Location Information .....	147
Searching for Messaging Users .....	147
Sending Messages .....	148
Sending a Message .....	148
Sending a Message to 2+ Users: Ad Hoc Groups .....	149
Sending a Message to an Existing Group .....	151
Group Management in InteleOrchestrator Messaging .....	153
Private Groups .....	153
Creating a Public or Private Group .....	153
Deleting a Public or Private Group .....	158
Leaving a Public or Private Group .....	158
Sending an Exam Link within Messaging .....	158
All Users are Credentialed .....	159
No Users are Credentialed .....	160
Users with Mixed Credentials .....	161
<b>Dictation and Reports .....</b>	<b>164</b>
Dictating Reports .....	165
Combining Reports and Reading Reports Consecutively .....	169
Changing a Dictation File After Submission for Transcription .....	172
InteleOrchestrator Speech Commands .....	176
Configuring Dictation Microphone Buttons to Perform Commands .....	179
<b>Protocol Workflow .....</b>	<b>181</b>
Opening Exams .....	182
Performing an Advanced Search .....	182
Understanding the Protocol Workflow .....	183

<b>Document Viewer</b> .....	<b>187</b>
Verify your Document Viewer Configuration Settings .....	188
To ensure that the following document viewer settings are enabled: .....	188
Document Viewer Preferences .....	188
General Configuration .....	189
Documents Category Configuration .....	190
Viewing Scanned Documents .....	190
<b>Peer Review Workflow</b> .....	<b>194</b>
Peer Review Workflow .....	195
About the Peer Review Worklist .....	195
Opening Peer Reviews .....	195
Peer Reviews to Evaluate Worklist .....	202
Peer Reviews to Arbitrate Worklist .....	205
Peer Review Scorecard .....	212
Peer Review Analysis Report .....	213
Creating Custom Reports .....	214
Manage Report Data .....	218
<b>InteleOrchestrator Analytics</b> .....	<b>222</b>
Analytics Overview .....	223
Opening Analytics .....	223
Understanding the Analytics Page .....	225
Dashboard Widgets .....	226
Viewing Widgets .....	226
Applying Widgets .....	232
Reports .....	240
Generating Reports .....	245
Creating a Report .....	245

Save and Schedule a Report .....	248
MIPS Reporting .....	265
<b>My Support Queue .....</b>	<b>266</b>
About My Support Queue .....	267
Understanding the My Support Queue Interface .....	267
Creating and Editing My Support Queue .....	274
To edit or create a new My Support Queue: .....	274
Changing the Properties of My Support Queue .....	279

# 1

## GETTING STARTED

This section details how to log in to InteleOrchestrator, select a shift, install the correct version of Extender for your dictation microphone, and describes the layout of the InteleOrchestrator user interface.

### In this chapter:

About This Help .....	10
Contacting Intelerad Technical Support .....	12
About InteleOrchestrator .....	13
Logging into InteleOrchestrator .....	13
Understanding the InteleOrchestrator Interface .....	19
Selecting a Shift or Work Hours .....	26
About Extender .....	31
Setting Your Status to Away .....	36
Viewing Your Current Productivity .....	39
Entering RVU and WU Credits Manually .....	40
Managing Your User Names and Passwords .....	41

# About This Help

This help provides instructions for end users, such as radiologists and technologists, on how to use InteleOrchestrator. We recommend that you browse through the entire Help and follow the procedures at your own pace. If you are already familiar with some aspects of InteleOrchestrator, you can use this Help for reference.

- Depending on your user role, geographical location, and assigned user permissions, you may not have access to all features and application interface elements described in this document.

## Document Conventions

Several conventions are used throughout this document. A list of these and examples of their use are provided below.

Convention	Example
Text that you enter in a field, or on a command line are in <b>courier font</b> .	In the Date field, enter 2003/04/04.
Keyboard commands are in <b>SMALL CAPS AND BOLD</b> .	Press <b>CTRL+C</b> to copy text.
New terminology or concepts are <i>italicized</i> .	The process of automatically distributing the images is referred to as <i>autorouting</i> .
Interface elements, such as menus, buttons, options, and preferences are <b>bold</b> .	From the <b>Font</b> list, choose the desired font.
Menu selections are separated by vertical lines.	Choose <b>File   Print</b> to print this page.
Information that is important for a user to know when performing a task, such as prerequisite information or restrictions, is represented with a note icon  .	<ul style="list-style-type: none"><li>To view reports, you must have the Report privilege enabled in your user account.</li></ul>

Convention	Example
Information that is helpful to a user, such as when describing an alternate or simpler way to perform a task, is represented with a tip icon  .	 You can also use the <b>CTRL+T</b> keyboard shortcut to show or hide thumbnail images.
Information that warns users to potential problems in the outcome of what they are doing, such as data loss or data breach, is represented with a warning icon  .	 Image measurements are saved for the current application session only. If you exit the application, all measurements are lost.

## Document Revisions

Each issue contains the features up to the version in the Software Version column.

Issue	Document Release Date	Software Version
001	Sep 19, 2023	4.4
002	May 21, 2024	4.5
003	Aug 01, 2024	4.5
004	Aug 30, 2024	4.5

## Obtaining Printed Documentation

Intelerad offers printed and bound versions of product documentation free of charge. To request printed copies of Intelerad documentation, contact your Client Success manager. The printed documents will be provided within 7 days or less.

## Comments and Questions

At Intelerad, we strive to create accurate and intuitive documentation that provides you with effective product training and troubleshooting support. To better help us develop documentation products that meet your needs, we encourage you to send your comments and questions to [documentation@intelerad.com](mailto:documentation@intelerad.com).

# Contacting Intelerad Technical Support

Your PACS administrator can assist you with any issues you may encounter. If you require additional assistance, you can contact Intelerad Technical Support, 24 hours a day, seven days a week.

To contact us:	Use:
On the Internet	<a href="https://serviceportal.intelerad.com/csm">https://serviceportal.intelerad.com/csm</a>
By telephone	Toll-free North America: 1-866-951-6222 Sans frais Amérique du Nord (français): 1 844-467-7227 Toll-free Australia: 1-800-286-418 Toll-free New Zealand: 0800-467-723 United Kingdom: 0113-360-2615 Other: +1-514-931-7127

These coordinates and a wealth of other information are also available on the Intelerad Service Portal.

<https://serviceportal.intelerad.com/csm>

You should regularly check the Intelerad knowledge base for the latest version of the documentation, as well as other product-specific resources such as TechNotes, downloads, and videos.

When you contact Intelerad Technical Support to report a problem, please have at hand the following information, as applicable:

- client code and location of your IntelePACS installation
- full error message and the steps required to reproduce the problem
- AE Titles of the affected devices
- operating systems of any affected machines
- description of the problem and when it first occurred

If the problem affects a particular study, please also provide the following:

- patient ID or patient number (M.R.N.)
- accession number/requisition number
- modality type and name

# About InteleOrchestrator


InteleOrchestrator offers a functionally rich, vendor-neutral, zero-footprint worklist that aims to achieve worklist perfection for each individual radiologist in your imaging group. Designed with flexibility in mind, InteleOrchestrator ensures that the appropriate radiologist is presented with the appropriate case at the right time, regardless of the complexity or disparate nature of your environment. With the ability to launch over 50 other applications through highly-tailored and dynamic worklists, InteleOrchestrator can adapt to the unique needs of any enterprise.

InteleOrchestrator offers:


- Support for multiple image viewers to integrate disparate medical image repositories.
- "Auto-Next" reading mode that radiologists can trust to present them with the most appropriate exam, without having to return to the worklist.
- Exams distributed across the practice automatically based on real-time variables such as availability, subspecialty, workload, location, and more.
- Dynamic worklist activations to automatically show and hide overflow worklists based on criteria-driven thresholds.
- AI-enabled worklist that can be augmented with actionable data through a vendor-neutral interface.

## Logging into InteleOrchestrator


1. Open a web browser session and navigate to the InteleOrchestrator URL that you were provided. If you don't know the URL contact your system administrator.  
 You must include the HTTPS:// portion of the URL for proper operation of the browser.
2. In the **Login** field, enter your username.
3. In the **Password** field, enter the password.



4. Select the **Domain** from the drop-down list.



5. Click **Login**.




6. InteleOrchestrator software allows for several different configurations, which may affect the screen that is displayed upon your initial login. If this is the case, from the Home Page, click **Worklist**.




7. If the application prompts you to install Extender Desktop Application, Click **OK** to download. If you do not have administrator privileges on your workstation, contact

your system administrator. For more information, see "[About Extender](#)"



8. The application may prompt you to select your current **Shift**, **Location** and **Room** if it is not already set.



- ⓘ Your practice may have integrated with a scheduling software, such as QGenda or Lightning Bolt or your InteleOrchestrator system administrator has already configured shifts for you. In this case, the InteleOrchestrator will already know the shift you are working when you log in and you will not be prompted to select your shift. For more information, see "[Selecting a Shift or Work Hours](#)."

 You may be prompted to select a shift even if you had previously selected and saved a shift. This can happen if you are assigned to more than one shift. For example, you are a Radiologist who is scheduled for a call shift and a day shift. You will need to select your primary shift when you log in to start your day shift. By selecting your primary shift, you will be assigned that shift's My Reading Queue and receive work unit credits for that shift. This is important if your organization is using the Workload Gauge.

- If your group is using a shift-based workflow, select the **Type** and the **Shift**.

Select Current Shift and Location

Select Shift       Set Working Hours

Type: Radiologist      Shift: Intel shift 1


Current Shift Details

Name:	Intel shift 1
Description:	
Location:	
Time/Date:	8:00 AM on 26-Nov-2024 to 11:40 PM on 26-Nov-2024
RVU target:	100
Work Unit target:	50

Reading Location

Location: Room:

**Save**



- If your group is not using a shift-based workflow, you can also set your working hours from this form. This may be useful to track how many hours you are working and your productivity during those hours.


Select Current Shift and Location

Select Shift       Set Working Hours

Time Range  
From: 07/02/2023  11:00 AM  To: 07/02/2023  4:00 PM

Reading Location  
Location:  Room:


**Save**



9. Once you have selected your shift or set your working hours, click **Save**.  
 Your practice may require you to provide certain information on this page, such as Shift, Location, and Room. If any required information is missing, an error message will prompt you to return and provide the necessary details.
10. The InteleOrchestrator interface is displayed.

## Understanding the InteleOrchestrator Interface


The InteleOrchestrator interface displays all of your worklists, the exams in your worklists, links to other InteleOrchestrator modules, and the means to communicate with other radiologists in your practice.



Section	Contents
Worklist Column	Displays your worklists, your exam history, and a quick search option. See <a href="#">"About the Worklist Column" (page 21)</a>
Patient View Tab	Displays exams that you opened for dictation but have not completed. This allows you to return to your work quickly if you are interrupted or need to change tasks quickly. See <a href="#">"Understanding the Patient View" (page 25)</a> <ul style="list-style-type: none"> <li>A Patient View tab will not open if you open an exam in View mode. It is only created if you open an exam in dictation or view and dictation mode. See <a href="#">"Opening Exams" (page 95)</a>.</li> </ul>
InteleOrchestrator Menu	Displays links to other InteleOrchestrator modules you may have access to.
Work Area	Displays the contents of your current worklist. If you click on a worklist in the Worklist Column, the exams in that worklist display in this area. See <a href="#">"Worklists and Exams" (page 90)</a> for more information. If you open an exam, the Patient View for that exam displays in this area. See <a href="#">"Understanding the Patient View" (page 25)</a> .
Communication Taskbar	Displays links to communicate with your colleagues. Includes messaging, alerts and notifications. See <a href="#">"Collaboration and Communication" (page 126)</a>

# About the Worklist Column

The Worklist column appears on the left side of the InteleOrchestrator interface, and contains panels from which you can access your worklists, view exam history, and search for exams.



## The Worklist Panel

The worklist panel contains all worklists that you have access to view your **My Reading Queue**.

**My Reading Queue** is considered your “master worklist”. If you are using a shift-based workflow, you will see all worklists associated with a shift under your **My Reading Queue**. These worklists will automatically populate as soon as your shift is selected.

Expand **My Worklist** to locate **My Reading Queue** in the worklist panel.

The screenshot shows the 'My Worklist' interface with the 'My Reading Queue' section expanded. The queue contains the following items:

- ED Overflow
- My Procedures 0
- Unread General 66
- My Subspecialty with exclusions 0

Below this, the 'My Support Queue' section is expanded, showing:

- Critical Results Worklists 14
- Test 12
- Today's Communications 0

The 'Unread Exams' section is also expanded, showing 306 unread exams. Other collapsed sections include 'EPIC', 'Time test 1', and 'Imaging Center 1'.

### To view the contents of a worklist:

1. Click the icon to expand the queue to display the available worklists.
2. Click a worklist in the worklist section to view the exams it contains in the work area.

The screenshot shows the 'Worklist' section expanded, displaying the 'Unread General Worklist'. The table lists the following exams:


Site	Accession	Actions	Status	Patient Name	Site Procedure	Priority	All
Imaging Center 1	131956		Unread	131956, 106664	NECK	Routine	
Imaging Center 1	CH27243		Unread	SAQUING, NICHOLAS	SKELETAL DIAGNOSIS	Routine	
Imaging Center 1	EGU81213		Unread	MASON, CHRISTIAN	ARM	Routine	
Imaging Center 1	EIT53133		Unread	BARBER, CHRISTIAN	NECK SCREENING	Routine	
Imaging Center 1	GKH63034		Unread	BROWN, SAMIR	LEG	Routine	
Imaging Center 1	HXJ69516		Unread	MARTINEZ, SARAH	SHOULDER	Routine	
Imaging Center 1	IHT80847		Unread	MCDONALD, NICOLE	ABDOMEN DIAGNOSIS	Routine	
Imaging Center 1	IMSQA29444		Unread	IMSAQA29444, IMSQA	PETCTC ONLY SKULL BASE TO MIDTHIGH	Routine	
Imaging Center 1	IMSQA29618		Unread	IMSAQA29618, IMSQA	KNEE ONE OR TWO VIEWS	Routine	
Imaging Center 1	IMSQA80406		Unread	NO, PATNAME INTIMAGE	ABDOMEN TEST	Routine	
Imaging Center 1	INV26920		Unread	RICHARDS, PASCAL	CAROTID	Routine	
Imaging Center 1	JDP73571		Unread	MCDONALD, CHRISTIAN	ABDOMEN DIAGNOSIS	Routine	Sub
Imaging Center 1	JTG45239		Unread	ST-JEAN, ISABELLE	BROKEN CAROTID	Routine	

### Auto-Next Functionality:

Auto-Next is a useful tool when considering Radiologist Workflow. The auto-next feature ensures that every exam in your worklist is read in the predetermined order. This order is

established according to how the worklists are set up, i.e., SLA rules, Shift management, subspecialty groups, etc.

Auto-Next controls are located above the worklists. Click the Auto-Next start icon,  to start the Auto-Next functionality. See “[Using Auto-Next](#)” for more information.



## History Panel

The **History** panel contains all the exams that you have dictated or taken action on the specified date range. This includes adding notes or changing exam information. You can search your exam history by selecting ‘Today’, ‘1wk’ or enter a custom date range by clicking the clock icon. See [Searching in Exam History](#) for additional details. To set which actions will appear in this list, see “[Setting History Indicators](#)” (page 47). The most common actions include:

- Exams you dictated
- Images you viewed
- Notes you created
- Information about exams you edited


History (12)											
	Today	1wk	Filter	Info							
<Accession, Procedure, Patient Name>											
Accession	Procedure	Patient	Draft	ED	Follow-up	New	New	Out	Draft	Waitlist	
7...	SHO...	W...	█					█			
N...	HAN...	B...	█					█			
T...	ABD...	C...	█					█			
J...	LUN...	S...	█					█			

💡 You can hover your mouse over the Report icon █ to view that patient's report, or you can click an exam to go to that exam's patient view.

History (12)											
	Today	1wk	Filter	Info							
<Accession, Procedure, Patient Name>											
Accession	Procedure	Patient	Draft	ED	Follow-up	New	New	Out	Draft	Waitlist	
7...	SHO...	W...	█					█			
N...	HAN...	B...	█					█			
T...	ABD...	C...	█					█			
J...	LUN...	S...	█					█			


## The Search Panel

You can use the Search panel to perform a quick search. The results of your search appear in the work area. See “[Quick Search](#)” (page 113) for instructions on how to use this search function.



## Understanding the Patient View

The Patient View displays all the patient and exam information in a single view for a specific patient. The Patient View opens when you open an exam to view the exam images, click on an exam in the worklist, or when you double-click an exam for dictation. There are five main sections in the Patient View dialog.



Callout	Section	Description
1	<b>Patient Information</b>	Contains patient information, demographics, and contact information. Your InteleOrchestrator system administrator configures the information displayed in this section.
2	<b>Report Area</b>	<ul style="list-style-type: none"> <li>Displays the order information for exams in the Order tab</li> <li>Displays the report information in the Report tab if a report has been created for this exam</li> <li>Displays the dictation area for you to dictate your</li> </ul>

Callout	Section	Description
		<p>reports</p> <ul style="list-style-type: none"> <li>• Displays any uploaded documents</li> </ul>
3	<b>Notes</b>	Displays all the notes associated with the exam. The notes are sorted with the most recent at the top. You can also create new notes or edit or complete existing notes.
4	<b>CRM</b>	The Customer Relationship Management section displays information pertaining to your customer sites and ordering physicians that is captured in the CRM module.
5	<b>Exam Information</b>	Lists the patient exams. The exam you opened from the worklist is selected in the exams list. You can view any other exam by clicking it or by hovering your mouse over the Report icon. If you are using InteleOrchestrator Master Patient Index (MPI), you will see all the reports for this patient that reside across multiple facilities. You can also view exam history and edit exam information.

## Selecting a Shift or Work Hours

If your practice does not use a scheduling application, such as QGenda or Lightning Bolt, or if your InteleOrchestrator system administrator did not configure shifts for you, you may need to manually select a shift when you log in to InteleOrchestrator. If this is the case, if you do not select a shift, you will have to manually create your own My Reading Queue or work from an individual worklist. See ["Managing My Reading Queue" \(page 91\)](#).

-  If your practice uses a scheduling application or your InteleOrchestrator system administrator created shifts for your practice, your shift is automatically set when you log in to InteleOrchestrator.

### To change your shift:

1. Log in to InteleOrchestrator. See ["Logging into InteleOrchestrator" \(page 13\)](#).
2. In the Select Current Shift dialog, click **Select Shift** or **Set Working Hours**.

Select Current Shift and Location

Select Shift       Set Working Hours

Type: Radiologist, Radiology Assistant,  Shift:

Shift Detail

Reading Location

Location:  Room:

3. Do one of the following:

- To select a shift, select one or more options in the **Type** list, and then select an option from the **Shift** list to set the shift values.

Select Current Shift and Location

Select Shift       Set Working Hours

Type: Radiologist      Shift: Intel shift 1

Current Shift Details

Name:	Intel shift 1
Description:	
Location:	
Time/Date:	8:00 AM on 26-Nov-2024 to 11:40 PM on 26-Nov-2024
RVU target:	100
Work Unit target:	50

Reading Location

Location: Room:

-  In some cases, there may be multiple instances of the shift that are available for selection. Multiple shift instances will be displayed when this happens. Select one of the instances displayed.

Select Current Shift and Location

Select Shift       Set Working Hours

Type: Radiologist, Administration, Radi Shift: Intel Shift

Instance:  1:30 PM on 18-Nov-2024 to 9:30 AM on 19-Nov-2024  
 1:30 PM on 19-Nov-2024 to 9:30 AM on 20-Nov-2024


Next Shift Details

Name:	Intel Shift
Description:	
Location:	
Time/Date:	1:30 PM on 19-Nov-2024 to 9:30 AM on 20-Nov-2024
RVU target:	-
Work Unit target:	-

Reading Location

Location:	New regression2
Room:	New_Test

Save



- To set your working hours, in the **From** and **To** fields, enter the date and time range

Select Current Shift and Location

Select Shift       Set Working Hours


Time Range

From: 07/02/2023 11:00 AM To: 07/02/2023 4:00 PM

Reading Location

Location:	
Room:	

Save



4. If prompted, select the **Location** and **Room**.

Select Current Shift and Location

Select Shift       Set Working Hours

Type: Radiologist      Shift: Intel shift 1

Current Shift Details

Name:	Intel shift 1
Description:	
Location:	
Time/Date:	8:00 AM on 26-Nov-2024 to 11:40 PM on 26-Nov-2024
RVU target:	100
Work Unit target:	50

Reading Location

Location:	Room:
-----------	-------


**Save**

5. Click **Save**.

- Your practice may require you to provide certain information on this page, such as Shift, Location, and Room. If any required information is missing, an error message will prompt you to return and provide the necessary details.

## Changing Shift or Working Hours

If you need to change your shift after the initial selection at login, you can do this by clicking the “Shift icon” in the top of the worklist panel.



## About Extender

The InteleOrchestrator uses the Extender to launch your Viewer and Dictation system. The Extender must be downloaded before these systems can be launched from InteleOrchestrator. You need to do this only once per workstation. After Extender is installed and configured, it will communicate with multiple viewer, dictation, and EHR/EMR client-side or web-based applications to update exams as you read them.

- If your practice is using InteleOrchestrator Desktop Monitoring Service, you will not need to install Extender. Ask your InteleOrchestrator system administrator if this step is required.


## Installing InteleOrchestrator's Extender Application

The InteleOrchestrator's Extender application should always be installed.

### To install InteleOrchestrator's Extender application:


- Log in to InteleOrchestrator.
- If this is the first time you are logging in to InteleOrchestrator, InteleOrchestrator may prompt you to download the Desktop Installer.

2. Click your user name in the upper-right corner of the InteleOrchestrator interface.



The screenshot shows the InteleOrchestrator software interface. At the top, there is a navigation bar with links for CRM, Protocoling, and Dictation. To the right of these links is a dropdown menu containing the text "Administrator, Initial". This dropdown menu is highlighted with a red rectangular box. Further to the right of the dropdown are two buttons: a blue "Logout" button and a blue circular button with a white icon. Below the navigation bar, there is a search bar with the placeholder text "<Accession, Procedure, Patient Name>". To the right of the search bar are several small icons: a magnifying glass, a printer, and a gear. The main workspace below the toolbar contains a table with columns labeled "Procedure", "Priority", "AI Status", "Procedure Code", and "Exam". There are two rows of data in the table. The first row shows "PROCEDURE 1 VIEW" in the Procedure column, "Stat 1" in the Priority column, "AI Status" in the AI Status column, and "DX.ABD1V" in the Procedure Code column. The second row shows "PROCEDURE 1 VIEW" in the Procedure column, "Stat" in the Priority column, "AI Status" in the AI Status column, and "DX.ABD1V" in the Procedure Code column. The word "Stat" in the second row's Priority column is displayed in red.

3. Click **About**.



4. The **About** box displays the InteleOrchestrator's Extender application installation choices. Click **Extender Desktop Installer**.



5. Run the InteleOrchestrator's Extender application file.

- Once the installation completes and InteleOrchestrator's Extender application is running, the InteleOrchestrator's Extender application  icon appears on your workstation taskbar.

## Additional Extender Actions

Occasionally, there may be unusual behavior when you are using InteleOrchestrator. To troubleshoot these issues, your InteleOrchestrator system administrator may need access to logs generated by your Extender application. Introduced in version 4.5 R3, log uploading occurs automatically in the background without any action by you. In versions 4.5 R2 and earlier, logs must be manually uploaded.

Here are some Extender actions to take that can make logs available to help your InteleOrchestrator system administrator troubleshoot issues.

Action	Reason	Instructions
Reinitialize	If your InteleOrchestrator system administrator asks you to reinitialize or restart Extender.	<ol style="list-style-type: none"><li>1. Right click the Extender  icon in your taskbar at the bottom of your screen.</li><li>2. Click <b>Reinitialize</b>.</li></ol>
View Initialization Logs	File Explorer opens to the location of the initialization log, allowing easy access to logs recorded before the completion of the Extender initialization. This log file contains any relevant errors encountered during this phase, including connection errors and log upload failures.	<ol style="list-style-type: none"><li>1. Right click the Extender  icon in your taskbar at the bottom of your screen.</li><li>2. Click <b>View Initialization Logs</b>.</li></ol>
Restart Log	Creates a new log file to help your InteleOrchestrator system administrator and Inteleraid Technical Support to troubleshoot issues.	<ol style="list-style-type: none"><li>1. Right click the Extender  icon in your taskbar at the bottom of your screen.</li><li>2. Click <b>Restart Log</b>.</li></ol>

Action	Reason	Instructions
Clear the Extender log  (Only version 4.5 R2 and earlier)	This will delete the existing Extender logs and restart logging. This action allows your InteleOrchestrator system administrator and Intelerad Technical Support to access the new log files.	<ol style="list-style-type: none"> <li>1. Right click the Extender  icon in your taskbar at the bottom of your screen.</li> <li>2. Click <b>Clear Log</b>.</li> </ol>
Upload Log  (Only version 4.5 R2 and earlier)	Uploads the logs to InteleOrchestrator to make the logs available to your InteleOrchestrator system administrator.	<ol style="list-style-type: none"> <li>1. Right click the Extender  icon in your taskbar at the bottom of your screen.</li> <li>2. Click <b>Upload log</b>.</li> </ol>

## Setting Your Status to Away

InteleOrchestrator provides you with the ability to set your status as “Away”. Setting a status as away locks your current InteleOrchestrator session.


When you are ready to resume working, InteleOrchestrator prompts you to select the task that you were performing while away, and credits you with the appropriate RVUs/WUs based on the amount of time you were away and the selected task.

When setting your status to away, you can also choose to:

- Suspend exam assignments
- Suspend the Auto-Next functionality
- Set your Messaging status to Away


## To activate an away status:

1. Click the **Lock System (Away)**  icon (next to **Logout**).
2. In the lock screen, you can choose to apply one or more of the following additional away functions:
  - **Pause Exam Assignment:** Exams are not automatically assigned to you by the InteleOrchestrator Assignment Engine while you are away. However, Operations users can still manually assign exams to you.
  - **Pause Auto-Next:** The Auto-Next functionality is paused while you are away. See "[Using Auto-Next](#)" (page 99) for more details.
  - **Set chat status to Away:** You will appear as "Away" to all other Messenger users. See ["About InteleOrchestrator Messaging"](#) for more information.




**To resume working after an away status:**

1. Enter your password to unlock InteleOrchestrator.




2. From the **Task** list, select the task that you were completing while you were away.

This list of possible tasks is managed by your practice administrators.



3. If required, enter the task details.

4. Click **Save**.



## Viewing Your Current Productivity

InteleOrchestrator offers two productivity gauges that you can use to view your current productivity: Productivity Gauge and Workload Gauge.

You can also view the exams for which you received RVU and work units for today.


If implemented by your practice, the gauges appear above the worklist menu area in InteleOrchestrator.

### To see the exams that you have received RVU credit for today:

- In the Worklist menu, click RVU/WU Credit Detail .


### About the Productivity Gauge

The Productivity gauge displays your current productivity for the shift you are working based on an productivity goal that is predetermined by your team. The Productivity gauge consists of two bars. The top bar shows your current productivity followed by the number of RVUs that you have dictated and the percentage in relation to your goal. The bottom bar shows where you should be in relation to your goal for the time of day. This gauge can be configured to show your productivity in RVUs, Work Units, or both. The below image is showing both units.



## About the Workload Gauge

The Workload gauge is configured by your practice, and it displays your current productivity in relation to your practice's current workload. The Workload gauge displays in a red/yellow/green format, and your practice may also choose a fourth color, magenta. The color indicates your productivity in relation to your practice's workload.



The practice factor (PF) indicates your practice's current workload.

PF Value	PF Value Meaning
1	Your practice has the exact amount of work it can handle
>1	Your practice is has more work than it can handle
<1	Your practice has less work than it can handle

Your gauge will also give you insight into your and your practice's workload.


Gauge Color	Meaning
Red	You are behind or your practice is very busy
Yellow	You are slightly behind or your practice is slightly busy
Green	You and your practice are right on target
Magenta	You are reading too fast or your practice is very slow

## Entering RVU and WU Credits Manually

Many of the radiologist's Relative Value Unit (RVU) or Work Unit (WU) credits come from the dictation of exams. However, you may complete non-reading tasks at times, such as a tumor board meeting. If your InteleOrchestrator system administrator has given you permission, you must manually enter RVU or WU credits for these tasks.

## To manually enter RVU or WU credits:

1. At the top of the Worklist column, click the **Enter RVU/WU Credits** icon.
2. In the RVU/WU Credits dialog, select the type of task you performed from the **Task** drop-down list.
3. In the **Duration** field, enter the amount of time you spent on that task.
4. If required, change the **Date** field to reflect when you performed the task.
5. If required, enter a comment in the **Message** field.
6. Click **Save**.




## Managing Your User Names and Passwords


Your group may manage all member passwords via LDAP. If not, you can use these instructions to change your password in InteleOrchestrator. If you launch an image viewer or a dictation engine that is not synchronized through LDAP, you also enter your user name and password for those systems here to log into them.

### To specify your user name and new password:


1. Log in to InteleOrchestrator. See "[Logging into InteleOrchestrator](#)" (page 13).
2. Click your user name in the upper-right corner of the interface.




3. Click **Password Management**.




4. In the Password Management dialog, enter your current password in the **Current** field.



5. Enter your new password in the **New** and **Verify** fields.



6. Click **Save**.



# 2

# APPLICATION SETTINGS AND USER PREFERENCES

While most application preferences and settings are set by your InteleOrchestrator system administrator, there are settings that you can change to customize the appearance of the application, ensure that other users can communicate with you, and control how the application notifies you of certain events.

## In this chapter:

Viewing Your User Settings .....	45
Setting History Indicators .....	47
Restoring Your User Settings Layout .....	49
Adding Your Contact Information .....	52
Profile Settings .....	54
Creating a Profile .....	54
Setting Your Display Mode .....	59
Setting Your Default Page .....	60
Filtering Your Reading Queue .....	62
Including All Shifts On Your Worklist .....	64
Exam Double Click .....	66
Start Viewer .....	67


Start Dictation .....	69
Enabling Combinable or Consecutive Reports .....	71
Selecting Sounds .....	75
Setting Search Result Page Size .....	77
Worklist Max Show Count .....	79
Open Peer Review in Window .....	81
Changing the InteleOrchestrator Theme .....	83
Setting the Default Image Viewer and Dictation System .....	85

## **Viewing Your User Settings**


You can view all user settings currently assigned to you, such as what groups you belong to, your subspecialties, your practices, and your roles. View your settings to confirm they are correct. If not, contact your InteleOrchestrator system administrator to modify these settings.

**To view your user settings:**


1. Click your username in the upper-right corner of the InteleOrchestrator interface.



2. Click **User Settings**.



3. User settings are displayed.




## Setting History Indicators

You can view, add, or remove the history indicator icons that appear in the History section of the Worklist tab. See [About the Worklist Column](#) for more details. The icons you make active will appear next to relevant exams and alert you to actions needed for an exam or an exam's status.


### To select history indicators:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **User Settings**.

3. Under the History Indicators, add or remove history indicators as required. To add a History Indicators icon, click a blank space in the desired History Indicators icon row under the Active column. To remove a History Indicators icon, click the checkmark in the History Indicators icon row under the Active column so that the checkmark no longer appears.



4. Click **Save**.




## Restoring Your User Settings Layout

If you make changes to your User Settings, you can reset those values to their original value. For example, if you changed the worklist layout, such as column width and order, and panel sizing, you could restore them to the default state. Additionally, if you clicked a


"Never show me this again" response to a notification but need to see them after all, you can reinstate any pop-up settings that you selected to never show.

## To restore your worklist layout:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **User Settings**.
3. At the bottom of the User Settings dialog, click **Restore Layout**.



4. A new pop-up window appears.




5. Click **Yes** to clear out all your application data to restore layout or Click **No** to cancel.

**To receive notifications again:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **User Settings**.

- At the bottom of the User Settings dialog, click **Clear Confirmation Settings**.



- The user settings are cleared. Click Close to close the User Settings dialog box.




If you click the Close button the User settings dialog box will be closed without saving the changes.

## Adding Your Contact Information

For Messaging and Email notification functions to work, you must add your contact information in InteleOrchestrator.

**To add your contact information:**

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **Contact Information**.



3. Enter your contact information in the corresponding fields. Any information you add will appear in the Messaging application.
- Tip:** To add multiple entries for a category, click **+**. To delete contact information, select the row you want to delete and click **trash** icon.

Contact Information

Contact Information																												
<b>Phone</b> <table border="1"> <tr> <td>Phone Number</td> <td>Description</td> <td>Type</td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td>970-000-1152</td> <td></td> <td>Home</td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td colspan="5"><input type="button" value="Add"/></td> </tr> </table>						Phone Number	Description	Type	<input type="button" value="Delete"/>	<input type="button" value="New"/>	970-000-1152		Home	<input type="button" value="Delete"/>	<input type="button" value="New"/>	<input type="button" value="Add"/>												
Phone Number	Description	Type	<input type="button" value="Delete"/>	<input type="button" value="New"/>																								
970-000-1152		Home	<input type="button" value="Delete"/>	<input type="button" value="New"/>																								
<input type="button" value="Add"/>																												
<b>Fax</b> <table border="1"> <tr> <td>Fax Number</td> <td>Description</td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td colspan="5"><input type="button" value="Add"/></td> </tr> </table>						Fax Number	Description	<input type="button" value="Delete"/>	<input type="button" value="New"/>	<input type="button" value="Add"/>																		
Fax Number	Description	<input type="button" value="Delete"/>	<input type="button" value="New"/>																									
<input type="button" value="Add"/>																												
<b>Address</b> <table border="1"> <tr> <td>Address Line 1</td> <td>Address Line 2</td> <td>City</td> <td>State</td> <td>Zip</td> <td>Type</td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td>House No 104</td> <td>CSK Prestige Villas</td> <td>Shamshabad</td> <td>CA</td> <td>501218</td> <td>Home</td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td colspan="7"><input type="button" value="Add"/></td> </tr> </table>						Address Line 1	Address Line 2	City	State	Zip	Type	<input type="button" value="Delete"/>	<input type="button" value="New"/>	House No 104	CSK Prestige Villas	Shamshabad	CA	501218	Home	<input type="button" value="Delete"/>	<input type="button" value="New"/>	<input type="button" value="Add"/>						
Address Line 1	Address Line 2	City	State	Zip	Type	<input type="button" value="Delete"/>	<input type="button" value="New"/>																					
House No 104	CSK Prestige Villas	Shamshabad	CA	501218	Home	<input type="button" value="Delete"/>	<input type="button" value="New"/>																					
<input type="button" value="Add"/>																												
<b>Email</b> <table border="1"> <tr> <td>Email Address</td> <td>Description</td> <td><input type="radio"/></td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td>bill.moore@intelerad.com</td> <td>test</td> <td><input type="radio"/></td> <td><input type="button" value="Delete"/></td> <td><input type="button" value="New"/></td> </tr> <tr> <td colspan="5"><input type="button" value="Add"/></td> </tr> </table>						Email Address	Description	<input type="radio"/>	<input type="button" value="Delete"/>	<input type="button" value="New"/>	bill.moore@intelerad.com	test	<input type="radio"/>	<input type="button" value="Delete"/>	<input type="button" value="New"/>	<input type="button" value="Add"/>												
Email Address	Description	<input type="radio"/>	<input type="button" value="Delete"/>	<input type="button" value="New"/>																								
bill.moore@intelerad.com	test	<input type="radio"/>	<input type="button" value="Delete"/>	<input type="button" value="New"/>																								
<input type="button" value="Add"/>																												

- Click **Save**.


## Profile Settings

You may have multiple profiles associated with your InteleOrchestrator user account. Different profiles allow you to apply different settings based on the location you are working at or the role you are working in that day. This can be useful for situations such as working in different locations that use a different computer or workstation setups. If you have multiple profiles, you can then select the profile that best fits your work needs for that location or workflow for the exams you are reading.


## Creating a Profile

To create a profile:


- Click your user name in the upper-right corner of the InteleOrchestrator interface.
- Click **Profile Management**.



3. Click Add next to the Profile field.




4. Enter a profile name.




5. Click **OK**.
6. Set the profile settings as needed.
7. Click **Save**. Your new profile now appears in the profile list.

#### To edit an existing profile:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Click the drop-down button from the **Profile** field.



4. Select the desired profile from the list of profile names displayed.
5. Make the desired changes to the selected profile.
6. Click Save.

#### To change to a different profile:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**. The different profiles will be displayed.


3. Click the profile that you will use during your current shift or work hours. You should now see the change to the behavior as defined in the selected profile.

## Setting Your Display Mode


One of the main preferences that you may want to change is the **Display Mode**. You can use this preference to change the report display from landscape mode to portrait mode or the reverse. This will allow you to support the monitor setup that may be different in the different locations that you work.

### To set your Display Mode:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, from the **Display Mode** list, select the mode that you desire.




4. Click **Save**.

## Setting Your Default Page


When you log in to InteleOrchestrator, the application opens on the page that is set as your default page. You can choose the **Home** page or any other page to which you have access. Most radiologists select **InteleOrchestrator** as their Default Page to save having to make an extra click at logon.

**To set your default page:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, from the **Default Page** list, select the page you want as the default.



4. Click **Save**.


## Filtering Your Reading Queue

Use the **Filter My Reading Queue** to exclude exams that have been assigned to another user or to groups that you are not a member of. It's a way to tidy up the worklist a little bit to keep those exams you don't need to see from displaying.


-  Most practices automatically remove all exams that are assigned to other users from your My Reading Queue.

**To set My Reading Queue to show or hide exams assigned to others:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, enable **Filter My Reading Queue**.




4. Click **Save**.

## Including All Shifts On Your Worklist


If your practice uses a shift-based workflow, you may be assigned to multiple shifts. If so, InteleOrchestrator prompts you to select your primary shift when you log in. However, you can choose to include all worklists from each shift in your My Reading Queue to see all exams from all shifts that you are reading for that day.

### To include all shift worklists in your My Reading Queue:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, enable **Include all shifts in Work Queues..**




4. Click **Save**.

## Exam Double Click


This setting determines the action that InteleOrchestrator will take (if anything) when a user double-clicks on a row within the Worklist/Advanced Search interface. The default behavior is to launch PACS and dictation while entering the patient view.

### To set exam double click:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, navigate to the **Exam double click** and choose the desired option from the drop-down menu.



4. Click **Save**.


## Start Viewer

By checking this box the default PACS system will automatically display when the user launches the Worklist page.


Not all PACS systems can be configured to support this.

To enable the start viewer:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, navigate to the **Start Viewer** and choose the desired option from the drop-down menu.



4. Click **Save**.


## Start Dictation

By checking this box the dictation/voice recognition system will launch based on the selected options when the Worklist page is launched.


Not all dictation systems can be configured to do this.

**To enable the start dictation:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, navigate to the **Start Dictation** and choose the desired option from the drop-down menu.



4. Click **Save**.

## Enabling Combinable or Consecutive Reports


You can configure InteleOrchestrator to prompt you when there are patient exams that could be combined or read consecutively. For example, if you would like to dictate a patient's CT Chest exam and their CT Abdomen/Pelvis exam on the same report, you can configure InteleOrchestrator to prompt you to combine the reports. Or you can configure InteleOrchestrator to prompt you when there are reports that can be read

consecutively. For example, your patient may have a screening mammogram that appears in your worklist because of your subspecialty, but also have a breast ultrasound that does not appear in your worklist. If you enable the Consecutive Reports prompt, after you read the screening mammogram, InteleOrchestrator will ask you if you want to read your patient's ultrasound next. This helps ensure all exams for a patient are dictated at the same time.

See [Combining Reports and Reading Reports Consecutively](#) to see how to interact with this prompt while reading exams.


**To enable the prompts for combinable or consecutive reports:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.




3. Under **Preferences**, do the following:

- Enable Check for Combinable Reports to configure InteleOrchestrator to prompt you to combine the reports.
- Enable Check for Consecutive Reports to configure InteleOrchestrator to prompt you when there are reports that can be read consecutively.




Once **Check for Combinable Reports** is enabled, additional fields such as the **Filter Combinable Reports by**, **Filter by Matching** dropdown and the **Match Event ID** checkbox will appear.



Functionality of each field is as follows:

- **Filter Combinable Reports by:** You can filter the exams to show only exams tagged with specified subspecialty using this field. This might be useful for a subspecialty radiologist who only reads or combines exams within their subspecialty.
- **Filter by Matching:** This field enables you to filter the list to show exams that precisely match the selected criteria, providing a more precise and targeted filtering approach. You can filter the exam list that match the selected criteria such as modality, subspecialty, body part, accession, and site accession.

- **Match Event ID:** This field will find any exams that match the Event ID of the selected exam.




4. Click **Save**.

## Selecting Sounds


InteleOrchestrator can play sounds when you receive an alert, notification, message, or collaboration in InteleOrchestrator. You can customize or disable the sound for each notification individually.

**To select a sound:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Select a sound option from the corresponding list in the **Preferences** section (alert, notification, reminder) or **Chat Preferences** section (Messaging).



💡 Click **Play** ⏪ to listen to a sample of the currently selected sound.

4. Click **Save**.


## Setting Search Result Page Size

When the Advanced Search functionality is used, InteleOrchestrator will limit the results to a maximum number of exams. Changing this value will set the maximum number of exams returned. This will not have any impact when doing a Count-type search or when doing an Export.


 Setting this value to a high number may have performance consequences, both on the user's local system and on the server (potentially impacting all users)

### To set up a search result page size:

1. Click your username in the upper-right corner of the InteleOrchestrator interface.
2. Click Profile Management.



3. Under **Preferences**, in the **Search Result Page Size** field, enter the number of exams to display on a single page.




4. Click **Save**.

## Worklist Max Show Count


In the Worklist tree on the Worklist page a user can select an option to show a count of how many exams are in a given worklist (without having to actually load the worklist). For performance reasons, the user is normally limited to only being able to turn on Show Count to a few worklists. Adjusting this value will set how many worklists can be set to Show Count.

To set up a worklist max show count:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, in the **Worklist Max Show Count** field, enter the number of exams to display on a single page.



💡 To avoid scrolling through an overly long list, a maximum of 200 exams is recommended.

4. Click **Save**.


## Open Peer Review in Window

Users can configure their profile to open the peer review window when initiating a peer review from a patient view. For more information about using Peer Review, see [About Peer Review Workflow](#).


💡 Popup blocker needs to be turned off in order to use this feature.

**To open peer review in window:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, enable **Open Peer Review in Window**.




4. Click **Save**.

## Changing the InteleOrchestrator Theme


You can change the color theme of the InteleOrchestrator interface to suit your working environment. You can choose between two color themes: Dark (default), or Classic.

**To edit the InteleOrchestrator theme:**

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Preferences**, from the **Theme** list, select the theme you want to use:



- Dark Theme: Darker backgrounds and lighter text.
- Classic Theme: Blue-gray backgrounds with dark text.


4. Click **Save**
5. Log out of InteleOrchestrator, and then log back in to apply the new theme. See "[Logging in to InteleOrchestrator](#)".

## Setting the Default Image Viewer and Dictation System


If your practice uses multiple image viewers or dictation systems, you can set the default image viewer and dictation system you want to use, or you can set the viewer and dictation systems individually for each site that you read for in InteleOrchestrator.

### To set the default image viewer and dictation system for all sites:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.




3. Under **Application Preferences**, select an option from the **Default Viewer** or **Default Dictation** list for the corresponding site.



4. Click the pencil icon to apply the selected preference to all sites.
5. Click **Save**.

#### To set the default image viewer and dictation system for a single site:

1. Click your user name in the upper-right corner of the InteleOrchestrator interface.
2. Click **Profile Management**.



3. Under **Application Preferences**, click the site in which you want to change the viewer or dictation system.

The screenshot shows the 'Profile Management' application window. At the top, it displays 'Profile: Default (active)' and 'Location:'. Below this, under 'Application Preferences', there is a checked checkbox 'Show only editable site application preferences'. Two dropdown menus are shown: '<Default Viewer>' and '<Default Dictation>'. A navigation bar below them includes icons for back, forward, search, and page numbers (Page 1 of 1). To the right, it says 'Displaying 1 - 13 of 13'. A table lists sites with their current viewer and dictation settings. The 'Fuji' row is selected, and its dropdown menu is open, showing options like 'Reporting' and 'Powerscribe 360'. Buttons for 'Update' and 'Cancel' are visible at the bottom of the dropdown.

Site	Default Viewer	Default Dictation
AIdoc		Reporting
CuraCloud	InteleViewer workflow4161	Reporting
Epic		Reporting
Fuji		Reporting
Hologic		Reporting
Imaging Center 2		Powerscribe 360
Imaging Center 3	InteleViewer workflow4181	Reporting
Lexmark		Reporting
McKesson		
McKesson2		

4. In the drop-down list of options, select the viewer or dictation system you want to use for this site.
5. Click **Update**.
6. Click **Save**.

# 3

# WORKLISTS AND EXAMS

InteleOrchestrator uses your practice's business logic to create worklists customized for you. Using these worklists, you can read exams, dictate reports, collaborate with fellow radiologists, and much more.

## In this chapter:


Managing My Reading Queue .....	91
Opening Exams .....	95
Managing Exam Assignments .....	96
Using Auto-Next .....	99
Using the Auto-Next Optimizer .....	106
Creating Your Own Worklist .....	107
Customizing the Work Area with Worklists .....	109
Searching in InteleOrchestrator .....	112
Using Folders .....	118

# Managing My Reading Queue

The My Reading Queue section of InteleOrchestrator is considered your master worklist. If you are using a shift-based workflow, you will see all worklists associated with a shift under your **My Reading Queue**. These worklists will automatically populate as soon as your shift is selected.

Expand **My Worklist** to locate **My Reading Queue** in the worklist panel.

- If your group is not using a shift-based workflow and you have the permissions, see [Manually Creating the My Reading Queue](#).



The ordering of the worklists within My Reading Queue is important, as this is a combined or stacked worklist based on the shifts you are working.


When you begin your shift, you should start at the top of the My Reading Queue and work down through the worklist.

## Customizing My Reading Queue


If you desire a different ordering of My Reading Queue and you have the permission to do so, you can drag and drop the worklists to set the reading priority. For example, you may want to read all the exams in the **Epic Exams** worklist first.

## To reorder My Reading Queue worklists:

1. Click and hold the worklist that you wish to move.



2. Drag the selected worklist to the new location within the My Reading Queue.




3. The changed worklist order will now take precedence over the default worklist order based on shift. This worklist order will not be saved automatically. You will see the default worklist order the next time you log on.

## Manually Creating the My Reading Queue


If your group is not using a shift-based workflow and you have the permissions, you can create your My Reading Queue manually.

### To manually create your My Reading Queue:


1. In the Worklist column, click the **Create Queue** icon.



2. Select the **Create a Reading Queue** from the drop-down menu.




3. Select the worklists that you would like to add to My Reading Queue.




4. Click **Create**.

#### To edit My Reading Queue after it is created manually:

1. Right-click **My Reading Queue**.
2. Select Edit from the drop-down menu.



3. Make the desired changes to My Reading Queue.



4. Click **Create**.

## Opening Exams

You can open an exam from the Worklist area. You can choose to open the exam for dictation only, to view exam images only, or both.

 You can use Auto-Next to automatically launch the next exam in your worklist when you finish dictating an exam. The exam will open in view and dictation mode. See "["Using Auto-Next" \(page 99\)](#)".

### To open an exam:

1. Access the desired worklist.
2. In the worklist area, find the exam you want to view and do one of the following:
  - In the Action Column for the required exam, you can single click on any of the below icons.

- Click to open the exam for dictation mode only
- Click to view the images only (view mode)
- Click for both dictation mode and view mode.
- Click to hide the exam from launching auto next. It acts as a way to pre-skip exams before you start auto next.
- If you do not want to use one of the single click actions icons, you can Double-click the exam to launch in dictation and view mode.

## Managing Exam Assignments

If you have the permissions to do so, you can assign exams to yourself or other users and groups in InteleOrchestrator.

When you assign an exam to yourself, you are locking it for yourself. Depending on how your group worklists are configured, this may remove that exam from everyone else's worklist.

### To assign an exam to yourself:

- In the Assignment column of the worklist, click the padlock icon for the exam you want to assign to yourself.
- Starting in InteleOrchestrator version 4.5, when you perform an Assign/Unassign, Lock/Unlock, or Subscribe action on an exam that has been combined with other exams, that action will also apply to all exams within the combination group.

My Reading Queue Worklist (47 RVU:1.50 WU:46.00)										ANDERSON, CALEIGH
Site	Accession	Actions			Assignment		Status			
<b>Unread General</b>										
Imaging Cente...	GKH63034								Un	
Imaging Center 1	HXJ69516						Lock Exam to Me			
Imaging Center 1	IMSQA29444								Un	
Imaging Center 1	IMSQA29618								Un	

You can unassign an exam that is assigned to you by clicking the padlock icon again.

### To assign an exam to someone else:

1. In the Assignment column of the worklist, click the Assign/Lock icon .
2. In the Assign dialog, select a **Radiologist**, **Group**, or **User** to assign the exam to.  
 This dialog also shows the assignment history for the exam.

**Assign**

Radiologist:	Group:	User:
Administrator, Initial		
<input type="checkbox"/> Make assignment persistent between login sessions		
Reason: Dictation		

**Assignment History**

Date	Assigned	Assignee	Reason	Persistent
09/21/2023 6:08 AM	-	Tobey, Amy	Auto assignment: Return from dictation	
09/21/2023 6:08 AM	Tobey, Amy	Tobey, Amy	Auto assignment: Dictation	Yes
09/11/2023 11:02 PM	-	Administrator, Initial	Auto assignment: Return from dictation	
09/11/2023	Administrator, Initial	Administrator, Initial	Auto assignment:	Yes

**Buttons:** Cancel, Unassign, Assign

3. If required, enter the reason you're assigning the exam.
4. Click **Assign**.


-  Starting in InteleOrchestrator version 4.5, when you perform an Assign/Unassign, Lock/Unlock, or Subscribe action on an exam that has been combined with other exams, that action will also apply to all exams within the combination group.

## Managing Exams Assigned to You when Logging Out


If you log out of InteleOrchestrator while you have exams still assigned to you, the Exam Assignment dialog will be displayed, prompting you to either keep your assigned exams or release your assigned exams back to the reading queue. The released exams can then be reassigned to other Radiologists.

### To keep your exams assigned to you when you log out:

1. Click **Logout** in the top-right corner of the InteleOrchestrator interface.



2. In the **Exam Assignment** dialog, leave the check-marks enabled next to the exams you want to retain.



3. At the bottom of the Exam Assignment dialog, click **Logout**.

## To release some of your exams to the reading queue when you log out:

1. Click **Logout** in the top-right corner of the InteleOrchestrator interface.
2. In the Exam Assignment dialog, deselect the check-marks next to the exams you want to return to the reading queue.
3. At the bottom of the Exam Assignment dialog, click **Logout**.

## To release all of your exams to the reading queue when you log out:

1. Click **Logout** in the top-right corner of the InteleOrchestrator interface.
2. At the top of the Exam Assignment dialog, click **Unassign All**.
3. At the bottom of the Exam Assignment dialog, click **Logout**.


# Using Auto-Next

InteleOrchestrator uses your practice's business logic to distribute exams to radiologists, based on their roles and specialties, prioritizing those exams for each user. The Auto-Next feature allows you to save time by automatically launching the next exam on your worklist when you finish dictating an exam.

## To start Auto-Next:

- Click the Auto-Next  icon above My Reading Queue.

This immediately launches the first exam on your worklist in your dictation and viewer applications. A mini version of your worklist, your Auto-Next queue, appears in the Worklist tab and the Patient View window for the current patient. For more details on the Patient View, see "[Understanding the Patient View](#)" (page 25).




As soon as you sign the dictation for the current exam, the next exam in your auto-next queue opens.

## To skip exams in your Auto-Next queue:

The radiologist can temporarily skip exams in the queue based on urgency or their preference not to report them. Radiologists may prioritize certain types of exams based on factors such as the complexity of the study, the need for consultation, or technical issues.

1. In the Worklist tab, select the exam you want to skip in your Auto-Next queue.
2. Click the Skip  icon for that exam. The exam will still appear in your Auto-Next queue, but Auto-Next will not open it for you.

 To return to the exam or "unskip" it, click the Skip icon again.



The screenshot shows the Radiologist Worklist Worklist interface. At the top, there's a navigation bar with Home, Worklist, and Scheduling. Below that is a toolbar with Auto-Next, Folder, and various icons. The main area is divided into two sections: 'Auto-Next' and 'Radiologist Worklist'. The 'Auto-Next' section shows a queue of exams with progress bars and current/suggested values. The 'Radiologist Worklist' section shows a list of exams with columns for Hospital, Patient Name, Status, and Type. One row in the worklist has a red box around its status column, which is currently set to 'Skip'. To the right of the worklist, there's a sidebar with patient details: Jinn, Qui-Gon (Female, MZ009966, DOB 12/20/1946, Age 77y, MPI 10339, Site MRN MZ009966 (Hospital A)).

## To set an exam to be next in your Auto-Next queue:

In situations where a radiologist receives a call from a referring provider requesting an urgent review or if a technologist asks them to prioritize an exam, marking it as 'Next' automatically moves urgent or critical exams to the top of the worklist. The case that is marked as "Next" will launch as soon as the radiologist is finished with their current exam.

1. In the Worklist tab, select the exam you want to read next in your Auto-Next queue.
2. Click the Next  icon for that exam. When you indicate that an exam should be next, that exam is assigned/locked to you and will load as soon as you finish dictating the current exam.

The screenshot shows the Radiologist Worklist Worklist interface. At the top, there are tabs for Home, Worklist (selected), and Scheduling. Below the tabs is a toolbar with icons for Auto-Next, Folder, Worklist (8 RVU:9), and various controls. The main area is divided into two sections: 'Auto-Next' on the left and 'Radiologist Worklist' on the right. The 'Auto-Next' section displays a list of exams with icons, patient names, and status. The 'Radiologist Worklist' section shows a detailed view for a patient named Jinn, Qui-Gon, with fields for Gender (Female), MRN (MZ009966), DOB (12/20/1946), Age (77y), MPI (10339), and Site MRN (MZ009966 (Hospital A)).


## To add an exam to your Auto-Next queue:

If you want to read an exam that does not appear in the Auto-Next queue, you can add that exam to the queue.

1. In the Worklist tab, click the Add **+** icon in the Auto-Next toolbar.

This screenshot is similar to the previous one, showing the Radiologist Worklist Worklist interface. The 'Auto-Next' tab is selected. A red box highlights the 'Add Exam to Auto-Next Worklist' button in the toolbar. A tooltip for this button is visible, showing the text 'Add Exam to Auto-Next Worklist'. The rest of the interface, including the patient details and the 'Radiologist Worklist' section, is identical to the first screenshot.

2. A new Quick Search window will appear on your screen, where you can search for the exam you want to add to your queue.



3. Search for the exam according to the fields and click the **Search** button.
4. Exams will appear in the search results. Click the lock icon to assign the exam to yourself.
5. After making your selection, the lock icon you clicked will turn orange to indicate that the exam has been added to your queue.

The screenshot shows the same 'Quick Search' window with search results. The results table includes columns for Priority, Accession, Procedure, Patient Name, Exam Time, Status, Site, and Actions. One of the rows in the Actions column contains a lock icon, which is highlighted with a red square. The data in the table is as follows:

Priority	Accession	Procedure	Patient Name	Exam Time	Status	Site	Actions
Stat	CT23905	CT HEAD WITHOUT	Roberts, Sta...	10/03/2021 ...	Final	Imaging Cen...	
Stat	DX7150	ANKLE RIGHT COMPLETE	Richards, Mo...	10/03/2021 ...	Final	Metro Hospital	
Routine	DX38547	Mammogram, Screening	Bean, April	09/03/2021 ...	Final	Imaging Cen...	
Stat	US5896	KNEE RIGHT 4 VIEWS	Diamond, John	10/03/2021 ...	Final	Imaging Cen...	
Routine	NM22939	C-arm Arthrogram Shoulder	Vang, Lan	10/03/2021 ...	Final	Metro Hospital	

6. To close the window, click on the 'x' button located in the upper-right corner of the pop-out window.

**To rearrange your Auto-Next columns:**

Your InteleOrchestrator system administrator enables the columns in Auto-Next, but you can rearrange them as required.

1. Click the gear icon above your Auto-Next queue.

The screenshot shows the Auto-Next interface with a toolbar at the top. A gear icon is highlighted with a red arrow pointing to it, labeled "Columns Sort Gear". Below the toolbar is a section titled "Worklist (43 RV)" containing two items: "RVU" and "WU", each with current and suggested values. Underneath is a section titled "Unread Exams" with three entries related to chest exams.

2. Once you click the gear icon, all the column headers for Auto-Next will appear.

The screenshot shows the Auto-Next interface after clicking the gear icon. The gear icon is now highlighted with a red box. Below the toolbar, a row of column headers is displayed, including "Body Part", "Time Req", "Assigned", "Acquisit", "Comb", and "Patient Na", all enclosed in a red box. The rest of the interface remains the same, showing the worklist and unread exams sections.

3. Click the drop-down arrow on one of the columns and then click **Columns**.

The screenshot shows the Auto-Next software interface. At the top, there's a toolbar with 'Auto-Next' and 'Folder' buttons, followed by a 'Worklist (8 RVU:0.16 WU:7.1)' header and various control icons. To the right, patient information is displayed: Segura, DOB: [redacted], Gender: [redacted], Age: [redacted], MPI: [redacted], MRN: [redacted], SSN: [redacted]. Below this are 'Report' and 'Protocol' buttons.


The main area is a 'Worklist' grid. It has columns for 'Patient Name', 'Body Part' (which is highlighted with a red box), 'Acquisition Ti', 'Time Remai...', 'Combined', and 'MRN'. A dropdown menu above the grid shows 'Sort Ascending' and 'Sort Descending' options. A 'Columns' button is also present, which opens a context menu on the right side of the screen. This menu lists various data items with checkboxes: Subspecialty, Status, Speech Engine, Site Procedure Code, Site Procedure, Signed Time, Priority, Modality, Location, Last Modified, Assignment (which is selected and highlighted in blue), Patient Name, Body Part, Acquisition Time, Time Remaining, Combined, and Site MRN. At the bottom of the context menu, there are 'FILE UPLOAD' and 'Insurance & Document' buttons.

4. Enable or disable a column type to show or hide the column, respectively.
5. Click and drag a column to its new position.


#### To exit Auto-Next:

When you exit Auto-Next, the main worklist page appears and none of your skipped exams are retained. If you pause Auto-Next, the main worklist appears and your skipped exams are retained for you when you resume Auto-Next.

1. Discard the exam that you are currently reading in your dictation system.




2. In the Discard Report dialog, click **Discard Report**.



3. In the Dictation Cancelled dialog, select whether you want to move to the next exam in your queue, pause Auto-Next, or Exit Auto-Next. Click one of the following:

- **Next Exam:** discards the current exam and moves to the next exam in your Auto-Next queue.
- **Pause:** Pauses Auto-Next, leaving the current report open on your screen.
- **Exit:** Exits Auto-Next and returns you to your My Reading Queue.

💡 You can exit or pause Auto-Next by clicking the **X** or **II** button above your Auto-Next queue.



## Using the Auto-Next Optimizer

You can set your Auto-Next queue to prioritize and automatically launch exams based on modality, body part, or application.

For example, if you prioritize by modality and you are currently reading a CT exam, Auto-Next launches the next CT exam in your worklist.


The feature does not re-sort your Auto-Next queue, it simply launches the next exam based on the efficiency settings you selected.

Optimizing by application applies only if your group uses multiple dictation systems or image viewers. For example, if you launch InteleViewer for some exams and Fuji for other exams, Auto-Next launches all InteleViewer exams before launching any Fuji exams.

STAT exams override any optimization settings that you have enabled.

### To set and enable Auto-Next optimizer:

1. In the Auto-Next toolbar, click the drop-down menu next to the Optimize Efficiency icon.



2. Choose an optimization method:


- By one sorting method only: Select the sorting method (Modality, Body Part, Application) in all three sorting options.
  - By multiple sorting methods: Select the sorting method to apply first in the **1st:** sorting option, and the second sorting method to apply as the **2nd:** sorting option, and the third sorting method as the **3rd:** sorting option. For example, sorting by Application, Modality, then by Body Part.
3. Click the Optimize Efficiency  icon to enable the optimizer. This icon turns green when Auto-Next Optimizer is enabled.

## Creating Your Own Worklist


If you have the permissions to do so, you can create your own worklist. Worklists that you create appear for you only.

## To create your own worklist:

1. Click the Advanced Search  icon.



2. Enter search criteria for the exams you want to add to this worklist. For example, if you want to find all unread exams from Community Hospital, select Community Hospital in the **Site** field and Unread in the **Status** field.
3. Click **Search**.
-  If you do not see the **Add to Worklist** button now, you do not have permission to create your own worklist.
4. Click **Add to Worklist** and then select **My Worklist** from the list.



5. In the Worklist Properties dialog, enter a **Name** for the worklist, and then select any other desired properties.

**Add to My Worklist Worklist**

**Options** Active Time 24/7

Name:	Test		
Refresh rate in seconds:	60		
Row coloring schema:	Status		
First sort order:	Time Remaining		
Second sort order:	Priority		
Third sort order:			
Threshold Counts:	<Lower>	<Upper>	
Active style for:	<input checked="" type="checkbox"/> Priority	<input type="checkbox"/> Status	<input type="checkbox"/> Time Remaining
Exclude from Auto Next:	<input type="checkbox"/>		
Wrap content:	<input type="checkbox"/>		

\*The limit on the number of worklists with activations has been reached.

**Cancel** **Clear** **Save**


- Click **Save**. The worklist appears under **My Worklists** in the Worklist tab.

## Customizing the Work Area with Worklists

You can customize the columns for the worklists in the work area and resize some of the sections. The changes that you make apply to all your worklists.

### To add or exclude columns in your worklist:

- Ensure that a worklist is currently displayed.
- Click the drop-down arrow on one of the columns and then click **Columns**.



3. Enable or disable a column type to show or hide the column, respectively.

	Accession	MRN	Actions
RA00021	MZ009966		
RA00019	MZ009966		
RA00020	MZ009966		
RA00056	MZ100036		
RA00089	MZ100057		
RA00076	MZ100050		
RA00067	MZ100043		
RA00111	MZ100071		

4. Click **Logout** to save your changes.

Your changes will persist during the current session but will not reappear unless you log out from InteleOrchestrator by using the Logout button. Clicking the (x) on your browser will close InteleOrchestrator without saving your changes

#### To reorder the columns in your worklist:

1. Ensure that a worklist is currently displayed.
2. Click and drag a column to the desired location.
3. Click **Logout** to save your changes.

Your changes will persist during the current session but will not reappear unless you log out from InteleOrchestrator by using the Logout button. Clicking the (x) on your browser will close InteleOrchestrator without saving your changes

#### To change the width of columns in your worklist:


1. Ensure that a worklist is currently displayed.
2. Hover the mouse over a column border until your mouse pointer changes to a dual arrow.
3. Click and drag the column to the desired width.
4. Click **Logout** to save your changes.

Your changes will persist during the current session but will not reappear unless you log out from InteleOrchestrator by using the Logout button. Clicking the (x) on your browser will close InteleOrchestrator without saving your changes

 You can double-click the column to resize it to the maximum width of the contents of the column.


### To resize the worklist panel:

1. Place the cursor on the worklist column boundary, that you wish to resize, until the cursor changes to .



2. Click and drag the boundary to increase or decrease the width of the worklist column.

 You can also click the double arrows  at the top of the worklist menu to completely collapse the worklist panel. Click the double arrows again to expand the worklist panel.



- Click **Logout** to save your changes.

**Tip:** Your changes will persist during the current session but will not reappear unless you log out from InteleOrchestrator by using the Logout button. Clicking the (x) on your browser will close InteleOrchestrator without saving your changes

## Searching in InteleOrchestrator


There are four search methods in InteleOrchestrator, and the one you choose to use depends on the exam you want to search for. The following table explains each search functionality, when you would want to use it, and links to instructions for how to conduct each type.

Search Type	Reason for Use	Instructions
Quick Search	When you know the name, procedure, accession number, site, and/or MRN of an active exam you want to find.	<a href="#">"Quick Search" (page 113)</a>
History Search	When you want to search for an exam you have recently opened or read.	<a href="#">"Searching in Exam History" (page 113)</a>
Worklist Search	When you want to search for an exam that is in your worklist	<a href="#">"Searching In A Worklist" (page</a>

Search Type	Reason for Use	Instructions
		<a href="#">115)</a>
Advanced Search	When you want to search for an exam that is outside your worklist.	<a href="#">"Advanced Search" (page 116)</a>

## Quick Search

You can use the Quick Search to find an exam in your worklist with a few key pieces of data. The results of your search appear in the work area.



 You can search for a partial string when using any of the search functions that are available in InteleOrchestrator. To do so, add an asterisk (\*) to the end of the partial string. For example, to search for all patients with a last name that start with the letters "Clar", enter "Clar\*" in the last name search field.

### To perform a quick search:

1. In the Search section of the worklist column, enter a patient's name (last and/or first), the procedure, the Accession number, the Site and/or the MRN. For information on where the Search section is, see ["About the Worklist Column" \(page 21\)](#).
2. Press **ENTER** or click the search  icon.


## Searching in Exam History

You can search for an exam from the list of exams you have recently opened or read. The results of your search appear in the work area.

 You can search for a partial string when using any of the search functions that are available in InteleOrchestrator. To do so, add an asterisk (\*) to the end of the partial string. For example, to search for all patients with a last name that start with the letters "Clar", enter "Clar\*" in the last name search field.

### To search in your History:

1. In the History section of the worklist column, enter a patient's name (last and/or first), the procedure, the Accession number, the Site and/or the MRN.



The screenshot shows the 'History (39)' worklist interface. At the top, there are buttons for 'Today' and '1wk', followed by three icons: a clock, a gear, and an information sign. Below these are buttons for 'Filter' and a refresh symbol. A search bar contains the text 'segura'. The main table has columns for 'Accession', 'Procedure', 'Patient', 'Draft', and 'Edit'. The data rows are:

Accession	Procedure	Patient	Draft	Edit
FFI-TEST3	XR Foot ...	Segura, Dora		
FFI-TEST3		HAMILTON, ALICE		
CL65205		MOORE, JUSTIN		
MR17442	BILAT DI...	Hamilton, Dennise		


2. If needed, click Filter to enable additional filters.

- Enable the **Needs Signature** check box to find exams you must sign.
- Enable the **Waiting Communication** check box to find exams that contain a note that is waiting for some type of action.

3. Press **ENTER**.

	Accession	Procedure	Patient	Draft	E
	FFI-TEST3	XR Foot ...	Segura, Dora		

## Searching In A Worklist



You can search for exams in the currently displayed worklist by using the Worklist Search. The results of your search appear in the work area.

You can search for a partial string when using any of the search functions that are available in InteleOrchestrator. To do so, add an asterisk (\*) to the end of the partial string. For example, to search for all patients with a last name that start with the letters "Clar", enter "Clar\*" in the last name search field.

### To search for exams in the current worklist:

1. Ensure that a worklist is currently displayed.
2. In the Worklist Search field, enter the search criteria to search by accession, procedure, or patient name.

The screenshot shows a table of patient records under the heading "Radiologist Worklist (8 RVU:9.50 WU:14.00)". The search bar at the top right contains the name "Jinn". The table columns include Status, Patient Name, Accession, MRN, Actions, Assignment, Modality, Site, Assigned, Radiologist, and Procedure. Several rows are highlighted in red, corresponding to the search results.

Status	Patient Name	Accession	MRN	Actions	Assignment	Modality	Site	Assigned	Radiologist	Procedure
Unread	Jinn, Qui-Gon	RA00021	MZ009966			XR	Hospital A			
Unread	Jinn, Qui-Gon	RA00019	MZ009966			XR	Hospital A			
Unread	Jinn, Qui-Gon	RA00020	MZ009966			XR	Hospital A			
Unread	Hutt, Jabba the	RA00056	MZ100036			XR	Hospital A			
Unread	Calrissian, Lando H	RA00089	MZ100057			MR	Hospital A			
Unread	Solo, Han	RA00076	MZ100050			XR	Hospital A			
Unread	Vader, Darth	RA00067	MZ100043			XR	Hospital A			DX.AAA0
Unread	Leia, Princess H	RA00111	MZ100071			XR	Hospital A			CF456

### 3. Press ENTER.

This screenshot shows the same worklist interface as above, but only the first row for "Jinn, Qui-Gon" is highlighted in red, indicating it is selected or the result of the search.

## Advanced Search

You can perform an advanced search to find an exam outside your worklist or with advanced search criteria.

- 💡 You can search for a partial string when using any of the search functions that are available in InteleOrchestrator. To do so, add an asterisk (\*) to the end of the partial string. For example, to search for all patients with a last name that start with the letters "Clar", enter "Clar\*" in the last name search field.

### To perform an advanced search:

1. Click the Advanced Search icon at the top of the InteleOrchestrator interface.


The Advanced Search dialog box is open, showing various search criteria sections. Key sections include:
 

- Report Section:** Includes fields for Section and Keywords.
- Enterprise:** Includes fields for Site, Location, and Site Group.
- Exam Date:** Includes fields for Exam Time, Dates Relative to Today, From, To, and Radiologist.
- Exclude Incomplete Notes:** Includes fields for Communication, ED Prelim, and Follow Up Notes.
- Report:** Includes fields for Patient, Personnel, and Exam Note.
- Exam:** Includes fields for Accession, Site Accession, Reason, Priority, Status, Modality, Type, Emergency, Billing Status, and Subspecialty.
- Patient:** Includes fields for MRN, MPI, Merge Status, First Name, Last Name, Gender, Age, Inconsistency, and Insurance.
- Personnel:** Includes fields for Radiologist, Resident, Resident Radiologist, Attending, Radiologist, Ordering, Ordering Group, Assigned (U), and Assigned (G).
- Exam Note:** Includes fields for Keywords, Created By, Type, and Created On.
- Patient Note:** Includes fields for Keywords, Created By, Type, and Created On.

2. In the Advanced Search dialog, enter the data you want to search for from the exams, such as Patient Note or Subspecialty.

-  If you do not add enough parameters to make a manageable search, InteleOrchestrator automatically limits the date range to within 3 days.

You can change the date range filters to search your data within a specified timeframe of your choosing. You can define a custom date range by choosing the start and end dates. Or, you can select from a list of predefined ranges, such as, Today, 3 days, 1 week, 1 month and 1 year.



### 3. Click **Search**.

 You can click **Count** to see how many results your search criteria returned.

### 4. The results of your search appear in the work area.


Referring Phy	Actions	Age at Exam	Site	Ordered Time	Begun Time	Final Time	Status	Prelim Time	Body Part	Signed Time	Communication	Radiologist	Comments
ALEX LEE		11y	Imaging Center 1	1992-10-...	2022-09-...		Dictating	<input type="checkbox"/>		2022-09-22 3...		Adminis...	
ALEX LEE		11y	Imaging Center 1	1992-10-...	2022-09-...		Final Addendum	<input checked="" type="checkbox"/>	abdomen,...	2022-09-22 3...		Adminis...	
ALEX LEE		11y	Imaging Center 1	1992-10-...	2022-01-...		Dictating	<input checked="" type="checkbox"/>	2021-11-...	abdomen	2021-11-05 9...	Unknow...	
Quin Knapp		0d	Unknown	2023-01-...			Dictating	<input checked="" type="checkbox"/>		2023-04-17 3...		Adminis...	
JOHN LEE		27y	Imaging Center 1	2021-09-...			Draft	<input checked="" type="checkbox"/>	lower arm...			Herbert...	
NICHOLA...		75y	Imaging Center 1	2021-08-...	2022-09-...		Final	<input checked="" type="checkbox"/>	chest	2022-08-05 9...		Hartma...	
REFER DE...		91y	Unknown				Draft	<input checked="" type="checkbox"/>				Adminis...	

### 5. To minimize advanced search, click on the "x" below the gear.

## Type of Advanced Searches:


There are different types of advanced searches. You can save worklists for all search types. The search type you use is dependent on what you are searching for and how you want your worklist to function.

1. Exam - This search focuses on exam related search filters, such as subspecialty, priority, and status.
2. Communication - This search type searches for communication notes.
3. Task - This search type searches for task notes.
4. Peer Review - This search type searches for peer review notes.
5. Tech QA - This search type searches for tech QA notes.
6. Teaching Note - This search type searched for teaching notes.
7. Tech Note (anonymized) - This search type searches for teaching notes, but anonymizes the results.




## Using Folders

You can track exams in folders in InteleOrchestrator. This feature is commonly used to track interesting cases, cases for a tumor board, or cases for peer learning. The Folder tab is situated next to the Worklist tab.




## To create a folder:


1. In the Worklist column, click the **Folder** tab.
2. Click the Add icon **+**.




3. Enable the **Remove Patient Identifiers** check box if you need to share the contents of this folder with others.




4. Select the sorting and appearance Properties for the exams in the folder, as needed.



5. Apply the style formatting needed by enabling the **Priority**, **Status**, and/or **Time Remaining** check boxes.




6. If you need to share the contents of this folder with other users, click the **Share Options** tab.



7. Select the level of access you to assign to these users.


- **View:** users can see the folder and the cases in the folder.
- **Contribute:** users can add cases to the folder.
- **Admin:** users can add cases to the folder and edit any of the properties.




8. Click **Save**.

## To share a folder:


1. In the Worklist column, click the **Folder** tab.
2. In the list of folders, right-click the folder you want to share, and then click **Share & Property**.




3. In the Edit Folder dialog, click the **Share Options** tab.
4. Select the users or groups with whom you want to share this folder.
5. Select the level of access you want to assign to these users.



6. Click Add  to add this permission set to the folder.




7. Click **Save**.




#### To add an exam to a folder:

1. Open an exam or a worklist.
2. Do one of the following:

- Right-click an exam on your worklist and select **Add to Folder**.



- From the Notes section of the Patient View, click the Teaching Note icon .
3. In the Assign Exam to Folder dialog, enter the subject of the note in the **Title** field and then enter a message in the **Note** field that describes why you are adding this exam to the folder.
  4. Enable the check box for the folder you want to add the exam to.
  5. Click **Add**.



#### To remove an exam from a folder:

1. In the Worklist column, click the **Folder** tab.
2. Click the folder that contains the exam you want to remove.
3. In the list of exams that opens, click the trash can icon next to the exam you want to remove.

#### To open an exam from a folder:

1. In the Worklist column, click the **Folder** tab.
2. Click the folder that contains the exam you want open.
3. Click the Launch in Viewer icon.

#### To edit a teaching note for an exam in a folder:

1. In the Worklist column, click the **Folder** tab.
2. Click the folder that contains the exam you want edit.
3. Click the pencil icon.

# 4

## COLLABORATION AND COMMUNICATION

InteleOrchestrator has several built-in features that allow you to communicate with other radiologists, collaborate with other radiologists on exams, and help administrators to broadcast information to the practice.

### In this chapter:

Collaborating With Other Radiologists .....	127
Sending Notifications and Alerts .....	127
Adding Notes to Exams .....	131
Assigning a Note to the Operations/Support Staff .....	140

# Collaborating With Other Radiologists

 The Collaboration tool has been deprecated for InteleOrchestrator version 4.4 and later. Please see [About Messaging](#) for more information about the replacement for Collaboration.

Collaboration is a tool, available in InteleOrchestrator version 4.3 or earlier, for radiologists to communicate with one another about a specific exam. This feature is only available after you open an exam. When starting or accepting a collaboration, both collaborating radiologists receive a pop-up window that contains information related to the exam. From this window, radiologists can see patient and exam information, preliminary or final reports, and any notes associated with the exam. Collaborating radiologists can chat with one another and the accepting radiologist can easily launch the images.

 You must turn off your web browser's pop-up blocker for the collaboration window to open.

## To start a collaboration:

1. Open the exam.
2. Click the Collaboration button in the Communication taskbar.
3. Select the **User** with whom you want to collaborate.

 This list is only populated with users who are currently online and using either the worklist or in a patient view.

4. Click **Submit**.

## To accept a collaboration:

- Click **Accept** when the collaboration invitation appears on your screen.

# Sending Notifications and Alerts

InteleOrchestrator allows users to send notifications and alerts to other InteleOrchestrator users for any type of communication you'd like to track, including critical results, critical tests, and standard call reports. These notifications can also track tasks that a radiologist needs completed before the radiologist is able to read an exam, such as prior reports or exam re-do. These types of notes could be used between the

radiologist and their support/operations team. See [My Support Queue](#) for how the support/operations team works with communications notifications.

The notification can be created by any InteleOrchestrator user as long as the user has the permission to create the notification and view the notification.


Notifications do not disrupt the radiologist's workflow and are often used for information that is not time sensitive. Notifications that you send to another user appear in the Notifications area on the toolbar.

An alert is a pop-up that the receiving user must acknowledge before resuming to work. As a result, alerts are typically used for time-sensitive information only. Alerts are stored in the lower-right corner of InteleOrchestrator.


 Users can access previous alerts from the Alerts area of the toolbar.

#### To create a notification or alert:

1. In the Communication taskbar, click **Notifications** or **Alerts**.




2. In the message dialog, click **Create**.




3. Enter the message.

A screenshot of a "Create Notifications" dialog box. It features a large text area labeled "<Message\*>" with a red border. Below it are fields for "Receivers:" (dropdowns for "User" and "Group") and "Delivery:" (date "14/02/2023" and time "12:30 AM"). A checkbox "Send only for currently logged in users." is present. At the bottom are "Cancel", "Clear", and "Create" buttons, with "Create" also highlighted with a red rectangle.


4. Select the users or groups that will receive the message.




5. Enter the date and time the message will be sent.



6. Enable **Send Only For Currently Logged In Users** if you want the message sent only to users who are logged in at the selected date and time.



7. Click **Create**.



## Adding Notes to Exams

You can add different types of notes to exams in InteleOrchestrator.

- The available note types depend on the note types your practice decided to deploy. Certain note types can trigger workflows that were configured by your organization. See your InteleOrchestrator administrator for more details.

You can add a note to an exam from:

- the worklist.
- the Patient View.

InteleOrchestrator offers the following types of notes for exams:

Note Type	Reason for Use	Icon
Communication note	Communicate different issues with exams. Used mostly to communicate between the radiologist and operations/support staff.  Example: Exam is missing images.	
Communication Log	Creates a record of issues with an exam, without assigning the follow-up action(s) to an individual or group.	
Patient Note	Communicate patient-specific information. The patient note will follow the patient throughout the InteleOrchestrator system.	
Exam Note	Communicate exam specific information. The exam note will follow the exam throughout the InteleOrchestrator system.  Example: Changing the priority of the exam for a particular reason.	
Peer Review Note	Peer Review specific notes. See <a href="#">Peer Review Workflow</a> for examples where peer review notes are used during the Peer Review workflow.	
Prelim Note	ED prelim notes or other notes from a preliminary read.	
Follow-Up Note	Communicate follow-up appointment needs, such as diagnostic breast follow-ups.	
File Upload	Records the filename and optional note, along with details of who uploaded the file and when.	
Teaching Note	Adds an exam to a Folder and prompts for a Title and a Note to describe why it is being added to that folder.	
Tech QA Note	Communicate to track technologist quality.	


**To add a note to an exam:**

1. Do one of the following:
  - In the worklist:
    - a. Right click on any exam on your worklist.
    - b. In the list that appears, hover the mouse over the **Add Note** option.
    - c. From the list that appears on the right, select the note type you would like


to enter.

	Read Exam	JING, NICHOLAS	SKELETA
	Open in Viewer	ON, CHRISTIAN	ARM
	Open in InteleViewer workflow4181	BER, CHRISTIAN	NECK SC
	Open in Dictation	WN, SAMIR	LEG
	Open in Reporting	TINEZ, SARAH	SHOULD
	Hide Study	DNALD, NICOLE	ABDOME
<input checked="" type="checkbox"/>	Assign	IA29444, IMSQA	PETCT-T0
	Add Note		
	Emergency Access		
	Assign to Peer Review		
	DICOM Route		
	DICOM Route and Read		
	Change Status		
	Change Priority		
	Change Attending Radiologist		
	Cancel Exam		
	Delete Exam	LIAM	HEAD SC
	Edit Subspecialty	HTS, ALICE	INJURED
	Add to Folder	ON, CHRISTIAN	HEAD DI
	Send Report	DNALD, NATALIE	BREAST
	Print Report	EUX, NICOLE	HEAD DI

- From the Patient View:
  - a. View or read an exam to access the Patient View
  - b. From the Notes section at the bottom of the Patient View, click the icon for the type of note you want to use.




2. For example, click on the Communication Note icon. Populate the fields in the note dialog that opens.
3. Click the **Type** dropdown menu and select the note type.




4. Assign to either a specific user(s) or to a group.

- Select the <User> drop-down to assign to specific users.



- Select the <Group> drop-down to the note to a group.



- Leave the note unassigned if you want to assign the note to the Operations team/support staff. See [Assigning a note to the Operations/Support Staff](#).
- ⌚ You can click the star next to a user name or group name for those contacts you use often. That will result in those favorites being listed at the top of the drop-down menu. Click the star to deselect the contact to no longer keep that user or group as a favorite in the Assign drop-down list. The user or group will then fall back to the original position in the list.
5. Add text to the message for example missing.
  6. Click **Submit**.

 **Communication Note**


Type: \* Issue with Exam

Assign: <User> Radiology Support Team

<Default Note Messages>

missing

**Cancel** **Clear** **Submit**



-  Indicating Public on a note means that the note is visible to anyone who can see the exam in a worklist. Indicating Private on a note means it is only visible to the users or groups indicated in the note.
-  If the type of note you want to create is not listed in the icons on top of the Notes section, click **Create** to access the other note types.

## To add a note to multiple exams from the worklist:

1. Select the check box next to the exams needing notes.

Radiologist Worklist Worklist (8 RVU:9.50 WU:14.00)							
Status	Patient Name	Accession	MRN	Actions	Assignment	Notes	Print
Unread	Jinn, Qui-Gon	RA00021	MZ009966	<input type="checkbox"/>			X
Unread	Jinn, Qui-Gon	RA00019	MZ009966	<input checked="" type="checkbox"/>			X
Unread	Jinn, Qui-Gon	RA00020	MZ009966	<input checked="" type="checkbox"/>			X
Unread	Hutt, Jabba the	RA00056	MZ100036	<input type="checkbox"/>			X
Unread	Calrissian, Lando H	RA00089	MZ100057	<input type="checkbox"/>			X
Unread	Solo, Han	RA00076	MZ100050	<input type="checkbox"/>			X
Unread	Vader, Darth	RA00067	MZ100043	<input type="checkbox"/>			X
Unread	Leia, Princess H	RA00111	MZ100071	<input type="checkbox"/>			X

2. Click on the Actions drop down.
3. In the list that appears, hover the mouse over the Add Note option.

4. From the list that appears on the right, select the note type you would like to enter.

The screenshot shows a software interface with a top navigation bar containing links: Home, Worklist (which is selected and highlighted in blue), Scheduling, Analytics, Messaging, Reporting, Management, Dashboard, CRM, and Help. Below this is a search bar with the placeholder ".00)". The main area displays a table titled "Worklist" with columns: MRN, Actions, Assignment, Modality, and Site. The table lists several patient records (e.g., MZ009966, MZ100036, MZ100057, MZ100050, MZ100043, MZ100071). For each record, there is a row of actions: Read Exam, Open in Viewer, Open in Dictation, Move to Reference, Hide Study, Add Note, Assign to Peer Review, DICOM Route, DICOM Route and Read, Change Status, Change Priority, Change Attending Radiologist, Cancel Exam, Delete Exam, Combine, Merge Exams, Edit Subspecialty, Add to Folder, Send Report, and Print Report. A context menu is open over the "Add Note" option for the first row, listing various note types: Communication Note (selected and highlighted in blue), Communication Log, Patient Note, Exam Note, ED Prelim Note, Follow-up Note, Tech QA Note, Teaching Note, and File Upload. The "Communication Note" option is currently active.

## Assigning a Note to the Operations/Support Staff

The previous section provided instructions for creating a note that can be assigned to specific groups or other users. You can leave the note unassigned, which would have the effect of sending the note to the whole operations team.

After the note is submitted, the exam containing the note will be removed from the radiologist's worklist and added to the My Support Queue. For more information about the My Support Queue, see [My Support Queue](#).

# 5

## INTELEORCHESTRATOR MESSAGING

InteleOrchestrator Messaging was created to offer improved performance, group management, message notification, and configurable permissions to users. With Messaging, you can quickly send secure messages to any user on the InteleOrchestrator platform - administrators, radiologist assistants, technologists, and other radiologists.

### **In this chapter:**

About InteleOrchestrator Messaging .....	142
Opening Messaging .....	142
Understanding the Messaging Page .....	143
Opening the Messaging App Window .....	144
Changing your Messaging Status .....	145
Setting Location Information .....	147
Searching for Messaging Users .....	147
Sending Messages .....	148
Group Management in InteleOrchestrator Messaging .....	153
Sending an Exam Link within Messaging .....	158

# About InteleOrchestrator Messaging

In InteleOrchestrator versions 4.4 and above, Messaging replaces the Chat and Collaboration features. Messaging was created to offer improved performance, group management, message notification, and configurable permissions to users. With Messaging, you can quickly send secure messages to any user on the InteleOrchestrator platform - administrators, radiologist assistants, technologists, and other radiologists.

With Messaging, you can:


- Find users by searching name, subspecialty, shift, location, and reading room
- Send messages to offline users
- Broadcast messages to multiple users
- View which users are online or offline
- Obtain a user's telephone number
- Send links of patient studies to provide or request a second opinion
- Set your status to available, busy, or offline
- Know where a user is currently situated

 To use this feature, ask your InteleOrchestrator administrator to enable the Messaging permission for your user account.

## Opening Messaging

The Messaging page can be opened from:

1. The Home screen.
2. The InteleOrchestrator Menu at the top of any InteleOrchestrator Application.
3. The Messaging icon in the Communication Taskbar.



When you have unread messages, the Messaging icon will be green and will show the number of unread messages.

 To use the Messaging feature, ensure that your InteleOrchestrator administrator has granted you the Access to User Messaging role within the User Role permissions. If you don't have access to the Messaging module, please reach out to your InteleOrchestrator administrator for assistance.

## Understanding the Messaging Page

The Messaging page provides a list of all Messaging users, including their current Status, Name, Specialties, Current Room Location, Department, Personal Phone Number, and Email Address.

This page allows you to initiate conversations, send messages to offline users, broadcast to multiple users, share study links with users, view telephone numbers, and copy email addresses.


You can also tell who is online by looking at the color of the messaging icon next to their name. In the below screenshot, all users are offline, as indicated by the gray messaging icon.

Group	User Type	Status	Name	Shift	Location	Reading Room	Reading Room Phone	Personal Phone	Email	
InteleOrchestrator Users	Radiologist Admin	Green	Administrator, Initial	Test Location	Room A					
		Gray	Hartmans, Jen	New regression2	New_Test	1231231234				
		Gray	Hebert, David	Unknown						
		Gray	Sampath, Kalaiavanan	New regression2	New_Test	1231231234				
		Gray	Shi, Dhi	Unknown						
		Gray	Two, Jen	New regression2	New_Test	1231231234				
		Gray	Sharma, Dhiraj H	New regression2	New_Test	1231231234				
		Gray	Nakade, Shyamal	New test test test	Test1	1231234123				
		Gray	test, test-Reg	Unknown						
		Gray	Tobey, Amy	test, test						
Portal Users	Hospital B	Gray	S. kalaiavanan	Unknown						
		Gray	uses, test	Unknown						
		Gray	Documentation, Admin	Unknown						
		Gray	Kawun, Paul	New regression2	New_Test	1231231234				
		Gray	Larsen, Sara	Unknown						
		Gray	zachtest, test	Unknown						
		Gray	Iyer, Revati	Unknown						
		Gray	tir, farshad	New test test test	Test1	1231234123				
		Gray	Kawun, paul2	New regression2	New_Test	1231231234				
		Gray	Sec, Dtest	Unknown						
Imaging Center 1	Gray	Moore, Bill	Unknown							
	Gray	hartmanns, jenitest	Unknown							
	Gray	coustells, kate	New regression2	New_Test	1231231234					
Imaging Center 2	Physician	Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
Imaging Center 3	Scheduler	Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								
		Gray								

## Opening the Messaging App Window

It's important to disable the pop-up blocker in your web browser to allow the InteleOrchestrator Messaging page to open successfully.

1. Click on the Messaging icon in the lower right corner of the Worklist or Messaging page.



- If you log in and you have unread messages, this Messaging icon will show green, you can click on it and it will take you to that conversation in the Messaging App window.

If you do not have any open messages, clicking this icon will open a new browser tab displaying the Messaging page.

The screenshot shows the InteleOrchestrator software interface. At the top, there is a navigation bar with links: Home, Worklist, Scheduling, Analytics, **Messaging**, Reporting, Projects, Management, Dashboard, CRM, Protocoling, Administrator, Initial, and Logout. Below the navigation bar, there are two main sections: "InteleOrchestrator Users" and "Portal Users".

**InteleOrchestrator Users:**

Role	Status	Name	Shift	Location	Reading Room	Reading Room Phone	Personal Phone	Email
Business Manager	0/2			Test Location	Room A			
Peer Review Admin	0/11			New regression2	New_Test	1231231234		
<b>Rad Assistant</b>	1/6	Sharma, Dhivaj H						
Radiologist	1/42	Tobey, Amy		New regression2	New_Test	1231231234		
Radiologist Admin	1/30	S. kalavanan		Unknown				
Resident	1/8	Iyer, Revati		Unknown				
Tech QA Reports	0/3	tir, farshad		New test test test	Test1	1231234123		
Transcriptionist	0/2							

**Portal Users:**

Hospital	Role	Status
Hospital B	Physician	0/15
	Scheduler	0/10
	Technologist	0/11
Imaging Center 1	Physician	0/16
	Scheduler	0/9
	Technologist	0/11
Imaging Center 2	Physician	0/16
	Scheduler	0/10
	Technologist	0/11
Imaging Center 3	Physician	0/7
	Scheduler	0/4

At the bottom right of the interface, it says "5:45 AM Mountain".

## Changing your Messaging Status

Your Messaging status determines your availability to communicate with colleagues via Messaging. Upon logging into InteleOrchestrator, your Messaging status icon is automatically set to the color green (Available). This means you are ready to send and receive messages and links related to ongoing orders.

The screenshot shows the InteleOrchestrator user dashboard. At the top, there are two tabs: "Users" and "Groups". Below the tabs, there is a section titled "InteleOrchestrator Users" which lists three users: "Business Manager", "Peer Review Admin", and "Rad Assistant". To the right of each user name is a green speech bubble icon with a white arrow pointing right, indicating they are "Available". Below the user names, their respective unread message counts are displayed: "0/12" for Business Manager and "0/0" for Peer Review Admin and Rad Assistant.


When you log out of the InteleOrchestrator your Messaging status icon automatically changes to the color gray (Offline).

The Messaging status icons are defined as:


Messaging status icon:	Indicates:
	Available: This status indicates that you are online and actively available to send and receive messages.
	Offline: When your status is set to "Offline," it means you are currently not logged into InteleOrchestrator. However, once you log in, you will receive any messages sent during your offline period.
	Busy: Choosing the "Busy" status signifies that you are online but occupied with other tasks. During this time, colleagues will not be able to send you messages.

## To change your Messaging status to Busy:

1. Click your Messaging status icon which is found in the header of the **InteleOrchestrator Users** list. Because you are online your Messaging status icon will show as green (Available).




2. A **Busy Time Limit** popup window is displayed. Set the desired time frame for the Busy status. The maximum Busy Time Limit is configured by the Administrator.



3. Click **Save**.
4. Your Messaging status icon will now be red with a line through it (Busy). Your colleagues will see a red Messaging status icon next to your name in their Messaging user list for the next 15 minutes.

## Setting Location Information

Location information is set when you open InteleOrchestrator. This pop up will appear at login if that permission is set by the administrator or you can click on the shift icon.



## Searching for Messaging Users

To include offline users in the list, check the "Show offline users" checkbox.

You can filter the list of InteleOrchestrator users by applying filters to display specific information and user groups.

To sort the columns in the Messaging list, simply click on a column heading. To reverse the sort order, click the column heading again. By default, the table is sorted by last name, first name in ascending order, while the Name column displays first name, last name.

# Sending Messages

## Sending a Message

1. Open Messaging. By default, all users online are shown.

The screenshot shows the InteleOrchestrator interface with the 'Messaging' tab selected. On the left, there are two sections: 'InteleOrchestrator Users' and 'Portal Users'. The 'InteleOrchestrator Users' section lists various roles and their counts: Business Manager (0/1), Peer Review Admin (1/1), Rad Assistant (1/6), Radiologist (2/42), Radiologist Admin (2/30), Resident (0/8), Tech QA Reports (0/3), and Transcriptionist (0/2). The 'Portal Users' section lists users from 'Hospital B' and three 'Imaging Center' locations. A search bar at the top right allows for searching by Name, Shift, Location, and Specialty. The bottom right corner shows the time as 1:13 PM Mountain.



The list only includes online InteleOrchestrator users. If you want to see offline users as well, make sure to select the **Show offline users** checkbox.

This screenshot shows the same messaging interface as above, but with a red box highlighting the 'Show offline users' checkbox in the search bar. The search bar also contains a placeholder text: '<Search Name, Shift, Location, Specialty>'. Below the search bar, there are columns for Location, Reading Room, Reading Room Phone, Personal Phone, and Email. A single user entry is visible: 'Headquarters Room A 206-555-4567'.

2. Double click on the user you want to send a message to.

InteleOrchestrator Users				
	Status	Name	Shift	Location
Business Manager	0/2	Two, Jen	New regression2	
Peer Review Admin	0/11	Sharma, Dhiraj H	New regression2	
Rad Assistant	1/6	Nakade, Shyamal	New test test test	
Radiologist	1/42	Yadav, Anjani	New test test test	
Radiologist Admin	1/30	Mart, Bobi	New regression2	
Resident	1/8	user, test	Unknown	
Tech QA Reports	0/3	McNamara, Jane	Unknown	
Transcriptionist	0/2	Iyer, Revati	Unknown	
		Automation, Testing	Test Location	
		Test, Login	Test Location	

3. This will open the Messaging App window.

In the screenshot, the browser title bar reads "InteleOrchestrator - Messaging App - Work - Microsoft Edge". The address bar shows the URL "https://devtest4-clario.intelerad.com/45/primary/swl/web/?talk#1724590562855". The main content area displays a list of messages. One message is visible, sent by "user, test" with the location "Unknown". At the bottom of the screen, there is a text input field with the placeholder "Type here".

4. In the dialog box, enter your text in the bottom field and press ENTER.

## Sending a Message to 2+ Users: Ad Hoc Groups

### Creating of Ad Hoc Groups

By using the ad hoc group functionality, you can efficiently engage with multiple users without creating a formal group for future messaging purposes.



When we close the Messaging App window, the Ad Hoc group will automatically go away.


1. On the left side of each user row, click the checkbox to select the desired users. You can select two or more users.

InteleOrchestrator Users			
	Status	Name	Shift
<b>Peer Review Admin</b>			
<input type="checkbox"/>		Two, Jen	
<input type="checkbox"/>		Sharma, Dhiraj H	
<input type="checkbox"/>		Nakade, Shyamal	
<input type="checkbox"/>		Yadav, Anjani	
<input type="checkbox"/>		Mart, Bobi	
<input checked="" type="checkbox"/>		user, test	
<input type="checkbox"/>		McNamara, Jane	
<input type="checkbox"/>		Iyer, Revati	
<input checked="" type="checkbox"/>		Automation, Testing	
<input checked="" type="checkbox"/>		Test, Login	
<input type="checkbox"/>		new, Create	

2. Once users are selected, the option '**Create ad hoc group**' will be active.
3. In the upper right corner of the window, click on '**Create ad hoc group**'.

InteleOrchestrator Users					
	Status	Name	Shift	Location	Reading Room
<b>Peer Review Admin</b>					
<input type="checkbox"/>		Two, Jen		New regression2	New_Test
<input type="checkbox"/>		Sharma, Dhiraj H		New regression2	New_Test
<input type="checkbox"/>		Nakade, Shyamal		New test test test	Test1
<input type="checkbox"/>		Yadav, Anjani		New test test test	Test1
<input type="checkbox"/>		Mart, Bobi		New regression2	New_Test
<input checked="" type="checkbox"/>		user, test		Unknown	
<input type="checkbox"/>		McNamara, Jane		Unknown	
<input type="checkbox"/>		Iyer, Revati		Unknown	
<input checked="" type="checkbox"/>		Automation, Testing		Test Location	Room A
<input checked="" type="checkbox"/>		Test, Login		Test Location	Room A
<input type="checkbox"/>		new, Create		Test Location	Room A

4. A Messaging App window will then appear, presenting an Ad Hoc group conversation that includes all the users selected in the previous step.




## Sending a Message to an Existing Group

1. On the Messaging page, click the Groups tab in order to display groups. For more information about creating private groups see [Creating a Private Group](#).

The screenshot shows a dark-themed user interface for managing groups. At the top, there are two tabs: 'Users' and 'Groups', with 'Groups' being the active tab and highlighted by a red box. Below the tabs is a search bar labeled 'Show All'. To the right of the search bar are two icons: a pencil for editing and a plus sign for creating new groups. A section titled 'Private groups (13)' follows, listing the names of the groups along with their status indicators. The groups listed are: 11 (1/1), AA (1/1), BB (1/5), CC (1/1), SS (1/1), Test (1/4), test private group (1/2), testgroup (1/2), testGroup (1/5), and TTTTTTTT (1/1). The 'testGroup' row is currently selected, indicated by a blue background.

Group <b>testGroup</b>	
	Status

2. Double click on the chosen group to initiate a conversation to open a Messaging App window.



3. In the dialog box, enter your text in the bottom field and press ENTER.

## Group Management in InteleOrchestrator Messaging

### Private Groups

Anyone with access to Messaging can build a Private Group. These will only be seen by members of the group and InteleOrchestrator administrators.

Public groups can be seen by anyone in Messaging. Permission to make them must be granted by the InteleOrchestrator administrator.

### Creating a Public or Private Group

As a user, you have the ability to create a Public or Private group that will remain in the lists of groups until removed. Certain User Role Permissions are required to create a Public group.

To create a group:

1. On the Messaging page, navigate to the Groups tab.

The screenshot shows the 'Groups' tab selected in the top navigation bar, indicated by a red box. Below the navigation bar, there is a 'Show All' checkbox and two icons: a pencil and a plus sign. A section titled 'Private groups (13)' is displayed, listing the following groups with their status:

Group Name	Status
11	1/1
AA	1/1
BB	1/5
CC	1/1
SS	1/1
Test	1/4
test private group	1/2
testgroup	1/2
testGroup	1/5
TTTTTTT	1/1

The 'testGroup' row is highlighted with a dark blue background.

2. Click on the plus icon to create a new group.

Users	Groups	
<input type="checkbox"/> Show All	<input type="button" value=""/>	<input type="button" value="+"/>
<input checked="" type="checkbox"/> Private groups (13)		
11	1/1	
AA	1/1	
BB	1/5	
CC	1/1	
SS	1/1	
Test	1/4	
test private group	1/2	
testgroup	1/2	
testGroup	1/5	

3. Provide a name for the new group.

Create new group

Group Name:	Example Group																																																																																														
Publish:	<input checked="" type="radio"/> Public <input type="radio"/> Private																																																																																														
InteleOrchestrator      Portal																																																																																															
Filter: <Search Name, Shift, Location, Specialty>																																																																																															
<table border="1"> <thead> <tr> <th>Name</th> <th>Shift</th> <th>Specialty</th> <th>Location</th> <th>Reading Room</th> <th>Reading Roo...</th> <th>Personal Phone</th> <th>Email</th> </tr> </thead> <tbody> <tr> <td>Kawun, paul2</td> <td></td> <td></td> <td>New regression2</td> <td>New_Test</td> <td>1231231234</td> <td></td> <td></td> </tr> <tr> <td>Kiser, Nora</td> <td></td> <td>Radiologists</td> <td>Test Location</td> <td>Room A</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Larsen, Sara</td> <td></td> <td>Radiologists</td> <td>Unknown</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Mart, Bobi</td> <td></td> <td>Bobi worklist group. CT Body, CT General, MSK, MSK CT Assignment Rule, Operations Group, Radiologist Reading Group A, Radiologists, Shr123, Shr12345, Shruti123, test 10 subspecialty, test 32 subspecialty</td> <td>New regression2</td> <td>New_Test</td> <td>1231231234</td> <td></td> <td></td> </tr> <tr> <td>Mart, Bobi_R...</td> <td></td> <td>Bobi worklist group. MSK, MSK CT Assignment Rule, Radiologist Reading Group A</td> <td>New regression2</td> <td>New_Test</td> <td>1231231234</td> <td></td> <td></td> </tr> <tr> <td>McNamara, J...</td> <td></td> <td>Radiologists</td> <td>Unknown</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Mel, Alex</td> <td></td> <td>Radiologists</td> <td>Unknown</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Mel, Alexandra</td> <td></td> <td>Radiologists</td> <td>Unknown</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Mirajkar, Amol</td> <td></td> <td>Cardiac</td> <td>Unknown</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td>Body CT- inactive20, radiotr1- inactive86</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>								Name	Shift	Specialty	Location	Reading Room	Reading Roo...	Personal Phone	Email	Kawun, paul2			New regression2	New_Test	1231231234			Kiser, Nora		Radiologists	Test Location	Room A				Larsen, Sara		Radiologists	Unknown					Mart, Bobi		Bobi worklist group. CT Body, CT General, MSK, MSK CT Assignment Rule, Operations Group, Radiologist Reading Group A, Radiologists, Shr123, Shr12345, Shruti123, test 10 subspecialty, test 32 subspecialty	New regression2	New_Test	1231231234			Mart, Bobi_R...		Bobi worklist group. MSK, MSK CT Assignment Rule, Radiologist Reading Group A	New regression2	New_Test	1231231234			McNamara, J...		Radiologists	Unknown					Mel, Alex		Radiologists	Unknown					Mel, Alexandra		Radiologists	Unknown					Mirajkar, Amol		Cardiac	Unknown							Body CT- inactive20, radiotr1- inactive86					
Name	Shift	Specialty	Location	Reading Room	Reading Roo...	Personal Phone	Email																																																																																								
Kawun, paul2			New regression2	New_Test	1231231234																																																																																										
Kiser, Nora		Radiologists	Test Location	Room A																																																																																											
Larsen, Sara		Radiologists	Unknown																																																																																												
Mart, Bobi		Bobi worklist group. CT Body, CT General, MSK, MSK CT Assignment Rule, Operations Group, Radiologist Reading Group A, Radiologists, Shr123, Shr12345, Shruti123, test 10 subspecialty, test 32 subspecialty	New regression2	New_Test	1231231234																																																																																										
Mart, Bobi_R...		Bobi worklist group. MSK, MSK CT Assignment Rule, Radiologist Reading Group A	New regression2	New_Test	1231231234																																																																																										
McNamara, J...		Radiologists	Unknown																																																																																												
Mel, Alex		Radiologists	Unknown																																																																																												
Mel, Alexandra		Radiologists	Unknown																																																																																												
Mirajkar, Amol		Cardiac	Unknown																																																																																												
		Body CT- inactive20, radiotr1- inactive86																																																																																													
<input type="button" value="Cancel"/> <input type="button" value="Submit"/>																																																																																															

4. If you have the permission, select either the Public or Private radio button. If you do not have the correct permission, the group will be Private.
5. Select the desired users.

Create new group

Group Name: Example Group

Publish:  Public  Private

InTeleOrchestrator Portal

Filter: <Search Name, Shift, Location, Specialty> Displaying 1 - 25 of 59 | Page 1 of 3

Name	Shift	Specialty	Location	Reading Room	Reading Roo...	Personal Phone	Email
Kawun, paul2			New regression2	New_Test	1231231234		
Kiser, Nora		Radiologists	Test Location	Room A			
Larsen, Sara		Radiologists	Unknown				
Mart, Bobi		Bobi worklist group, CT Body, CT General, MSK, MSK CT Assignment Rule, Operations Group, Radiologist Reading Group A, Radiologists, Shr123, Shr1245, Shruti123, test 10 subspecialty, test 32 subspecialty	New regression2	New_Test	1231231234		
<input checked="" type="checkbox"/> Mart, Bobi_R...		Bobi worklist group, MSK, MSK CT Assignment Rule, Radiologist Reading Group A	New regression2	New_Test	1231231234		
<input checked="" type="checkbox"/> McNamara, J...		Radiologists	Unknown				
<input checked="" type="checkbox"/> Mel, Alex		Radiologists	Unknown				
Mel, Alexandra		Radiologists	Unknown				
Mirajkar, Amol		Cardiac	Unknown				
		Body CT- inactive20, Cardiac1- inactive86					

Cancel Submit

6. Click on the submit button.

Create new group

Group Name: Example Group

Publish:  Public  Private

InTeleOrchestrator Portal

Filter: <Search Name, Shift, Location, Specialty> Displaying 1 - 25 of 59 | Page 1 of 3

Name	Shift	Specialty	Location	Reading Room	Reading Roo...	Personal Phone	Email
Kawun, paul2			New regression2	New_Test	1231231234		
Kiser, Nora		Radiologists	Test Location	Room A			
Larsen, Sara		Radiologists	Unknown				
Mart, Bobi		Bobi worklist group, CT Body, CT General, MSK, MSK CT Assignment Rule, Operations Group, Radiologist Reading Group A, Radiologists, Shr123, Shr1245, Shruti123, test 10 subspecialty, test 32 subspecialty	New regression2	New_Test	1231231234		
<input checked="" type="checkbox"/> Mart, Bobi_R...		Bobi worklist group, MSK, MSK CT Assignment Rule, Radiologist Reading Group A	New regression2	New_Test	1231231234		
<input checked="" type="checkbox"/> McNamara, J...		Radiologists	Unknown				
<input checked="" type="checkbox"/> Mel, Alex		Radiologists	Unknown				
Mel, Alexandra		Radiologists	Unknown				
Mirajkar, Amol		Cardiac	Unknown				
		Body CT- inactive20, Cardiac1- inactive86					

Cancel Submit

7. The group is created and is displayed in the private group list.

- Both Public or Private groups will remain in the lists of groups until removed.

## Deleting a Public or Private Group

Right click the group name in the Groups tab and click Delete.

If you do not see this option, you do not have the Permission set for your role. If you feel this is a mistake, speak to your InteleOrchestrator administrator.

## Leaving a Public or Private Group

Right click the group name in the Groups tab and click Leave.

If you do not see this option, you do not have the Permission set for your role. If you feel this is a mistake, speak to your InteleOrchestrator administrator.


## Sending an Exam Link within Messaging

With InteleOrchestrator Messaging, you can send a link to your current exam from the Messaging App. Your exam link can be shared in conversations both with individuals or among a group.

 When sharing the exam link with other individuals or a group, the individuals or group members must be credentialed in order to view the exam.

### To send an exam link within Messaging:


- Open the exam you want to send a link to.
- Click the Messaging icon. Select the people you want to send the exam link to. The Messaging App window will open.
- Click the **Send link of current exam** button.



Whether the exam link will be sent to the users or groups you've added to the message will depend on if those users and groups are credentialed to view the exam.


- **All Users are Credentialed**

If every user in the message is credentialed to access the exam, the message containing the exam link will be sent to all users.




- **No Users are Credentialed**

In cases where none of the users in the message are credentialed to access the exam, an error message will be displayed, which states: **The link to the current exam was not sent. These users are not credentialed for this exam.**




- **Users with Mixed Credentials**

When some of the users are not credentialed to view the exam, you can send the exam link to the group. However, a warning message will be displayed to the sender and all users within the message, stating: **One or more users in this conversation might not be credentialed.**



- Any users without the necessary credentials, attempting to open the exam link will receive an error message: **You are not credentialed for this exam.**



# 6

## DICTATION AND REPORTS

InteleOrchestrator Reporting offers a powerful way to use InteleOrchestrator virtually hands-free. Using only your microphone and the speech-to-text functionality, you can open and read exams, then dictate and sign your reports with minimal interaction with your keyboard or mouse. Whether you use this voice recognition functionality or dictate files that are sent to a transcriptionist, InteleOrchestrator makes it easy to dictate reports with report templates and the ability to combine reports or read reports consecutively.

### **In this chapter:**


Dictating Reports .....	165
Combining Reports and Reading Reports Consecutively .....	169
Changing a Dictation File After Submission for Transcription .....	172
InteleOrchestrator Speech Commands .....	176
Configuring Dictation Microphone Buttons to Perform Commands .....	179

# Dictating Reports


InteleOrchestrator provides users with two methods to dictate reports: dictating an audio file for a transcriptionist, or using InteleOrchestrator "virtually hands-free" navigation that allows you to use speech commands to navigate through InteleOrchestrator and your report, and dictate your findings directly into the fields of the report template you select.

## To dictate a report:


1. Open an exam in dictation mode. See "[Opening Exams](#)" (page 95)
2. In the Reports area, click the **Dictation** tab.



3. Under the Dictation tab, select either Voice Recognition (VR) or Transcription (TX) as your Dictation mode.




4. From the list of templates in the Template column, select a mapped template or a template structure to provide the framework for your report.



Whether you choose a mapped template or one of the template structures is a personal preference. The template structures can provide some or no structure, guiding you to dictate "free-form" if you wish, while the templates will have distinct fields to enter data into.


5. Do one of the following:

- Press the button on your microphone that you configured to start dictating the report in Voice Recognition mode. See ["Configuring Dictation Microphone Buttons to Perform Commands" \(page 179\)](#).
- In the Reports section, click **Start Recording** to start dictating the report in Transcription mode.



If you are using the Transcription mode, the background of the dictation area will change to a dark red background with a light red border to indicate that the microphone is active and audio is being recorded. The recording

 time indicator updates after you stop recording.




6. Dictate your findings into either the report fields or your transcription file.

 If you discover a critical finding that requires immediate action at any time during your dictation, you can click **Sign Prelim** button to create a preliminary report that can be acted on immediately.


7. Do one of the following to finish your report:

- For Voice Recognition Mode:
  - a. Press the button on your microphone that you configured to stop recording.
  - b. Press the button on your microphone that you configured to sign the report, use the voice command to sign the report, or click **Sign Report**.
- For Transcription Mode:
  - a. Click **Stop Recording** .

 The recording time indicator appears or updates. You can now review and add to the dictation, and change the playback speed and volume.



- b. When you are ready, click **Send to Transcription**.



There is a series of icons across the top of the Report Area you can use to stop your dictation, or interact with your reports in other ways. The following table explains the icons available and when to use each one.

Action	Icon	Reason for use
Save Report		This saves all of the data you have added to the report, and closes the exam. The exam stays assigned to you. You can open it to continue your reporting by clicking the <b>Dictation</b> module in the InteleOrchestrator menu and opening the report from your list of exams.
Cancel Report		This removes all of the changes you have made to the report since it was last saved and closes the exam. The exam stays assigned to you. You can open it to continue your reporting by clicking the <b>Dictation</b> module in the InteleOrchestrator menu and opening the report from your list of exams.
Discard Report		This discards all changes you have made to the report, closes the exam, and releases the exam to the queue for someone else to read.
Print Report		This prints the report.
Pop out Report Editor		This expands the Report Editor to its own window. This can work in conjunction with the Profile Management setting to keep the InteleOrchestrator Dictation Window Topmost. Setting this preference will keep the dictation window on top of any other applications you have open, and should prevent your system from losing focus while you are dictating your findings.

# Combining Reports and Reading Reports Consecutively

If a patient has multiple unread exams, you can configure InteleOrchestrator to prompt you to combine those exams or read them consecutively. While you can combine exams or choose to read them consecutively from the Patient View at any time, configuring this prompt can help you find exams that may not be visible on your worklist. Here are some examples why you might choose to enable these prompts:

- If your patient has an exam that could be combined with your current exam, InteleOrchestrator will prompt you when you open the exam to determine if you want to combine your current exam with one or more exams that share the same subspecialty or body part, for example. Additionally, your practice may have InteleOrchestrator configured to automatically combine exams based on specific criteria. This configuration will be done by the system administrator. For more information on the configuration, please contact your system administrator
- If your patient has multiple exams that could logically be read consecutively, InteleOrchestrator can ask you if you want another of these patient exams after you finish reading your current exam.

## To combine exams from the prompt:

1. If you have not done so already, verify that you have enabled combinable and/or consecutive reports in Profile Management. See [“Enabling Combinable or Consecutive Reports” \(page 71\)](#).
2. Open an exam. See [“Opening Exams” \(page 95\)](#)
3. In the Multiple Unread Exams Found dialog, click the exams you wish to combine or read consecutively.
4. Do one of the following:
  - Click **Read Consecutively** to open the selected exams one at a time after you finish reading the current exam.  
 To read exams consecutively, you must use the Auto-Next feature. See [“Using Auto-Next” \(page 99\)](#)
  - Click **Combine Exams** to combine the selected exams and open them together in an image viewer and reporting application.

**Multiple Unread Exams Found**

Additional exam(s) were found for this patient. Select the exam(s) you wish to open:


	Accession	Procedure	Modality	Date/Time	Priority	Assigned to	Status
<input checked="" type="checkbox"/>	CL7253	HEAD ANALYSIS	US	19/12/2022 07:...	Stat 1	-	Unread
<input type="checkbox"/>	CL7253	HEAD ANALYSIS	CT	19/12/2022 07:...	Routine	-	Unread
<input checked="" type="checkbox"/>	CL7253	HEAD ANALYSIS	US	19/12/2022 06:...	Stat 1	-	Unread

**Cancel** **Read Consecutively** **Combine Exams**

 You will not have the option to deselect the exams with the same modality once they are automatically combined.

## To combine exams from the Worklist:

1. Open an exam. See "[Opening Exams](#)" (page 95)
2. In the Worklist area, select the exams you want to combine.




The screenshot shows the RadAI Testing Worklist interface. A red box highlights the top row of the worklist table, which contains three exams. The first two exams have the 'Assignment' column dropdown open, showing options like 'Call PE', 'ASAP', 'Routine', and 'Medium'. The third exam's dropdown is closed. The worklist table has columns for Assignment, Priority, Status, Modality, Final Time, Exam Time, Site Procedure, and Accession. The worklist sidebar on the left shows various patient lists and a history section.

3. Click the drop-down arrow in the Actions column at the top of the Worklist area.
4. Select **Combine**.

The screenshot shows a software interface with a top navigation bar containing links for Home, Worklist, Scheduling, Analytics, Messaging, Reporting, Management, and Dashboard. Below this is a header with the text '(4.00)'. The main area is a table with columns for MRN, Actions, Assignment, Modality, and Site. There are eight rows of data, each representing an exam entry. A context menu is open over the fourth row (MRN MZ009966). The menu items are: Read Exam, Open in Viewer, Open in Dictation, Move to Reference, Hide Study, Add Note, Assign to Peer Review, DICOM Route, DICOM Route and Read, Change Status, Change Priority, Change Attending Radiologist, Cancel Exam, Delete Exam, and Combine. The 'Combine' option is highlighted with a red box.

5. In the Combine Exams dialog, click **Generate** to generate a combination ID number.




6. Click **Save**.

## To combine exams from the Patient View:

When the exam is already open in the Patient View and the user wants to add more exams to the report, the combined report icon is used.

1. Select the exams you want to combine in the Exams area of the Patient View. See "[Understanding the Patient View](#)" (page 25) for information about the Exams area.
2. Click the **Add for Combined Report** icon in the actions column to combine the selected exams.



3. If you haven't already, launch the viewer to see the exams for the combined report.
- Tip:** The exams will be combined only if they are opened in dictation mode.


## Changing a Dictation File After Submission for Transcription

If you need to make a change to a dictation after you submitted it for transcription, you open the report and append new audio to the transcription file. If your InteleOrchestrator administrators have allowed access to the Transcription worklist, you may also retrieve exams from the Exams Sent to Transcription worklist.


- Tip:** Once the transcriptionist opens an exam and begins transcribing, you can no longer edit the dictation.

## To change a dictation you have already submitted for transcription:

1. Search in your exam History or Exams Sent to Transcription worklist to find the exam you want to edit. See "[About the Worklist Column](#)" (page 21) for more details on the exam History section.



2. Open the exam in dictation mode. See "[Opening Exams](#)" (page 95).
3. Add dictation to the file.
4. Click **Send to Transcription** to resubmit the appended file to the transcriptionist.



## Template Structure Levels

Reporting offers three levels of report structure. The icons for these are located at the top of the report window just below the template header. The three options are Structured Template, Less Structured Template, and Narrative Mode.



These options can be utilized on any template once it has been loaded into the report editor. Select the appropriate structure level by either selecting it manually or using the respective voice command to select it.

- **Structured Template** - Utilizes the entire original template as it was created. It included fields, field options, auto fields, and default text.


The screenshot shows a dark-themed report editor interface with various sections and their default values:

- STUDY DESCRIPTION**: MR Abdomen w/o Contrast
- COMPLETED DATE/TIME**: 04/01/2020 6:20 PM
- REASON FOR STUDY**: reason
- COMPARISON**:
- EXAM PARAMETERS**:
  - TECHNIQUE:
  - LIMITATIONS: None
- FINDINGS**:
  - ESOPHAGUS: No structural or mechanical abnormality.
  - GALLBLADDER: No stones. Normal wall thickness. No pericholecystic fluid.
  - APPENDIX: Normal .
- IMPRESSION**:

- **Less Structured Template** - This option removes the default field option text from the report editor. However, the field options remain available when the respective field is selected

<b>STUDY DESCRIPTION</b>
MR Abdomen w/o Contrast
<b>COMPLETED DATE/TIME</b>
04/01/2020 6:20 PM
<b>REASON FOR STUDY</b>
reason
<b>COMPARISON</b>
<input type="checkbox"/>
<b>EXAM PARAMETERS</b>
TECHNIQUE: <input type="checkbox"/>
LIMITATIONS: <input type="checkbox"/>
<b>FINDINGS</b>
ESOPHAGUS: <input type="checkbox"/>
GALLBLADDER: <input checked="" type="checkbox"/>
APPENDIX: <input checked="" type="checkbox"/>
<b>IMPRESSION</b>
<input type="checkbox"/>

- **Narrative Template** - This option removes all field options and default text from the report editor. It only has auto fields and a single blank field within each section for navigation and to provide dictation.



## InteleOrchestrator Speech Commands

If your organization purchased InteleOrchestrator, you can use speech-to-text technology to navigate within InteleOrchestrator and dictate your reports. This allows you to work virtually hands-free for the entire time you are reading your exams. The following commands are configured by default.

- These speech commands can be configured by your InteleOrchestrator system administrator, so the commands you use may be different. Contact your InteleOrchestrator system administrator to get new speech commands.

### Dictating Commands

The Dictating commands are used when dictating a report. Simply speak the command run the associated action.

Action	Spoken Command
Select a template	"Load [name of template]"
Insert structured template	"Structured template"

Action	Spoken Command
Insert less structured template	"Less structured template"
Insert narrative template	"Narrative template"
Insert a macro	"Insert [name of macro]"
Select option from pick list	"Select [option]"
Begin a new line in the report	"New line"
Go to a different field	"Go to 'field name'" (ex. "Go to 'Reason for Study")
Save the report as a draft	"Save report"
Finalize the report	"Sign report"
Sign a prelim	"Sign prelim"
Cancel any edits you were making to a draft or prelim	"Cancel report"
Delete or discard your report	"Discard report"
Delete your last statement	"Delete that"
Go to previous field	"Go to previous"
Go to next field	"Go to next"
Agree to default text	"Agree"
Agree to default text and go to next field	"Agree and next"
Cancel the dialog	"Cancel"
Clear the dialog	"Clear"
Highlight your last statement	"Select that"
Highlight the entire report	"Select all"
Make text bold	"Format bold"
Make text italicized	"Format italics"
Make text underlined	"Format underline"
Add a superscript	"Format superscript"
Add a subscript	"Format subscript"
Strikethrough	"Format strikethrough"
Align left	"Format left"
Align right	"Format right"

Action	Spoken Command
Align center	"Format center"
Fit text to window	"Format justify"
Create bullet points	"Format bullets"
Create numbers	"Format numbering"
Indent a line	"Format indent"
Outdent a line	"Format outdent"
Add a template remark	"Add template remark"
Insert dates for today, yesterday, tomorrow, now	"Date [date option]"
Save	Save
Save report	"Save report"
Submit	"Submit"

## Worklist Commands

There are commands to navigate worklists. Simply speak the command run the associated action.

Action	Spoken Command
Launch Auto-Next	"Start reading"
Pause Auto-Next	"Pause reading"
Resume Auto-Next	"Resume reading"
Stop/Exit Auto-Next	"Stop reading"
Read first exam on the worklist	"Read first exam"
Refresh worklist	"Refresh my worklist"
Search accession	"Search accession"
Search first name	"Search first name"
Search last name	"Search last name"
Search history	"Search history"
Search MRN	"Search M R N"
Search site	"Search site"
Search site procedure	"Search site procedure"

Action	Spoken Command
Search submit	"Search submit"
Search worklist	"Search worklist"

## Notes Commands

There are commands to create the most common notes in exams. Simply speak the command to run the associated action.

Action	Spoken Command
Add a communication note	"Create comm note"
Add an exam note	"Create exam note"
Add a patient note	"Create patient note"

# Configuring Dictation Microphone Buttons to Perform Commands

InteleOrchestrator supports the use of several dictation microphones and USB headsets. See the *Intelerad Workstation and Accessories Guide* on the Intelerad Service Portal for a full list of supported microphones.


 The following configuration options apply to Phillips, Olympus, and Nuance dictation microphones only.

You can configure the handset buttons on your microphone to perform different commands, based on your preference. This feature enables you to read your exams, navigate between auto-text fields, and dictate and sign your reports with minimal interaction with your keyboard or mouse.

### To configure handset buttons:

1. Click your user name in the top right corner of the interface, and then select **Profile Management**.
2. In the **Reporting Preferences** section, click **Configure**.

3. In the **Microphone Button Configuration** dialog, for each button, select the associated command.



4. Click **Save**.

# 7

## PROTOCOL WORKFLOW

InteleOrchestrator provides the protocol requests as a worklist of orders requiring protocoring directly in your worklist. These exams are identified from order status from the HL7 message or by special handling rule setup in InteleOrchestrator.

### In this chapter:

Opening Exams .....	182
Performing an Advanced Search .....	182
Understanding the Protocol Workflow .....	183

# Opening Exams

From the worklist of orders requiring prototyping, you can open an exam. The exam will open in prototyping mode.

## To open an exam:

1. Access the worklist containing protocol requests.

The screenshot shows the 'Prototyping Needed' worklist. On the left, there is a sidebar with 'My Worklist' sections: 'My Reading Queue' (137 items), 'My Support Queue' (0 unread exams), 'EPIC' (Time test 1), 'Imaging Center 1', and 'Prototyping Needed' (1 item). The main area displays a grid of patient records with columns for Ordered Time, Referring Physician, Accession, Actions, Assignment, Patient Name, Begun Time, Final Time, Prelim Time, Signed Time, and Status. A red box highlights the 'Actions' column for the first row, which has a checkmark icon.

2. In the Action Column for the required exam, click to open the exam for prototyping only.

# Performing an Advanced Search

If the desired exam is not in your worklist, you can perform an advanced search to find an exam outside your worklist or with advanced search criteria. The results of your search appear in the work area.

## To perform an advanced search:

1. Click the Advanced Search icon at the top of the InteleOrchestrator interface.

The screenshot shows the Advanced Search dialog. It includes fields for 'Priority', 'Status', 'Modality', 'Final Time', 'Exam Time', 'Site Procedure', and 'Accession'. A red box highlights the 'Status' dropdown menu, which is set to 'Protocol Needed'. Below the dialog, the worklist shows several exams listed with their status as 'Protocol'.

2. In the Advanced Search dialog, enter the data you want to search for from the exams.

For example, if you want to find all protocols that need exams, select Protocol Needed from the drop-down menu in the **Status** field.

The screenshot shows the InteleOrchestrator search interface with the following details:

- Search Results (0)**: The search term "Protocol Needed" is entered in the search bar.
- Advanced Search**: The "Exam" filter is selected.
- Keywords**: "Protocol Needed" is listed under the "Status" dropdown.
- Count**: The count is set to "Unlimited".
- Search**: The "Search" button is highlighted with a red box.
- Results**: A table displays search results with columns: First Name, Last Name, Status, Modality, Final Time, Exam Time, Site Procedure, and Accession.
- Filters**: Various filters are applied, including "Protocol Needed" in the Status field and "Protocol Needed" in the "Assigned To" dropdown.
- Time Range**: The time range is set to "Today" and "3d".
- Buttons**: "Add to Worklist", "Print", "Email", and "Copy" buttons are visible.

InteleOrchestrator automatically limits the date range to within 3 days unless you specify a date range.

3. Click **Count** to see how many results your search criteria returned. [Optional]
4. Click **Search**.
5. Select the desired exam in the list of exams returned from your search.

## Understanding the Protocol Workflow

This section provides the different steps to be followed to protocol an exam.


### To protocol an exam:

1. If an exam is not opened then open an exam in prototyping mode using the steps as provided in the "[Opening Exams](#)" section.


The screenshot shows the "Protocol Needed Worklist" interface with the following details:

- Worklist**: The "Protocol Needed" folder is selected.
- Protos**: A progress bar indicates 100% completion.
- Actions**: A column of icons for each exam entry.
- Columns**: The table has columns for Patient Name, Begun Time, Final Time, Prelim Time, and Signed Time.
- Data**: A list of patients with their corresponding physician names and accession numbers.
- Filters**: "Protocol Needed" is selected in the "Assignment" dropdown.
- Buttons**: "Print", "Email", and "Copy" buttons are visible.

- The Protocols tab in patient view will be displayed. From the list of templates in the Template column, select a template to provide the framework for your Protocol. Templates will be shown based on system administrator configuration.



- Click in a field to edit the text or you can use Voice Recognition mode. Click **Start Recording** to start dictating. For more information, see [Configuring Dictation Microphone Buttons to Perform Commands](#).




4. Select an icon to interact with your protocol in other ways. The following table explains the icons available and when to use each one.

Action	Icon	Description
Save Report		This saves all of the data you have added to the protocol, and closes the exam. The exam stays assigned to you. You can open it to continue your protociling by opening the Protocol from your list of exams.
Cancel Report		This removes all of the changes you have made to the protocol since it was last saved and closes the exam. The exam stays assigned to you. You can open it to continue your protociling by opening the Protocol from your list of exams.
Discard Report		This discards all changes you have made to the protocol, closes the exam, and releases the exam to the queue for someone else to read.
Print Report		This prints the protocol.
Pop out Report Editor		This expands the Protocol Editor to its own window.

5. When you click the save icon for an edited protocol, the protocol will be saved as a draft. Similarly, if you start to edit a protocol and click out to a worklist, the protocol will be saved as a draft.

If you are not sure which worklist contains your draft protocol, contact your system administrator.

6. When you are finished with the protocol, Click **Submit Protocol**



7. The protocol will be sent to your RIS or EMR based on rules defined by your system administrator.

# 8

## DOCUMENT VIEWER

InteleOrchestrator offers a powerful way to view documents directly within the InteleOrchestrator Patient View. This feature allows InteleOrchestrator users to view documents easily when they may be reading from disparate systems where documents are not reliably stored in the same place.

### **In this chapter:**

Verify your Document Viewer Configuration Settings .....	188
Document Viewer Preferences .....	188
Viewing Scanned Documents .....	190

# Verify your Document Viewer Configuration Settings

To ensure that the following document viewer settings are enabled:

1. Open a web browser session, and go to the InteleOrchestrator URL.
2. Log in to InteleOrchestrator.
3. From the InteleOrchestrator menu at the top of the interface or from the list of icons, click **Management**.
4. Click **Configuration**.
5. Set the Document Viewer preferences to customize the behaviour of the Document Viewer.

The screenshot shows the InteleOrchestrator configuration interface. The left sidebar has a 'Tables' section with various options like Active Directory, Configuration, Document Viewer (with Vendor selected), Exam Workflow, Notes, Organization, Portal, Productivity Gauge, Reporting, Security, Services, and Worklist. The main panel is titled 'Configuration of Vendor' and shows 'Document Vendor: Kailo'. Under 'General Configuration', there are fields for 'Source' (http://52.243.70.225:5145/), 'Default Document' (Requisition), 'Document Display' (Display Automatically checked), and 'Default Display' (Split View selected). Below this is a 'Documents Category Config' section with a table:

Category	Name	Icon	Role (permission)	Inverted
Requisition	Requisition	file		
Report	Report	file		
Clinical	Clinical	file		
Patient	Patient	file		

## Document Viewer Preferences

Set the Document Viewer preferences to customize the behaviour of the Document Viewer.

The screenshot shows the 'Configuration of Vendor' screen in the YOURRad Practice application. The left sidebar lists various configuration categories like Tables, Active Directory, Configuration, Document Viewer, Vendor, Exam Workflow, Notes, Organization, Portal, Productivity Gauge, Reporting, Security, Services, and Worklist. The 'Vendor' category is selected. The main panel has two sections: 'General Configuration' and 'Documents Category Config'. In 'General Configuration', the 'Document Vendor' is set to 'Kailo'. Under 'Source', the URL is 'http://52.243.70.225:5145/'. The 'Default Document' is set to 'Requisition'. Under 'Document Display', 'Display Automatically' is checked. Under 'Default Display', 'Split View' is checked. In the 'Documents Category Config' section, there is a table titled 'Adjust Category settings for scanned document management' with columns for Category, Name, Icon, Role (permission), and Inverted. The categories listed are Requisition, Report, Clinical, and Patient, each with their respective names and icons.

Category	Name	Icon	Role (permission)	Inverted
Requisition	Requisition	document icon		
Report	Report	document icon		
Clinical	Clinical	document icon		
Patient	Patient	document icon		


## General Configuration

This section is about the configurations specific to vendor.

**Default Document** – The document type that is displayed upon entering the Document tab in Patient View.

**Document Display** – When you enable this option and you open an order, the Document Viewer automatically opens and displays the scanned documents that are associated with that order.

**Default Display** - When you enable this option and you open an order, the display will be in split view and the document specified in drop down will be displayed in the second window.



## Documents Category Configuration

This section is about the configurations of each individual document category applies the preferences described below only when you open new orders. These preferences do not affect orders that you have already opened.

**Category** – Type of documents which are supported.

**Name** – Ability to rename document type in application.

**Icon** – Upload different icon for the document type.

**Role** – Users associated to specific roles are configured in order to have access to that document type.


**Inverted** – Reverses the colors of a scanned document. For example, black text appears white, and a white background appears black.

## Viewing Scanned Documents

In the Document Viewer, you can view scanned documents that are associated with an order, such as requisitions and reports stored in RIS. You can also view the scanned documents for all priors associated with the current order.

To use this feature, ask your administrator to enable the Documents Display privilege in your user account.

You can set the Document Viewer to display only documents of specific types, such as REQUISITION and REPORT. If you would like to limit the type of documents that appear in the Document Viewer, have your administrator contact Intelerad Technical Support to filter out the unwanted document types.



In the Document Viewer, you can view Adobe Portable Document Format (PDF) files and scanned documents saved as the following image files:

- PNG
- JPEG
- BMP
- GIF
- TIFF (including multipage TIFF)

For unsupported document types it will show up as a download link.

You can use the following tools to manipulate scanned document:


Click :	To:
	Pan around the document when it is larger (zoomed in) than the viewport available.
	Rotate the current page left or right in 90-degree increments.
	Inverts the text and background colours for the document.
	Switch between the single and split screen views for the scanned docs.

Click :	To:
	Choose the magnification of the documents: Fit Page (default), Fit Width, Fit Height, Actual Size, 50%, 150%, or 200%.
	Auto fit the document to the width of the screen.
	Auto fit so the entire document is visible in the screen.
	Revert the documents to their original magnification, orientation, and colour settings.
	Retain the function of uploading scanned docs to an exam, however we'll move that function from the 'file uploads' note to the Scanned Docs tab toolbar. The actual 'note' or upload window that pops up should remain the same.
	Pop the scanned document viewer into a separate browser window (existing functionality).
	Switch between documents in the scanned document viewer.
	Pager if document contains more than one page.
	Information tooltip to provide more information on document.

### To view scanned documents:

1. From the worklist or the Search tool, open the order that contains the scanned documents that you want to view. If you have your Document Viewer preferences set to open the Document Viewer automatically, the Document Viewer appears when you open the order.
2. To view the scanned document, click a document thumbnail.
3. To navigate between pages in a multi-page document, use the navigation arrows underneath the scanned document. If the Zoom mode is set to Fit Page, you can use the mouse wheel to scroll between documents or between pages in a multi-page document.
4. To manipulate the scanned document, use the tools in the Document Viewer.
5. To view scanned documents into a separate browser window, click the Pop-out Document icon in the top right corner of the document viewer window. The scanned

documents which open in a separate window helps the users logout sessions.



6. To close the Document Viewer window, click the “x” button in the top right corner of the window. To reopen the Document Viewer, click Show Scanned Documents in the Patient Record tab.

# 9

## PEER REVIEW WORKFLOW

The Peer Review Workflow is an essential component of a comprehensive radiology department quality assurance program. The Peer Review Workflow provides a useful mechanism for evaluating the diagnostic accuracy of the radiologist's reports. InteleOrchestrator can be configured for Peer Review workflows, providing greater efficiency and accuracy in the interpretation of reports.

### In this chapter:

Peer Review Workflow .....	195
Opening Peer Reviews .....	195
Peer Reviews to Evaluate Worklist .....	202
Peer Reviews to Arbitrate Worklist .....	205
Peer Review Scorecard .....	212
Peer Review Analysis Report .....	213

# Peer Review Workflow

The Peer Review Workflow is an essential component of a comprehensive radiology department quality assurance program. The Peer Review Workflow provides a useful mechanism for evaluating the diagnostic accuracy of the radiologist's reports. InteleOrchestrator can be configured for Peer Review workflows, providing greater efficiency and accuracy in the interpretation of reports.

## About the Peer Review Worklist

The Peer Review Worklist displays cases assigned to the radiologist for Peer Review. The Peer Reviews are assigned based on the rules setup by the Peer Review Administrators.

These are the peer review worklists that the system will show by default if Peer Review is configured. These worklists can be customized by your system administrator:

- Peer Reviews Assigned to Me – Peer reviews in a pending status that are assigned to the current user.
- Peer Reviews to Evaluate – Reports that have been peer reviewed by another colleague that the current user needs to review and potentially take action.
- Peer Reviews to Arbitrate – Reports that have been peer reviewed that require the Peer Review Administrator to determine the outcome of a discrepancy between peer review ratings.

Actions	Site	Location	Accession	Site Procedure	MRN	Patient Name	Radiologist
[Actions]	Hospital A	UNKNOWN	AG1440653...	NULL	CL-192501E	Tavares, ANTONIO	
[Actions]	Imaging Center 1	UNKNOWN	CL-203261-IR...	IR REPLACEMENT NON TUNNEL...	CL-192437R	Derdzinski, BENITO	
[Actions]	Imaging Center 1	UNKNOWN	CL-203231-...	CT TAC W/O CM	CL-192356H	Cherry, FLORENTINO	Nakade, Shyan
[Actions]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaivanan	Administrator, [REDACTED]
[Actions]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaivanan	Administrator, [REDACTED]
[Actions]	Hospital A	UNKNOWN	00021	Cervical Spine X-Ray	00002	Dhi, Raj	Hartmann, Jen

## Opening Peer Reviews

To perform a peer review of an exam:

1. Open InteleOrchestrator.
2. If required, set your shift hours and click **Save**.
3. In the **Worklist** column select **Peer Review**. The Peer Review worklist will expand to show multiple choices.

The screenshot shows the InteleOrchestrator interface with the 'Worklist' column expanded to show the 'Peer Review' section. The 'Peer Review' section contains three items: 'Peer Reviews Assigned to Me' (6), 'Peer Reviews to Arbitrate' (1), and 'Peer Review to Evaluate' (1). A red box highlights this section. The main work area displays a grid of peer review assignments with columns for Actions, Site, Location, Accession, Site Procedure, MRN, Patient Name, Radiologist, Ordering, and Reason. The grid shows six rows of data.

Actions	Site	Location	Accession	Site Procedure	MRN	Patient Name	Radiologist	Ordering	Reason
[checkbox]	Hospital A	UNKNOWN	AG1440653...	NULL	CL-192501E	Tavares, ANTONIO	James Tournas	Manual	
[checkbox]	Imaging Center 1	UNKNOWN	CL-203261-IR	IR REPLACEMENT NON TUNNEL...	CL-192437R	Derdzinski, BENITO	Trish Smith	Manual	
[checkbox]	Imaging Center 1	UNKNOWN	CL-203231-...	CT IAC W/O CM	CL-192356H	Cherry, FLORENTINO	Nakade, Shyamal	Trish Smith	Manual
[checkbox]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaisaran	Administrator, Initial	BOB BECHARA	Manual
[checkbox]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaisaran	Administrator, Initial	BOB BECHARA	Manual
[checkbox]	Hospital A	UNKNOWN	00021	Cervical Spine X-Ray	00002	Dhi, Raj	Hartmann, Jen	BOB BECHARA	Manual

4. Click **Peer Reviews Assigned to Me** from the Worklist column. The Peer Review notes are displayed in the work area.

The screenshot shows the InteleOrchestrator interface with the 'Worklist' column expanded to show the 'Peer Review' section. A red arrow points to the 'Peer Reviews Assigned to Me' item, which has a value of 6. The main work area displays a grid of peer review assignments with columns for Actions, Site, Location, Accession, Site Procedure, MRN, Patient Name, Radiologist, Ordering, and Reason. The grid shows six rows of data.

Actions	Site	Location	Accession	Site Procedure	MRN	Patient Name	Radiologist	Ordering	Reason
[checkbox]	Hospital A	UNKNOWN	AG1440653...	NULL	CL-192501E	Tavares, ANTONIO	James Tournas	Manual	
[checkbox]	Imaging Center 1	UNKNOWN	CL-203261-IR	IR REPLACEMENT NON TUNNEL...	CL-192437R	Derdzinski, BENITO	Trish Smith	Manual	
[checkbox]	Imaging Center 1	UNKNOWN	CL-203231-...	CT IAC W/O CM	CL-192356H	Cherry, FLORENTINO	Nakade, Shyamal	Trish Smith	Manual
[checkbox]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaisaran	Administrator, Initial	BOB BECHARA	Manual
[checkbox]	Hospital A	UNKNOWN	00018	Cervical Spine X-Ray	1111222334	s, kalaisaran	Administrator, Initial	BOB BECHARA	Manual
[checkbox]	Hospital A	UNKNOWN	00021	Cervical Spine X-Ray	00002	Dhi, Raj	Hartmann, Jen	BOB BECHARA	Manual

5. Click a Peer Review note in the work area.

The screenshot shows the InteleOrchestrator Worklist interface. On the left, there's a sidebar with navigation links like Home, Worklist, Scheduling, Analytics, Messaging, Reporting, Management, Dashboard, CRM, Protocoling, and Dictation. The main area displays a table titled "Peer Reviews Assigned to Me Worklist (6)". The columns are Actions, Site, Location, Accession, Site Procedure, MRN, Patient Name, and Radiologist. One row is highlighted with a red border, representing a cervical spine X-ray from Hospital A.

6. The Peer Review note is displayed. Here you can view the order, read the report and launch images in the viewer. Click located in the top right corner to open images in the viewer.

The screenshot shows a detailed view of a peer review note for a CXR study. At the top, it says "Smith, Ramon" and "DX33717 - CHEST SINGLE VIEW". Below this is a table with patient and exam details. Further down, there's a section for "EXAM DESCRIPTION", "REASON FOR STUDY", "COMPARISON", "NUMBER OF VIEWS", "TECHNIQUE", and "Note". At the bottom, there are fields for "Reason", "Rating", "Peer Learning", and a message input field. The "Launch Viewer" button in the top right corner is highlighted with a red box.

7. You can choose to either dismiss a peer review or perform a peer review. Click **Dismiss Peer Review** if you choose to dismiss the peer review. The peer review will be removed from the assigned worklist.
8. To perform a Peer Review, select a peer review **Rating** from the drop down menu.

Peer Review

Smith, Ramon  
DX33717 - CHEST SINGLE VIEW

Site	Community Hospital	DOB	03/29/2019	Ordered Time	03/09/2021 11:03 PM
Priority	Stat	Gender	Male	Begun Time	03/18/2021 12:59 AM
Exam Reason	TRAUMA	Pre-Authorization		Unread Time	03/07/2021 7:49 PM
Location	Ortho	Age	3y	Acquisition Time	
Ordering	T. William Fisher	MPI		Prelim Time	12/11/2020 7:59 PM
Radiologist	N/A	MRN	6900	Final Time	12/29/2020 12:49 PM

Report templates courtesy of EIDETICO. Visit eidetico.com for license details.

EXAM DESCRIPTION: CHEST SINGLE VIEW

REASON FOR STUDY: rm tr2 mvc

COMPARISON: None.

NUMBER OF VIEWS: One view.

TECHNIQUE: single frontal radio-graphic view of the chest acquired.

Note

Reason:  Manual

Rating:

Peer Learning:  Test

<Message>

1 - Agreement  
2a - Subtle abnormality (not likely to be clinically significant)  
2b - Subtle abnormality (likely to be clinically significant)

Cancel Clear Auto-next Dismiss Peer Review Submit Peer Review

9. If required, select a **Peer Learning** option from the drop down menu.

Peer Review

Always Launch Viewer

**Smith, Ramon**  
**DX33717 - CHEST SINGLE VIEW**

Site	Community Hospital	DOB	03/29/2019	Ordered Time	03/09/2021 11:03 PM
Priority	Stat	Gender	Male	Begun Time	03/18/2021 12:59 AM
Exam Reason	TRAUMA	Pre-Authorization		Unread Time	03/07/2021 7:49 PM
Location	Ortho	Age	3y	Acquisition Time	
Ordering	T. William Fisher	MPI		Prelim Time	12/11/2020 7:59 PM
Radiologist	N/A	MRN	6900	Final Time	12/29/2020 12:49 PM

Report templates courtesy of Eidecico. Visit [eidecico.com](http://eidecico.com) for license details.

EXAM DESCRIPTION: CHEST SINGLE VIEW

REASON FOR STUDY: rm tr2 mvc

COMPARISON: None.

NUMBER OF VIEWS: One view.

TECHNIQUE: Single frontal radiographic view of the chest acquired.

Note

Reason:<sup>\*</sup> Manual

Rating:<sup>\*</sup>

Peer Learning:<sup>\*</sup>

<Message>

- Feedback 1.0
- Feedback 2
- Peer Learning Worklist popout Test
- Peer Learning Worklist Test

Cancel Clear Auto-next Dismiss Peer Review Submit Peer Review

10. Enter a message (Optional).

Peer Review

Always Launch Viewer

**Smith, Ramon**  
**DX33717 - CHEST SINGLE VIEW**

Site	Community Hospital	DOB	03/29/2019	Ordered Time	03/09/2021 11:03 PM
Priority	Stat	Gender	Male	Begun Time	03/18/2021 12:59 AM
Exam Reason	TRAUMA	Pre-Authorization		Unread Time	03/07/2021 7:49 PM
Location	Ortho	Age	3y	Acquisition Time	
Ordering	T. William Fisher	MPI		Prelim Time	12/11/2020 7:59 PM
Radiologist	N/A	MRN	6900	Final Time	12/29/2020 12:49 PM

Report templates courtesy of Eidecico. Visit [eidecico.com](http://eidecico.com) for license details.

**EXAM DESCRIPTION:** CHEST SINGLE VIEW

**REASON FOR STUDY:** rm tr2 mvc

**COMPARISON:** None.

**NUMBER OF VIEWS:** One view.

**TECHNIQUE:** Single frontal radiographic view of the chest acquired.

**Note**

Reason: \* Manual

Rating: \* 1 - Agreement

Peer Learning: \* Peer Learning Worklist popout Test

<Message>

11. You also have the option to first click **Auto-next** and then click **Submit Peer Review** to automatically go to the next exam to be peer reviewed.

Peer Review

Always Launch Viewer

**Smith, Ramon**  
**DX33717 - CHEST SINGLE VIEW**

Site	Community Hospital	DOB	03/29/2019	Ordered Time	03/09/2021 11:03 PM
Priority	Stat	Gender	Male	Begun Time	03/18/2021 12:59 AM
Exam Reason	TRAUMA	Pre-Authorization		Unread Time	03/07/2021 7:49 PM
Location	Ortho	Age	3y	Acquisition Time	
Ordering	T. William Fisher	MPI		Prelim Time	12/11/2020 7:59 PM
Radiologist	N/A	MRN	6900	Final Time	12/29/2020 12:49 PM

Report templates courtesy of Eeidetico. Visit eidetico.com for license details.

EXAM DESCRIPTION: CHEST SINGLE VIEW

REASON FOR STUDY: rm tr2 mvc

COMPARISON: None.

NUMBER OF VIEWS: One view.

TECHNIQUE: Single frontal radiographic view of the chest acquired.

Note

Reason:  Manual

Rating:  1 - Agreement

Peer Learning:  Peer Learning Worklist popout Test

<Message>

Cancel Clear  Auto-next Dismiss Peer Review Submit Peer Review

12. Click **Submit Peer Review** to save your changes.

Peer Review

Always Launch Viewer

**Smith, Ramon**  
**DX33717 - CHEST SINGLE VIEW**

Site	Community Hospital	DOB	03/29/2019	Ordered Time	03/09/2021 11:03 PM
Priority	Stat	Gender	Male	Begun Time	03/18/2021 12:59 AM
Exam Reason	TRAUMA	Pre-Authorization		Unread Time	03/07/2021 7:49 PM
Location	Ortho	Age	3y	Acquisition Time	
Ordering	T. William Fisher	MPI		Prelim Time	12/11/2020 7:59 PM
Radiologist	N/A	MRN	6900	Final Time	12/29/2020 12:49 PM

Report templates courtesy of Eidecico. Visit [eidecico.com](http://eidecico.com) for license details.

**EXAM DESCRIPTION:** CHEST SINGLE VIEW

**REASON FOR STUDY:** rm tr2 mvc

**COMPARISON:** None.

**NUMBER OF VIEWS:** One view.

**TECHNIQUE:** Single frontal radiographic view of the chest acquired.

**Note**

Reason:<sup>\*</sup> Manual

Rating:<sup>\*</sup> 1 - Agreement

Peer Learning:<sup>\*</sup> Peer Learning Worklist popout Test

<Message>

Cancel Clear Auto-next Dismiss Peer Review **Submit Peer Review**

## Peer Reviews to Evaluate Worklist

The Peer Reviews to Evaluate Worklist contains reports that have been peer reviewed by another colleague that the current user needs to review and potentially take action.

**To evaluate reports that received positive feedback or an unfavorable peer review rating:**

1. Click **Peer Reviews to Evaluate** from the worklist column. The details of the **Peer Reviews to Evaluate** are displayed in the work area.

InteleOrchestrator™

Home | **Worklist** | Scheduling | Analytics | Messaging | Reporting | Management | Dashboard | CRM | Protocoling | Dictation | <Access

**Worklist** | Folder

Peer Review

**RVU** | **1** Hospital A Worklist


Current: 0.00 % Suggested: 24.49%

- My unread exams: 8
- Unread: 338
- CombineWL: 0
- Imaging Center 1: 102
- Test WL: 429
- My Support Queue**
- \* Resident Que ICL
- \* Resident Que HA
- Radiologist Worklists
- Operations
- Peer Review
  - Peer Reviews Assigned to Me: 6
  - Peer Review to Evaluate: 1**
  - Peer Reviews to Arbitrate: 1


Actions Site Location Accession Site Procedure MRN Patient Name Radiologist

... Imaging Center 1 UNKNOWN CL-203187... CT UPPER EXT RIGHT W CL-192014A Bartolero, HARLAN Larsen, Sara

- Click the Peer Review Note in the work area.



- This will direct you to the Patient View where you can see the report and view the received Peer Review Notes.



- Click the Peer Review Note in the Evaluation status of the Notes panel.

5. A new window opens, where you can agree or disagree with the Peer Review rating.

**To Agree:** Click the **Agree** button. The Peer Review will move to final status.

Peer Review

Rating: 3a - Finding should usually be made (not likely to be clinically significant)

Classification: EQUIP 3

Reason: Requested

Category: Vascular

Peer Learning: Peer Learning Worklist popout Test

Created: Larsen, Sara, 03/20/2021 9:49 AM

Evaluation

Classification: EQUIP 3


Rating: 4

Comment (optional)

**Close** **Disagree** **Agree** **Agree and Addend** **Comment**

**To Disagree:** Select a new score from the **Rating** drop-down menu. You have the option to enter a message or reason for your disagreement.

Click **Disagree**, if configured, the peer review will move to an arbitration status for the Peer Review Administrator to review and assign the final rating.



## Peer Reviews to Arbitrate Worklist

The Peer Reviews to Arbitrate Worklist contains peer reviewed reports that require the Peer Review Administrator to determine the outcome of a discrepancy between peer review ratings. Only the Peer Review Administrator will see this worklist.

**To arbitrate exams that received a disagreement for peer review rating during the evaluation:**

1. Click Peer Reviews to Arbitrate from the Peer Review worklist in the Worklist column. The list of Peer Review Notes are displayed in the work area.

The screenshot shows the InteleOrchestrator software interface. On the left, there is a navigation sidebar with sections like 'Worklist', 'Peer Review', 'Operations', and 'Peer Review'. Under 'Peer Review', the 'Peer Reviews to Arbitrate' item is highlighted with a red box. The main work area displays a table with columns: Actions, Site, Location, Accession, Site Procedure, MRN, Patient Name, Radiologist, Ordering, and Reason. One row is visible, showing Hospital A, UNKNOWN, 14009, IR PLACEMENT NEPHRO URETE..., CL-192116J, Orlow, DONNELL, Hartmann, Jen, Alexandre Messier. At the bottom, there is a search bar and a footer with various links and a timestamp of 2:25 PM Mountain.

2. Click on a Peer Review Note displayed in the work area.

This screenshot is similar to the previous one, showing the InteleOrchestrator interface. The 'Peer Reviews to Arbitrate' item in the navigation sidebar is highlighted with a red box. The main work area shows the same table with one row of data. The search bar and footer are also present.

3. This will direct you to the Patient View where you can see the report. You can view the Peer Review Notes received from both the interpreting radiologist and the peer reviewer.

The screenshot shows the InteleOrchestrator software interface. At the top, there's a navigation bar with links like Home, Worklist, Scheduling, Analytics, Messaging, Reporting, Management, Dashboard, CRM, Protocols, Dictation, Shares, and Logout. Below the navigation bar is a search bar and a 'Worklist' tab. The main area displays a 'Peer Review' section for 'Orolva, DONNELL'. It includes fields for Gender (Female), MRN (CL-192116), DOB (10/12/1992), Age (31y), and MPI (10056). The 'Exam' section shows 'IR PLACEMENT NEPHRO U' with a 'Final Time' of 10/27/2023 at 3:00 PM. The 'Site' section is 'Hospital A' with 'Accounting' details: Phone 123-456-7894, Address 654 W Alton Seattle, WA 97214. The 'Ordering' section is 'Alexandre Messier'. Below this, there's a 'Report' section, an 'Exams' grid, and a 'Notes' section where a note from Hartmann, Jen is shown. A red box highlights the 'Arbitration' section of this note.

#### 4. Click the note marked Arbitration to begin reviewing the Peer Review Note.

This screenshot is similar to the previous one but shows a different peer review note. The 'Notes' section now displays a note from Hartmann, Jen dated 10/19/2023 at 11:47 AM. The note content is identical to the previous screenshot: 'Assigned to: Hartmann, Jen, 10/19/2023 11:47 AM', 'Regression, peerreview, 10/30/2023 7:41 AM', and 'Message: queue, Regression, peerreview, 10/30/2023 7:41 AM'. A red box highlights the 'Arbitration' section of this note.

#### 5. A new window opens where you can accept the peer review rating, overrule the peer review rating, or assign the arbitration to a different Peer Review Administrator.

**Peer Review** **Arbitration**

*Classification: EQUIP 1  
Reason: Ongoing  
Category: (No Category)  
Assigned To: Kiser, Nora at 01/08/2021 6:46 AM  
Created: Wood, Chris, 12/09/2020 2:33 AM*

*Rating: 3b - Finding should usually be made (likely to be clinically significant)  
Message: messed up  
Created: Wood, Chris, 12/09/2020 6:17 AM*

*Rating: 2a - Subtle abnormality (not likely to be clinically significant)  
Message: test  
Reason for Change: test  
Disagreed: Larsen, Sara, 12/24/2020 9:56 AM*

**Close**  **Needs Addendum** **Assign** **Overrule** **Accept** **Comment**

6. To **Accept**: Select the desired Peer Review rating by clicking the checkbox next to the rating you agree with and click **Accept**. The peer review will no longer display in the Worklist.

**Peer Review** **Arbitration**


*Classification: EQUIP 1  
Reason: Ongoing  
Category: (No Category)  
Assigned To: Kiser, Nora at 01/08/2021 6:46 AM  
Created: Wood, Chris, 12/09/2020 2:33 AM*

*Rating: 3b - Finding should usually be made (likely to be clinically significant)  
Message: messed up  
Created: Wood, Chris, 12/09/2020 6:17 AM*


*Rating: 2a - Subtle abnormality (not likely to be clinically significant)  
Message: test  
Reason for Change: test  
Disagreed: Larsen, Sara, 12/24/2020 9:56 AM*

**Close**  **Needs Addendum** **Assign** **Overrule** **Accept** **Comment**


## 7. To Overrule: Click Overrule.



Enter a new rating, enter a comment and click **Overrule** to submit. The peer review will no longer display in the Worklist.



8. To Assign: Click **Assign**.



Select a Peer Review Administrator from the **Assign to Radiologist** drop-down menu. Enter a message and click **Assign** to submit. The peer review will remain in arbitration status and will be assigned to the new Peer Review Administrator.

**Peer Review** **Arbitration**

Classification: EQUIP 1  
Reason: Ongoing  
Category: (No Category)  
Assigned To: Kiser, Nora at 01/08/2021 6:46 AM  
Created: Wood, Chris, 12/09/2020 2:33 AM

Rating: 3b - Finding should usually be made (likely to be clinically significant)  
Message: messed up  
Created: Wood, Chris, 12/09/2020 6:17 AM


Rating: 2a - Subtle abnormality (not likely to be clinically significant)  
Message: test  
Reason for Change: test  
Disagreed: Larsen, Sara, 12/24/2020 9:56 AM

Classification: EQUIP 1

Assign to Radiologist: |

Enter an assignment message

**Close** **Needs Addendum** **Assign**




- When the Peer Review Administrator finalizes a peer review, they have an option to require/request an addendum.

**Peer Review** **Arbitration**


Rating: 1 - Agreement  
Classification: EQUIP4  
Reason: Variance  
Category: Abdomen  
Created: Test123, test123, 03/06/2021 9:16 PM

Rating: 2a - Subtle abnormality (not likely to be clinically significant)  
Message: test  
Reason for Change: test  
Disagreed: Larsen, Sara, 03/19/2021 9:29 AM

**Close** **Needs Addendum** **Assign** **OVERRULE** **ACCEPT** **COMMENT**



- A pop-up box will appear. Including in the Reason for the Addendum is required.



11. Click **OK**.

- Note: Peer review ratings are applied against the originating radiologist, even when an addendum is written.

## Peer Review Scorecard

Radiologists can view their own results via the Peer Review Scorecard widget.

### To view the Peer Review Scorecard:

1. Select Analytics on the InteleOrchestrator home screen.

The screenshot shows the InteleOrchestrator interface. The top navigation bar has tabs: Home, Worklist, Scheduling, **Analytics** (which is highlighted with a red box), and Messaging. The left sidebar is titled 'Worklist' and shows 'My Worklists' with items like 'My Reading Queue' (429 items) and 'My Support Queue' (59 items). The main content area is titled 'Peer Reviews Assigned to Me Worklist (0)' and contains an advanced search form with fields for Keywords, Must Have Text, Created By, Status (set to Pending), Reason, and Category. To the right of the search form, there are filters for 'Assigned On' (Today, 3d, 1w, 08/29/2023) and 'Created On'.

2. Click **Dashboard | Peer Review | Peer Review Scorecard**. The peer review scorecard is displayed in the work area.

The screenshot shows the InteleOrchestrator dashboard with the 'Analytics' tab selected. On the left, there's a sidebar with various widgets like Turnaround, Planning, Productivity, Results Communication, and Peer Review. The 'Peer Review Scorecard' item is highlighted with a red box. The main area displays three charts: 'Stacked Turnaround Time' (No results returned), 'Peer Review Scorecard' (Administrator, Initial) for February 2022, and 'Exam Turnaround Time'. The 'Peer Review Scorecard' chart shows data for the current month and year-to-date across four rating categories (T1, 1, 2a, 2b).

Rating	Reviews You Received	Reviews You Performed	Group at Large
T1	0 0%	0 0%	0 0%
1	0 0%	0 0%	0 0%
2a	0 0%	0 0%	0 0%
2b	0 0%	0 0%	0 0%
Total	0 0	0 0	0 0

Rating	Reviews You Received	Reviews You Performed	Group at Large
T1	0 0%	0 0%	0 0%
1	0 0%	0 0%	0 0%
2a	0 0%	0 0%	0 0%
2b	0 0%	0 0%	0 0%
Total	0 0	0 0	0 0

3. The scorecard shows the reviews received during the current month or year to the current month.

This is a detailed view of the Peer Review Scorecard for the 'Administrator, Initial' role in February 2022. It includes two tables: 'Current Month' and 'Year to This Month', each showing data for four rating categories (T1, 1, 2a, 2b). Below each table is a summary of reviews received and performance metrics. A 'View Messages' button is located at the bottom right.

Rating	Reviews You Received	Reviews You Performed	Group at Large
T1	0 0%	0 0%	0 0%
1	0 0%	0 0%	0 0%
2a	0 0%	0 0%	0 0%
2b	0 0%	0 0%	0 0%
Total	0 0	0 0	0 0

Rating	Reviews You Received	Reviews You Performed	Group at Large
T1	0 0%	0 0%	0 0%
1	0 0%	0 0%	0 0%
2a	0 0%	0 0%	0 0%
2b	0 0%	0 0%	0 0%
Total	0 0	0 0	0 0

*Last updated: 09:04:59 PM*

## Peer Review Analysis Report

The Peer Review Analysis report includes a summary page summarizing peer review results by receiving radiologist, performing radiologist, and modality. It also includes details of each peer review. Note: If you cannot see the report contact your InteleOrchestrator administrator.

**To view the Peer Review Analysis report:**

1. Click **Reports | Peer Review Analysis**. The Peer Review Analysis report is displayed in the Archive section.

Report	Name	Created
Communications Turnaround Time		08/25/2024 1:06 PM
Resident Productivity Report		08/25/2024 12:59 PM
Peer Review Export for eRadPeer		08/25/2024 12:58 PM
<b>Peer Review Analysis</b>		07/02/2024 11:39 AM
Out of Ordering Reading		07/02/2024 11:38 AM
Exam Assignment		07/02/2024 11:37 AM
Manual Rcv History Report		07/01/2024 10:49 AM
Tasks & Report		04/16/2024 9:59 PM
MIPS Reporting		04/15/2024 3:22 PM
GRID Report		04/15/2024 1:58 PM
ACR QCDR Measures Report		04/09/2024 11:49 AM
Peer Review Analysis		02/18/2024 10:37 AM
Verbose Data Export		10/25/2023 11:59 AM
Daily Physician Productivity Summary		

## Creating Custom Reports

You can create custom reports by changing the filters as described below. See [Manage Report Data](#) section for information about the filters.

1. To create custom reports, set the filters as needed.

Report-Peer Review Analysis

Create Report | Save Clear

Options

Anonymize radiologist  
 Include all radiologists

Peer Review

Reason:

Category:

Rating:

Created By:

Assigned:

Enterprise

Site Group:

Site:

Location:

Exam

Site Procedure:

Priority:

Modality:

Exam Date

Exam Time

Dates Relative to Today

Today | 3d | 1w | 1m | 1y

From:

To:

Personnel

Radiologist:

Ordering:

Practice:

Peer Review Date

Created On:


Today | 3d | 1w | 1m | 1y | Absolute

to

Final On:

Today | 3d | 1w | 1m | 1y | Absolute

to



2. Click **Save** to apply the custom filter(s).

Report-Peer Review Analysis-Amol

Create Report **Save** Clear

Options

Anonymize radiologist  
 Include all radiologists

Peer Review

Reason:

Category:

Rating:

Created By:

Assigned:

Enterprise

Site Group:

Site:

Location:

Exam

Site Procedure:  

Priority:   

Modality:   

Exam Date

Exam Time

Dates Relative to Today

Today | 3d | 1w | 1m | 1y 

From: 07/31/2022 

To: 08/24/2022 


Personnel

Radiologist:   

Ordering:  

Practice:  

3. Enter a name for the custom filter and click **Save**.



4. Click **Create Report** to generate the Peer Review Analysis report.

**Report-Peer Review Analysis-Amol**

<b>Create Report</b>	Save	Clear
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Options</b>   <input checked="" type="checkbox"/> Anonymize radiologist             <input type="checkbox"/> Include all radiologists         </div>		
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Peer Review</b>             Reason: <input type="text"/>             Category: <input type="text"/>             Rating: <input type="text"/>             Created By: <input type="text"/>             Assigned: <input type="text"/> </div>		
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Exam</b>             Site Procedure: <input type="text"/>             Priority: <input type="text"/>             Modality: <input type="text"/> </div>		
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Exam Date</b>             Exam Time: <input type="text"/>             Dates Relative to Today: <input type="text"/>             Today   3d   1w   1m   1y   <input type="button" value="X"/>             From: <input type="text" value="07/31/2022"/>             To: <input type="text" value="08/24/2022"/> </div>		
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Enterprise</b>             Site Group: <input type="text"/>             Site: <input type="text"/>             Location: <input type="text"/> </div>		
<div style="border: 1px solid #ccc; padding: 5px;"> <b>Personnel</b>             Radiologist: <input type="text"/>             Ordering: <input type="text"/>             Practice: <input type="text"/> </div>		

## Manage Report Data

To customize Peer Review Analysis report data, you can:

- Anonymize Radiologist
- Select Peer Review type
- Select Enterprise
- Select type of exam
- Set the Exam date and time range
- Select user
- Set Peer Review Date

### Anonymize Radiologist:

You can filter the report to anonymize radiologists or Include all radiologists (include identifying radiologist information).

**Options**

<input type="checkbox"/> Anonymize radiologist
<input type="checkbox"/> Include all radiologists

### Select Peer Review type:

You can filter the Peer Review Reports by selecting peer review details from the drop down menus.

Peer Review

Reason:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Category:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Rating:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Created By:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Assigned:	<input type="button" value="▼"/>

#### Select Enterprise:

You can filter the Peer Review Report by selecting the Site Group, Site, and Location.

Enterprise

Site Group:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Site:	<input type="button" value="▼"/>
Location:	<input type="button" value="▼"/>

#### Select type of exam:

You can filter the Peer Review Report by selecting the type of exams from the drop down menu.

Exam

Site Procedure:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Priority:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>
Modality:	<input type="button" value="▼"/> <input checked="" type="checkbox"/> <input type="button" value="▼"/>

#### Set the Exam date and time range:

You can filter the Peer Review Report by setting the required date range for the report.

The interface shows a dropdown menu labeled "Exam Date". Below it is a dropdown menu labeled "Exam Time". Under "Exam Date", there is a section titled "Dates Relative to Today" with a list of time intervals: "Today", "3d", "1w", "1m", and "1y". To the right of this list is a small "X" button. Below these are two input fields: "From:" and "To:", each with a calendar icon to its right.

#### Select user:

You can filter the Peer Review Report by selecting the type of user or practice.

The interface shows a dropdown menu labeled "Personnel". Below it are three filter sections: "Radiologist:" with a dropdown arrow and a checkmark icon; "Ordering:" with a search icon; and "Practice:" with a dropdown arrow.

#### Set Peer Review Date:

You can filter the Peer Review Report by setting the peer review date range.

Peer Review Date

Created On:

Today | 3d | 1w | 1m | 1y | Absolute

Final On:

Today | 3d | 1w | 1m | 1y | Absolute

# 10

## INTELEORCHESTRATOR ANALYTICS

InteleOrchestrator Analytics provides enterprise-wide operational data and real-time practice information. InteleOrchestrator creates clear visualizations using data from your worklist.

### **In this chapter:**

Analytics Overview .....	223
Opening Analytics .....	223
Understanding the Analytics Page .....	225
Dashboard Widgets .....	226
Reports .....	240
Generating Reports .....	245
MIPS Reporting .....	265

# **Analytics Overview**

InteleOrchestrator Analytics provides enterprise-wide operational data and real-time practice information. InteleOrchestrator creates clear visualizations using data from your worklist. If desired, you can generate reports to share with the rest of your organization.

## **Analyze Your Practice Data**

Instead of relying solely on radiologist's perceptions of where productivity is being lost, InteleOrchestrator software analyzes the practice data. While radiologists can guide the analysis, data always surpasses individual opinions.

Some of the data included in the analysis are:

- Shift and radiologist productivity
- Interruptions
- Interfaces with Operations Staff and Radiology Assistants
- Exam assignment patterns
- Modality, site, subspecialty patterns
- Turnaround Time (TAT) and Service Level Agreement (SLA) compliance patterns
- Inflow of exams and rejection patterns for incomplete exams
- Peer review data as a function of time of day radiologist, and exam type

# **Opening Analytics**


There are two ways to open the Analytics page.

Opening the Analytics page from the Home page:

1. From the Home page, click **Analytics**.



2. The Analytics page is displayed.



Opening the Analytics page from another InteleOrchestrator module:

1. From the list of applications click **Analytics**.

Assignment	AI Priority	Actions	Priority	Status	Modality	Final Time
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Stat	✉️ Unread	XR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Call PE	✉️ Unread	XR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	ASAP	✉️ Unread	XR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Routine	✉️ Unread	MR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Medium	✉️ Unread	XR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Routine	✉️ Unread	XR	
■	🔒 ↴ ⏪	🔗 🔍 🔁 ➔	Routine	✉️ Unread	XR	

2. The Analytics page is displayed.

- If you don't have access to the Analytics Module, please reach out to your InteleOrchestrator administrator for assistance.

## Understanding the Analytics Page

The screenshot shows the InteleOrchestrator Analytics page. On the left, there's a sidebar with two main sections: 'Dashboard' and 'Reports'. Under 'Dashboard', there are several collapsed categories like Turnaround, Planning, Productivity, and Communications. Under 'Reports', there are also collapsed categories like Peer Review and ED Prelim. The main content area is a dark dashboard with various widgets, though they are not clearly legible. At the bottom right, there's a footer with the InteleOrchestrator logo, a notifications icon, and a reminder for '6:17 AM Pacific'.

The Analytics page provides the following tools to view data in InteleOrchestrator:

- Dashboard Widgets - provides a real-time view of exam data.
- Reports - built-in analytic reports with the capability of creating or customizing reports.

# Dashboard Widgets

A widget provides a real-time view of data. Accessing these widgets can assist your administration team in understanding the current operation of your practice, aiding in resource planning. Each widget can be filtered by many variables, including shift. In all cases, the data used to generate the graph can be exported to an Excel-readable format.

Each widget operates on a permission-based system, with permissions configured in the Management module under User Role Management. Contact your system administrator for more information.


## Viewing Widgets

To view available widgets:


1. From the Home page, click **Analytics**.



2. Click **Dashboard**.



3. The list of widgets is displayed.



Each widget is described below:

Widgets	Description
Exam Turnaround Time	<p>Measure exam turnaround time by Modality, Radiologist, Priority, Site and Shift. The turnaround time bubble chart displays three dimensions of data:</p> <ul style="list-style-type: none"><li>• Y-axis: Turnaround time = Final Timestamp - Unread Timestamp (Default) Note: Users can define Exam Turnaround time in the widget.</li><li>• X-axis: Modality, Radiologist, Site, Shift, Priority</li><li>• Bubble size: Exam count</li></ul> <p>Statuses available: Acquired, Begin, Completed, Escalated, Exam Time, Final, Modified, OnPACS, Ordered, Prelim, Signed or</p>

Widgets	Description
	Unread.
Stacked Turnaround Time	Measure average exam turnaround time for multiple segments of the exam life cycle. The stacked turnaround time widget plots average turnaround time by shift, modality, radiologist, site or priority. You can select multiple turnaround time segments. For example. Exam Time to Unread Time and Unread Time to Final Time.
Current Practice	<p>Shows the current number of exams, open communications, and open tasks for a defined set of sites or site groups:</p> <ul style="list-style-type: none"> <li>For exam count, you must specify exam status. Optionally, you can filter by Modality, Site, Site group or Priority.</li> <li>For open communication count, you can select type. If no selection is made all types will be included in the count.</li> <li>For open tasks count, you can select type. If no selection is made all types will be included in the count.</li> </ul>
Exam Backlog	<p>Plot the number of exams in the Unread status at each interval throughout a day or averaged across multiple days.</p> <p> This graph performs a simple count of exams. The export will not include specific exams and the graph cannot be filtered.</p>
Resource Planning	Assess exam count trends throughout the day, week, or year by averaging longer term data. The resource planning histogram displays the distribution of exam count by hour of day, grouped by Radiologist, Modality, Site, Shift, or Subspecialty. This widget is useful when determining staffing needs.
Unread Exam	<p>Plot exam count by Radiologist, Modality, Site, Shift, or Subspecialty with a stacked third dimension showing time spent in a selected status. The unread exam count chart displays two dimensions of data:</p> <ul style="list-style-type: none"> <li>Y-axis: Exam Count = number of exams in unread status.</li> <li>X-axis: Modality, Radiologist, Site, Shift or Subspecialty.</li> </ul>
Exam Count	Plot exam count by Modality, Radiologist, Site Shift or Subspecialty. The exam count chart displays two dimensions of data:

Widgets	Description
	<ul style="list-style-type: none"> <li>• Y-axis: Exam Count = number of exams in Final status.</li> <li>• X-axis: Modality, Radiologist, Site, Shift or Subspecialty.</li> </ul>
Load Balancing Queue Count	<p>Track radiologist's current exam assignment queue for each load balanced assignment rule. The queue count table displays each load balanced assignment rule followed by the radiologists who currently have exams assigned to them by that rule. For each rule, the count of currently unassigned exams is displayed. For each radiologist, the count of exams assigned by that rule and the total count of all assigned exams is displayed.</p>
Practice Productivity	<p>Plot professional RVUs or Work units and exam count by day and month. The Practice Productivity mixed series chart displays three dimensions of data:</p> <ul style="list-style-type: none"> <li>• Left Y-axis: Exam count = number of exams in a Final status.</li> <li>• Right Y-axis: Professional RVUs or Work units = sum of professional relative units for the exams in a Final status.</li> <li>• X-axis: Time = Day of Week or Month of Year.</li> </ul>
Radiologist Productivity	<p>Track radiologist's productivity including exam count, professional RVU or WU and turnaround time. The multi-panel productivity chart displays radiologist productivity as measured by exam count, professional RVUs or Work units count and turnaround time.</p> <ul style="list-style-type: none"> <li>• Exam count = number of exams in a Final status</li> <li>• Total RVUs or Work units = sum of professional relative value units for the exams in a Final status.</li> <li>• Avg TAT = <math>(\text{Final Timestamp} - \text{Unread Timestamp}) / (\text{Exam Count})</math></li> </ul>
SLA Compliance	<p>Shows the percent of exams that were read within the defined service level agreement (configured in SLA Management). SLA Compliance can be plotted as a function of modality, radiologist, shift, site, or priority.</p>
Stacked Exam Count	<p>Plot exam count by day or month with a stacked third dimension by Modality, Radiologist, Site, Shift or Subspecialty. The stacked exam count chart displays three dimensions of data:</p>


Widgets	Description
	<ul style="list-style-type: none"> <li>● Y-axis: Exam Count = number of exams in a Final status.</li> <li>● X-axis: Time = Day of Week or Month of Year.</li> <li>● Stacked by: Modality, Radiologist, Site, Shift or Subspecialty.</li> </ul>
Subspecialty Reading Rate	<p>Shows the percent of exams tagged with a subspecialty that were read by a radiologist within the same subspecialty.</p> <p> Exams with multiple subspecialty tags may be counted more than once. Exams with an unknown radiologist will not be counted.</p>
Auto-next Productivity	<p>Plot monthly radiologist productivity measured by RVU or WU as a percent of exams read in auto-next mode. Each data point shows the total monthly RVU/WU for one radiologist.</p>
Non-reading Task	<p>Plot the number of hours reported for each non-reading task or interruption type. Export to see the breakdown by radiologist, shift, and reading location.</p> <p> The widget is applicable for InteleOrchestrator version 4.3 and earlier.</p>
Communications Turnaround	<p>Measure communication turnaround time by communication type. The communication turnaround time bubble chart displays three dimensions of data:</p> <ul style="list-style-type: none"> <li>● Y-axis: Communications Turnaround time = Completed Timestamp - Created Timestamp.</li> <li>● X-axis: Communication type.</li> <li>● Bubble size: Communication note count</li> </ul>
Peer Review Scorecard	<p>Displays peer review activity for a selected radiologist, for a selected month, and for the year through the selected month. The peer review activity displays peer reviews received and performed by the selected radiologist and compares these to the group at large. You can filter by modality, site and by priority. If no selection is made on these filters, the scorecard will include all modalities, sites and priorities. From the scorecard you can view messages entered on peer reviews for the selected month. The scorecard can be exported to the Excel spreadsheet.</p>

Widgets	Description
	 Users with the 'All radiologist access' user permission can select any user from the radiologist dropdown menu. Without this permission, they will only have current user available for selection.
Peer Review Error Rate by Time of Day	Plot the number of unfavorable peer reviews and corresponding error rate as a function of time of day.
ED Prelim Discrepancy	<p>Plot the number of ED Prelims and the percentage resulting in discrepancy. The ED discrepancy mixed series chart displays three dimensions of data:</p> <ul style="list-style-type: none"> <li>• Left Y-axis: ED Prelim Note count = number of ED Prelim Notes in a Complete status.</li> <li>• Right Y-axis: Percentage with Discrepancy = Percent of ED Prelim Notes in a Complete status that resulted in a discrepancy.</li> <li>• X-axis: Time = Hour of Day, Day of Week, or Month of Year.</li> </ul>


## Applying Widgets

### To add widget to dashboard:

1. Click on a widget name in the widgets panel on the left side of the dashboard.



- The selected widget is displayed on the main section of the dashboard.




- █ The same widget can be selected and added to your dashboard more than once.
- 💡 Add the same widget multiple times with different properties and filters applied.

### To edit widget:

- Click the Edit icon  in the upper right-hand corner of the widget.




- Enter the criteria you would like to view in the widget.



█ Users with only view permissions are currently unable to select options from the site dropdown menu.


3. Click **Save**.




⚠ Caution: Adding excessive widgets that are set to refresh can negatively impact performance. This is especially true when the refresh frequency is high.

## To export widget data:

1. Click the Export  icon.




2. A pop-up will appear.




3. Select an export option:

- Click **Export** to export to an Excel file, a quick view of the information seen in the widget. This view only shows the Summary and the Options tabs.




- The exported data is presented in the UTC time zone.
- Click **Export with Data sheet** to export to an Excel file all information, including the exam details.



- Users must have the Add Exams sheet to widget export permission.

## To delete a widget from the dashboard:

1. Click the Edit icon  in the upper right-hand corner of the widget.



2. Click **Delete**.

**Edit-Exam Turnaround Time**

Properties

Widget: Exam Turnaround Time 

Name: Exam Turnaround Time

Refresh: None

Graph Options

Plot By: Modality

Turnaround Time Definition

Beginning Status: Unread

Ending Status: Final

Exclude Comm Note:

Filter



Modality:  

Site:  

Radiologist:  

Priority:  

Shift: 

Time Range:  to 


Date:

Today | 3d | 1w | 1m | 1y | Absolute  


12/31/2023  to 01/30/2024 

**Cancel** | **Delete**  **Save**

3. A pop-up will appear.



4. Click **Yes**.



- ☞ The Delete Widget is useful to remove any extra copies of a widget when you have added multiple copies of a widget on the dashboard by mistake.

## Reports

The Reports module provides enterprise-wide reporting capabilities, offering a comprehensive view of data across the organization.

Reports can be accessed through the Analytics module. Reports are downloadable in Excel format, can be scheduled for auto-generation, and can be sent via email to recipients in a password-protected format.

Each report is permission-based, and permissions are configured in the Management module under **User Role Management**.

To view available reports:

1. From the Home page, click **Analytics**.



2. Click **Reports**.

The image shows the 'Reports' section of the InteleOrchestrator interface. On the left, a sidebar lists various report categories such as Communications Turnaround Time, Physician Productivity Report, Daily Physician Productivity Summary, and MIPS Reporting. The 'Communications Turnaround Time' report is currently selected and displayed on the right. This screen allows users to configure the report with filters for Site, Location, Laterality, Image Count, Subspecialty, Time Remaining, Accession, Site Accession, Reason, and Priority. It also includes fields for Communication Note, Created By, Type, Status, Exam Date, Exam Time, and Dates Relative to Today. A date range from 'Today' to '1y' is set, and the 'From' and 'To' fields are empty.

3. The list of reports is displayed.

## Reports



- ▶ **Communications Turnaround Time**
- Physician Productivity Report**
- Resident Productivity Report**
- ▶ **Daily Physician Productivity Summary**
- Work Unit Reporting**
- Workload Balancing Report**
- Peer Review Export for eRadPeer**
- Verbose Data Export**
- Peer Review Analysis**
- Out of Ordering Reading**
- Exam Assignment**
- Manual Rvu History Report**
- Tech QA Report**

 The list of reports displayed here is for InteleOrchestrator version 4.5.

The following reports are available for InteleOrchestrator version 4.5, 4.4 and earlier:

Reports	Description
SLA Compliance	This report summarizes exam turnaround time into predefined time brackets and by priority, modality or patient type. The report includes a summary page showing the number of exams that meet each turnaround time bracket. Individual exam details are included on subsequent tabs for all exams and for each exam priority. Available in InteleOrchestrator version 4.4 and earlier.
Communications Turnaround Time	This report provides a detailed breakdown of communication turnaround times for each exam. Turnaround time is presented in

Reports	Description
	the format Days:Hours:Minutes.
Physician Productivity Report	This report provides a detailed breakdown of physician productivity, both overall and categorized by exam and facility. The report includes measurements of RVUs or Work Units and turnaround time.
Resident Productivity Report	This report provides a detailed breakdown of resident productivity, encompassing overall, by exam, and by facility. The report incorporates measurements of RVUs or Work Units and turnaround time.
Daily Physician Productivity Summary	<p>This report offers a summary for each radiologist assigned to a shift, presenting shift details, login session information, the number of exams read, total RVUs generated, and the number of peer reviews performed. The summary table calculates practice totals and averages.</p> <p> A radiologist must be assigned to a shift to be included in this report, and only exams read by a radiologist assigned to a shift are included in the summary.</p>
Work Unit Reporting	This report summarizes procedure work unit totals. Work units are calculated per modality and totaled by site and radiologist. The report includes measurements of WU, RVU, and turnaround time per exam. The values for WU and RVU can be configured in your management pages.
Workload Balancing Report	This report summarizes work units and load balancing gauge results. Results are calculated per radiologist based on the time interval start time. The Work Unit goal per time interval can be configured in your management pages.
Peer Review Export for eRadPeer	<p>This report generates an XML file for upload to eRadPeer™. Only peer reviews meeting the following requirements will be included in the export:</p> <ul style="list-style-type: none"> <li>• Reviewer and reviewee must have a RadPeer ID saved in their profile.</li> <li>• Group ID and Site Code(s) must be defined in the configuration file.</li> <li>• Exam modality must be one of the following: DX, US, CT,</li> </ul>

Reports	Description
Verbose Data Export	<p>MR, NM, MG, IR, PT.</p> <p>This report summarizes all exams within a specified date/time range.</p> <p> Caution: Users should select a limited range to avoid performance impact.</p>
Peer Review Analysis	<p>This report provides detailed information on peer reviews for each exam. The report includes a summary page summarizing peer review results by the receiving radiologist, performing radiologist, and modality. Additionally, it contains a detailed listing of each peer review. Radiologist names can be anonymized.</p>
Out of Order Reading	<p>If users are configured to enter a reason for reading an exam out of order, the Out of Order Reading report will provide a summary showing the number of exams read out of order and the number of exams skipped in auto-next for each radiologist. Exam details are provided on additional tabs. One tab displays the exams read out of order and up to 10 of the exams skipped to read that exam. The other tab shows exams that were marked as skipped in auto-next.</p>
Exam Assignment Summary	<p>This report provides a summary of exam assignments made manually and by the assignment engine.</p>
Manual RVU History Report	<p>This report provides a summary of the RVU and Work Unit credits assigned to RVU tasks.</p>
Tech QA Report	<p>Generates a report with Summary page and an Exam Details page.</p> <p>The Summary page shows a count of Tech QA Notes for each technologist, in total and broken down by modality.</p> <p>The Exam Details page shows every Tech QA Note generated during the report time frame. The details include Location, Site Procedure, Modality, Technologist, Radiologist, Referring Physician, Rating, Message, Created By, Created On, Completed On, Note Completed By, Status, and MRN.</p>

# Generating Reports

You can view results from existing reports, create new reports or create custom reports, based on existing reports.

## Creating a Report

Users can create, save, or schedule reports for automatic execution or email delivery. When a report is saved, the selected filters are stored under the report name in the report table. This allows users to access and download the report at any time.

To create a report:

1. From the Home page, click **Analytics**.



2. Click **Reports**.

The screenshot shows the Analytics dashboard interface. At the top, there is a navigation bar with links: Home, Worklist, Scheduling, Analytics (which is highlighted in blue), Messaging, Reporting, Projects, Management, Dashboard, and Help. Below the navigation bar is a sidebar titled "Widgets" containing sections like Turnaround, Planning, Productivity, Results Communication, Peer Review, and ED Prelim. The main content area is titled "Dashboard" and contains various data visualizations and reports.

3. Select a report from the available list.

The screenshot shows the configuration screen for the "Communications Turnaround Time" report. The left sidebar lists various report options under the "Reports" category, with "Communications Turnaround Time" being the selected item. The main panel displays the report configuration form with fields for Site, Location, Exam, and Communication Note, along with filters for Laterality, Image Count, Subspecialty, Time Remaining, and various date/time parameters.

4. Apply the necessary filters to specify the criteria for data extraction.


This screenshot shows the same report configuration screen as the previous one, but with specific filters applied. The "Laterality" dropdown is set to "Head/Neck". The "Image Count >=" dropdown is set to "10". The "Subspecialty" dropdown is set to "Neurology". The "Time Remaining <=" dropdown is set to "1 min". The "Time Remaining >=" dropdown is set to "0 min". The "Accession" dropdown is set to "All". The "Site Accession" dropdown is set to "All".

## 5. Click **Create Report**.

The screenshot shows the 'Report - Communications Turnaround Time' creation interface. On the left, a sidebar lists various report types. The 'Create Report' button at the top of the main form is highlighted with a red box. The form includes sections for Site, Location, Exam, Communication Note, and Date/Time parameters.

For larger reports, users could receive a pop-up notifying them that the report will take some time to generate. In such cases, they are advised to check back later.

## 6. Upon successful report creation, a confirmation message will appear at the center of the Home page.



## 7. In the confirmation message, click on the InteleOrchestrator **Report** to start the download.

The screenshot shows the 'Analytics' tab selected in the top navigation bar. On the left, a sidebar menu under 'Reports' lists various report types, with 'Communications Turnaround Time' currently selected. The main panel displays the 'Report - Communications Turnaround Time' configuration screen. A modal window titled 'Create Report - Communications Turnaround Time' is open, showing a download link for 'Communications Turnaround Time 08-25-2024 1306.xls' (17 KB in 1 sec). The entire modal window is highlighted with a red box.

8. Once the report is downloaded, you can access the report in your **Downloads** folder.
9. If you did not download the report at the time it was generated, you can select the report from the Archive section.

The screenshot shows the same 'Analytics' tab and 'Communications Turnaround Time' report configuration screen as the previous image. However, the right side of the interface now displays the 'Archive' section. This section contains a table listing three reports: 'Communications Turnaround Time' (created 08/25/2024 1:06 PM), 'Resident Productivity Report' (created 08/25/2024 12:59 PM), and 'Peer Review Export for eRadPeer' (created 08/25/2024 12:58 PM). The entire 'Archive' table is highlighted with a red box.

The archive section displays the last 20 reports generated in total.

10. Click the Excel icon to download the report.

Archive			
	Report	Name	Created
	Communications Turnaround Time		08/25/2024 1:06 PM
	Resident Productivity Report		08/25/2024 12:59 PM
	Peer Review Export for eRadPeer		08/25/2024 12:58 PM

## Save and Schedule a Report

In InteleOrchestrator, you can save reports to be regenerated later. You can also configure InteleOrchestrator to generate reports at set intervals. After being generated,

InteleOrchestrator can send the reports via email, with optional password protection. For extra security, you can choose to save a copy of the report in your Archive section.

**To save a report:**

1. From the Home page, click **Analytics**.



2. Click **Reports**.

The screenshot shows the Analytics dashboard interface. At the top, there is a navigation bar with links: Home, Worklist, Scheduling, Analytics (which is highlighted in blue), Messaging, Reporting, Projects, Management, Dashboard, and Help. Below the navigation bar is a sidebar titled "Widgets" containing sections like Turnaround, Planning, Productivity, Results Communication, Peer Review, and ED Prelim. The main content area is titled "Dashboard" and contains various report cards and summary statistics.

3. Select a report from the available list.

The screenshot shows the configuration screen for the "Communications Turnaround Time" report. The left sidebar lists various reports under the "Reports" category, with "Communications Turnaround Time" being the selected item, indicated by a red box. The main panel displays the report's configuration fields, including dropdown menus for Site, Location, Exam, Laterality, Image Count, Subspecialty, and Time Remaining, along with date and time filters for Exam Date, Exam Time, and Dates Relative to Today.

4. Apply the necessary filters to specify the criteria for data extraction.

This screenshot shows the same report configuration screen as the previous one, but with specific filters applied. The "Site" and "Location" fields are populated with values. The "Laterality" dropdown is set to "Left". The "Image Count >=" dropdown has "1" selected. The "Subspecialty" dropdown has "All" selected. The "Time Remaining <=" dropdown has "0 min" selected. The "Time Remaining >=" dropdown has "0 min" selected. The "Accession" and "Site Accession" fields are empty. The "Communication Note" section includes fields for "Created By" (set to "<Type User>"), "Type" (dropdown menu), and "Status" (dropdown menu). The "Exam Date" section includes "Exam Time" (dropdown menu) and "Dates Relative to Today" (dropdown menu) with "Today" selected. The "From" and "To" date pickers are empty.


5. Click **Save**.

The screenshot shows the 'Report - Communications Turnaround Time' interface. On the left, there's a sidebar with various report options like 'Physician Productivity Report' and 'MIPS Reporting'. In the center, there are search filters for 'Enterprise', 'Site', 'Location', 'Exam', and 'Communication Note'. A large red box highlights the 'Save' button at the top right of the main content area.

6. A pop-up screen will appear.


This is a 'Save Report - Communications Turnaround Time' dialog box. It has sections for 'Name' (with a required field indicator), 'Auto-Generate' (with fields for 'Repeat', 'Days on', 'Every', and 'Active' date), 'Email Options' (with 'Address' and 'Password' fields and a 'Show' checkbox), and an 'Archive' checkbox. At the bottom are 'Cancel' and 'Save' buttons. The 'Save' button is highlighted with a red box.

7. Add a name for the report.




8. Set the frequency for report generation in the repeat field. The frequency can be set to:
  - **No Schedule:** To save the report to be manually generated when desired, in the Auto-Generate section, set Repeat to No Schedule. The report will not be emailed.
  - **Single:** The report is generated only once on the specified date.
    - ⌚ If you save a report with a report interval to today (single repeat), setting the date to today, and this action occurs after the configured Report Start Time, the report may not generate as expected.
  - **Daily:** The report generation recurs daily. The **Every** parameter denotes the number of days. For example, Repeat: Daily + Every: 1 means the report is sent every 1 day.
  - **Every Weekday:** The report is generated every Monday through Friday, weekly.
  - **Weekly:** The report generation recurs weekly.


- **Monthly:** The report generation recurs monthly. The Every parameter denotes the number of months. For example, Repeat: Daily + Every: 1 means the report is sent every 1 month.




9. In the **Email Options** section, add the email address for each recipient of the report.



10. Enter a password to add password protection to the report.



11. Choose **Save sent report in my Archive list** to save a copy of the report to your archive list since the report is only sent to the recipients in the address field.



12. Click **Save** to complete scheduling the report.

**Save Report-Communications Turnaround Time**

Name: \* **Test Report**

Auto-Generate

Repeat: **Weekly**

Days on:  S  M  T  W  T  F  S

Every: **1**

Active: \* **09/25/2024**

Email Options

Address: **abc@gmail.com  
xyz@gmail.com**

Password: **\*\*\*\*\***  Show

Archive:  Save sent report in my Archive list.

**Cancel** **Save**

The **Save** button is highlighted with a red box.

 The same widget can be selected and added to your dashboard more than once.

## Viewing Report History

You can view the history of scheduled reports and all the previously generated reports.

To view the history of reports:

1. In the reports section, expand the report which was saved earlier.

The screenshot shows a software application window with a dark theme. At the top, there are two tabs: "Dashboard" and "Reports". Below these, the word "Reports" is displayed in a larger, bold font. To the right of "Reports" is a small icon. The main content area is titled "Communications Turnaround Time". Under this title, there is a list of items. The first item, "Test Report", is highlighted with a red rectangular box around it. The other items listed are "Test234", "Physician Productivity Report", "Resident Productivity Report", "Daily Physician Productivity Summary", "Work Unit Reporting", "Workload Balancing Report", "Peer Review Export for eRadPeer", "Verbose Data Export", "Peer Review Analysis", and "Out of Ordering Reading".

2. Right click the scheduled report.
3. Click **History**.

Reports



► Communications Turnaround Time

Test Report

Test234



History

Physician Pr.



Edit

Resident Pro



Delete

► Daily Physician Turnaround

Work Unit Reporting

Workload Balancing Report

Peer Review Export for eRadPeer

Verbose Data Export

Peer Review Analysis

► Out of Ordering Reading

4. Details of the scheduled reports are displayed.

Report History - Test Report																																																																				
Generated	Sent	Runtime	Size	Schedule Dates																																																																
02/26/2024	abc@gmail.com xyz@gmail.com	0h 0m 35s	2 KB	<table border="1"> <thead> <tr> <th></th><th>Date</th><th>Day</th><th>ID</th></tr> </thead> <tbody> <tr><td>1</td><td>03/04/2024</td><td>Monday</td><td>5</td></tr> <tr><td>2</td><td>03/11/2024</td><td>Monday</td><td>12</td></tr> <tr><td>3</td><td>03/18/2024</td><td>Monday</td><td>19</td></tr> <tr><td>4</td><td>03/25/2024</td><td>Monday</td><td>26</td></tr> <tr><td>5</td><td>04/01/2024</td><td>Monday</td><td>33</td></tr> <tr><td>6</td><td>04/08/2024</td><td>Monday</td><td>40</td></tr> <tr><td>7</td><td>04/15/2024</td><td>Monday</td><td>47</td></tr> <tr><td>8</td><td>04/22/2024</td><td>Monday</td><td>54</td></tr> <tr><td>9</td><td>04/29/2024</td><td>Monday</td><td>61</td></tr> <tr><td>10</td><td>05/06/2024</td><td>Monday</td><td>68</td></tr> <tr><td>11</td><td>05/13/2024</td><td>Monday</td><td>75</td></tr> <tr><td>12</td><td>05/20/2024</td><td>Monday</td><td>82</td></tr> <tr><td>13</td><td>05/27/2024</td><td>Monday</td><td>89</td></tr> <tr><td>14</td><td>06/03/2024</td><td>Monday</td><td>96</td></tr> <tr><td>15</td><td>06/10/2024</td><td>Monday</td><td>103</td></tr> </tbody> </table>		Date	Day	ID	1	03/04/2024	Monday	5	2	03/11/2024	Monday	12	3	03/18/2024	Monday	19	4	03/25/2024	Monday	26	5	04/01/2024	Monday	33	6	04/08/2024	Monday	40	7	04/15/2024	Monday	47	8	04/22/2024	Monday	54	9	04/29/2024	Monday	61	10	05/06/2024	Monday	68	11	05/13/2024	Monday	75	12	05/20/2024	Monday	82	13	05/27/2024	Monday	89	14	06/03/2024	Monday	96	15	06/10/2024	Monday	103
	Date	Day	ID																																																																	
1	03/04/2024	Monday	5																																																																	
2	03/11/2024	Monday	12																																																																	
3	03/18/2024	Monday	19																																																																	
4	03/25/2024	Monday	26																																																																	
5	04/01/2024	Monday	33																																																																	
6	04/08/2024	Monday	40																																																																	
7	04/15/2024	Monday	47																																																																	
8	04/22/2024	Monday	54																																																																	
9	04/29/2024	Monday	61																																																																	
10	05/06/2024	Monday	68																																																																	
11	05/13/2024	Monday	75																																																																	
12	05/20/2024	Monday	82																																																																	
13	05/27/2024	Monday	89																																																																	
14	06/03/2024	Monday	96																																																																	
15	06/10/2024	Monday	103																																																																	

## Editing the report scheduling

After viewing the report history, you may decide to change the frequency of the report, or the day report is emailed. You can make these changes to report scheduling without having to delete the old schedule and create a new schedule.

To edit the report's schedule:

1. In the reports section, expand the report which you want to edit.

The screenshot shows a software application window with a dark theme. At the top, there are two tabs: "Dashboard" and "Reports". Below these, the word "Reports" is displayed in a larger, bold font. To the right of "Reports" is a small icon. The main content area is titled "Communications Turnaround Time". Under this title, there is a list of items. The first item, "Test Report", is highlighted with a red rectangular box. The other items listed are "Test234", "Physician Productivity Report", "Resident Productivity Report", "Daily Physician Productivity Summary", "Work Unit Reporting", "Workload Balancing Report", "Peer Review Export for eRadPeer", "Verbose Data Export", "Peer Review Analysis", and "Out of Ordering Reading".

2. Right click the scheduled report.
3. Click **Edit**.

The screenshot shows a software interface for managing scheduled reports. At the top, there are links for 'Dashboard' and 'Reports'. Below this, a section titled 'Communications Turnaround Time' is displayed. Underneath, a 'Test Report' is listed with several options: 'Test234' (with a history icon), 'History', 'Physician Pro' (with a red box around it), 'Resident Pro' (with a red box around it), 'Edit' (highlighted with a red box and a pencil icon), 'Daily Physician' (with a trash bin icon), and 'Delete'. Below this section, there are other report types: 'Work Unit Reporting', 'Workload Balancing Report', 'Peer Review Export for eRadPeer', 'Verbose Data Export', 'Peer Review Analysis', and 'Out of Ordering Reading'.

4. Make the desired changes to the schedule.
5. Click **Save**.

#### **Deleting a Scheduled Report:**

If you no longer need to send out a report, you can delete the report schedule from the **Reports** section.

To delete the scheduled report:

1. In the reports section, expand the report which was saved earlier.

The screenshot shows a software application window with a dark theme. At the top, there are two tabs: "Dashboard" and "Reports". Below these, the word "Reports" is displayed in a larger, bold font. To the right of "Reports" is a small icon. The main content area is titled "Communications Turnaround Time". Under this title, there is a list of items. The first item, "Test Report", is highlighted with a red rectangular box. The other items listed are "Test234", "Physician Productivity Report", "Resident Productivity Report", "Daily Physician Productivity Summary", "Work Unit Reporting", "Workload Balancing Report", "Peer Review Export for eRadPeer", "Verbose Data Export", "Peer Review Analysis", and "Out of Ordering Reading".

2. Right click the scheduled report.
3. Click **Delete**.

Reports



Communications Turnaround Time

Test Report

Test234



History

Physician Pro



Resident Pro



Edit

► Daily Physician



Delete

Work Unit Reporting

Workload Balancing Report


Peer Review Export for eRadPeer

Verbose Data Export

Peer Review Analysis

► Out of Ordering Reading

4. In the pop-up dialog, select **Yes** to confirm the deletion or select **No** if you want to keep the report.



## Creating a Custom Report

Contact your Customer Success Manager for more information.

# MIPS Reporting

Merit-based Incentive Payment System (MIPS) reporting involves submitting data on various measures related to patient care, health information technology, practice improvement, and cost with the goal of encouraging high-quality care, interoperability, and overall improvement in healthcare delivery.

The following reports are available:

Reports	Description
ACR Grid Report	Generates the PQRS Grid Report for MIPS Reporting.
ACR QCDR Measures Report	Generates the ACR QCDR Measures Report for MIPS Reporting.

The screenshot shows the InteleOrchestrator software interface. On the left, there's a sidebar with a tree view of reports. The 'MIPS Reporting' node is expanded, and its children, 'GRID Report' and 'ACR QCDR Measures Report', are both highlighted with a red box. The main central area is titled 'Report - GRID Report' and contains various filter and search fields. To the right, there's a preview window with the text 'This report generates the PQRS GRID Report.' and an archive list on the far right showing a history of reports with columns for Report, Name, and Created date.

# 11

## MY SUPPORT QUEUE

My Support Queue is a worklist designed to handle radiologist issues which helps to improve the radiologist productivity. The My Support Queue contains exams with communication Notes assigned to the operations staff. If you are using shift-based worklists, the My Support Queue will be created for you when you log in.

### In this chapter:

About My Support Queue .....	267
Understanding the My Support Queue Interface .....	267
Creating and Editing My Support Queue .....	274

# About My Support Queue

My Support Queue is a worklist designed to handle radiologist issues which helps to improve the radiologist productivity. The My Support Queue contains exams with communication Notes assigned to the operations staff. If you are using shift-based worklists, the My Support Queue will be created for you when you log in.

Some typical communication Notes that you may receive from a radiologist are:

- Prior exams missing
- Needing additional information
- Critical findings
- Contacting the referring physician for more information or follow-up needed

All of these notes are customizable.

All of the worklists contained within My Support Queue are based on communications Notes and not exams. This means that any exams included in the My Support Queue worklists would have to include a communications Note, while the radiologist's My Reading Queue would be exam based and only include exams with no communications Notes or exams with completed communications Notes.

# Understanding the My Support Queue Interface

## To open the My Support Queue:

1. Log into the InteleOrchestrator. For more information, see [Logging in to InteleOrchestrator](#).
2. InteleOrchestrator allows for several different configurations, which may affect the screen that is displayed upon your initial login. If this is the case, from the Home Page, click **Worklist**.



3. Set your shift or set your working hours, click **Save**.
4. The InteleOrchestrator interface displays the My Support Queue.

The screenshot shows the InteleOrchestrator Worklist interface. The top navigation bar includes Home, Worklist (selected), Scheduling, Analytics, Messaging, Reporting, Projects, Management, Dashboard, CRM, and Protocols. The main area is titled "My Support Queue Worklist (59 RVU:119.50 WU:120.00)". On the left, there's a sidebar with "My Suppo" and "My Worklists" sections, including a "Waiting communication notes" section with a count of 59. The main table lists patient information with columns for MRN, Patient Name, Ordering, Type, and Message. The table contains 10 rows of data.

MRN	Patient Name	Ordering	Type	Message
CL-192296P	Terry, LYMAN	Theodore Jones	Missing Images	Missing Coronals
146914S	Lohrer, HARLAN	Alexandre Messier	Standard Call Report	
CL-19239...	Ardalan, CODY	James Tournas	Missing Images	
CL-192203R	Barrios, FELIX	Alexandre Messier	Critical Result	
mmr29144203	Ardalan, CODY	Theodore Jones	Missing Images	Missing Sagittals Missing Coronals
CL-192061D	Lopez, WHITNEY	James Caplan	Missing Images	Missing Sagittals Missing Coronals
mmr30121011	Dhiraj, Dhiraj	PRINCE DIANA	Missing Images	common notes for two exams
mmr30121011	Dhiraj, Dhiraj	PRINCE DIANA	Missing Images	common notes for two exams

Some columns that provide information about the communication notes are:

- The **Type** column displays the type of communication Note attached to the exam.
- The **Message** column displays the message that is included in the communication Note.

My Support Queue Worklist (17 RVU:6.06 WU:12.50)					, BENJAMIN	, JOHN	..EDA..ERLY,	131956, 106664
Patient Name	Radiologist	Referring Physician	Type	Message				
JOHNSON, GEORGE	Administrator, Ini...	BENJAMIN TURNER	Critical Issue with Exam	-asdfsdfasdfsadfasdf				
Fletcher, Sami		Sara Larsen	Critical Result	- Call result into referring physicia...				
RICKERSON, EMILY		MERLIN LEMIEUX	Critical Issue with Exam	-				
TAYLOR, JAMES	Hartmann, Jen	FREDERIC MAURI...	Alert to ED	-				
Smith, Bernice	Larsen, Sara	Mei Li	Critical Result-High Pri...	-test change				
Vigil, Andrew		ADEOLA DARDEN	Alert to ED	-test				

Expand **My Worklist** to locate **My Support Queue** in the worklist panel. You can expand My Support Queue to further see how My Support Queue is broken out into different worklists.

The screenshot shows the 'My Support Queue' interface. On the left, there is a 'Worklist' panel with a tree view. The 'My Worklist' node is expanded, showing the 'My Support Queue' node, which is also expanded. Under 'My Support Queue', there are several categories: 'All Waiting Comm Notes', 'All Critical Results', 'Critical Exam Issues', 'Comm Notes - Completed Today', and 'Incomplete exams'. To the right of the worklist panel is a large table titled 'My Support Queue Worklist (28 RVU:20.79 WU:26.79)'. The table has columns for Site, Location, and Accession. It contains several rows of data, each with a checkbox, a small icon, and text indicating the site (e.g., Clinic C, Hospital A), location (e.g., ED, 3W, NEURO), and accession number (e.g., 10128, MERGED-CT7090, 10765). Below the main table, there are three additional sections: 'All Critical Results', 'Critical Exam Issues', and another section that is partially visible.

	Site	Location	Accession
All Waiting Comm Notes	Clinic C	ED	10128
	Clinic C	3W	MERGED-
	Imaging Center 1		CT7090
	Clinic C	3W	10765
	Hospital A	NEURO	3546972
	Hospital A	ED	10214
	Community Hospital		DX38290
	Hospital A	6NW	21103-1
All Critical Results	Clinic A	ORTH	17010
	Clinic A	ORTH	17010
	Clinic A	ORTH	17010
Critical Exam Issues	Hospital A	ED	24693074
	Clinic A	ORTH	17010

## To open the communication Notes directly from the My Support Queue worklist:

1. Click the Open Notes icon.

Search Results (3 RVU:1.00 WU:3:00)							
	Site	Location	Accession	Site Procedure	MRN	Actions	Patient Name
<input type="checkbox"/>	Hospital B	UNKNOWN	acn0920165...	HAND	mrn0920165...		Study, Frank
<input type="checkbox"/>	Hospital B	UNKNOWN	acn0920165...	HAND	mrn0920165...		Study, Frank
<input type="checkbox"/>	Hospital B	UNKNOWN	CL-203156...	ABD 3 VIEWS NO CXR	CL-192203R		Barrios, FELIX

2. The communication Note is displayed.

Communication-Critical Result

**Waiting**

**Study, Frank**      **mrn0920165149**

**acn0920165149**      **HAND**

Testing critical message

Waiting: Administrator Initial, 03/15/2024 1:54 PM

Type: **Critical Result**

<Note Message>

**Close**      **Comment**      **Attempted**      **Complete**

## To open a communication Note from Patient View:

1. Click on an exam in the worklist.

My Support Queue Worklist (59 RVU:119.50 WU:120.00)				
Worklist	Folder	MRN	Patient Name	Ordering
My Suppo		CL-192296P	Terry, LYMAN	Theodore Jones
		1469145	Lohrer, HARLAN	Alexandre Messier
		CL-19239...	Ardalan, CODY	James Tournas
		CL-192203R	Barrios, FELIX	Alexandre Messier
		mmr29144203	Ardalan, CODY	Theodore Jones
		CL-192061D	Lopez, WHITNEY	James Caplan
		mmr30121011	Dhiraj, Dhiraj	PRINCE DIANA
		mmr30121011	Dhiraj, Dhiraj	PRINCE DIANA


2. Click on the communication note in the notes section of the patient view.

The screenshot shows the 'My Support Queue Worklist' interface. At the top, it displays 'My Support Queue Worklist (17 RVU:6.06 WU:12.50)' and a list of users: BENJAMIN, JOHN, EDA.ERLY, 131956, 106664, Aardvark, Amanda, Ashley, Em, and BECH. Below this, the 'Patient' section shows details for 'JOHNSON, GEORGE': DOB 10/20/1914, Gender Male, Age 108y, MPI 10659, MRN PUI21946. The 'Exam' section is titled 'INJURED BONE' with details: Accession PUI50791, Modality CR - Computed Radiography, Procedure INJURED BONE, Diagnostic Code (empty), and Priority Routine. Below these sections are tabs for 'Report', 'Protocol', 'Dictation', and 'Document'. The 'FINDINGS' section contains a table for 'INJURED BONE' with rows for Radius, Ulna, 1st Metacarpal, 3rd and 5th metacarpals, and a summary row for 'BONE' with columns for STAGE (A - G) and SCORE. The 'Notes' section at the bottom lists three communication notes:

- COMMUNICATION NOTE-Case Issue - Already Dictated  
ghjkjk  
Larsen, Sara, 11/09/2021 00:28
- COMMUNICATION NOTE-Case Issue - Addendum Request  
-high  
Larsen, Sara, 11/09/2021 00:27
- COMMUNICATION NOTE-Case Issue - Addendum Request  
-tytutyuftfyu  
Assigned: Larsen, Sara  
Larsen, Sara, 11/09/2021 00:25


The first note is highlighted with a red box.

3. The communication Note is displayed.




## To work with a communication Note:

1. With the communication Note opened you have the following options:
  - a. Comment - you can let the radiologist know that you are working on this communication Note and have acknowledged it. When you add a comment, this Note will remain in your My Support Queue worklist because it is not completed.
  - b. Attempted - you attempted to find or perform the task requested by the radiologist, but were not able to. When you set the Note to Attempted, this communication Note will remain in your My Support Queue worklist because it is not completed.
  - c. Complete - you have now completed the work.



2. You do have an option to type a message back to the radiologist.
3. With each step (comment, attempted, and complete) you can click close to close the note until you have done more work on this Note.



- Once you set the communication Note to Complete, the communication Note will no longer display as part of your My Support Queue and the exam will automatically display back on the radiologist's worklist.


## Creating and Editing My Support Queue

If you are using shift-based worklists, the My Support Queue will be created for you when you log in. If you are not using shift-based worklists and have not had a My Support Queue created for you, you can create a My Support Queue for yourself.


The My Support Queue is also customizable if you have the permissions to do so.

### To edit or create a new My Support Queue:

- Click the Q icon in the top right of the Worklist panel.




- Select Create a Support Queue.




3. Dialog opens where you can create or change your My Support Queue by selecting specific worklists or you can quickly reload your My Support Queue with the default worklists from another shift.

 Only communication note worklists can be included in My Support Queue.




4. You can do the following:


- Expand the worklist groupings or click **Expand All** to view all worklists and individually add or change your My Support Queue.




- Select **My Shift** to reload the default worklists associated with your current shift.



- Select another shift from the **Shift Reading Queue** drop-down to load the default worklists associated with a different shift.




5. Click Create.




6. Worklists will be added to the My Support Queue section.

## Changing the Properties of My Support Queue

1. Right click on **My Support Queue**.
2. Select **Properties**.




3. Set the Sort Orders, Active Style and optionally you can choose to wrap the My Support Queue content:
  - a. Sort Order - Sort order is the selection to define how you want the communications notes to be sorted on the worklist. For example, if it is set to **Accession**, the worklist would be sorted by number order of accession.



3. Set the Sort Orders, Active Style and optionally you can choose to wrap the My Support Queue content:
  - b. Active style -Active style defines the colors used on the worklist. When **Status** is checked, different communication note statuses are displayed in various colors. If **Status** is not checked, the different communication note statuses are

displayed in white.



4. Click **Change**.

