

PERMISSION GROUP	PERMISSION	SELECT CHECKBOX TO GIVE USER
ANALYTICS (4)	<b>Access to analytics</b>	Access to the Analytics module. The user can click on the Analytics module but there will be no widgets or reports available to select. Permission must be given separately to make each widget and report available for selection.
	<b>Add Exams sheet to widget export</b>	Ability to export a spreadsheet with exam information from widgets. When the user clicks the widget Export icon, they are presented with the <b>Export Option</b> window where they can select Export with Data sheet. The exported spreadsheet contains an extra worksheet detailing the exams that contribute to the widget data.
	<b>Add widgets to dashboard</b>	
	<b>View Custom Reports</b>	Access to the Custom Reports panel in the Reports tab. The panel is not visible without this permission.
ANALYTICS REPORTS (15)	<b>ACR QCDR Measures Report</b>	Ability to generate an ACR QCDR Measures Report from the MIPS Reporting panel.
	<b>Communications Turnaround Time</b>	Ability to generate a Communications Turnaround Time report.
	<b>Daily Physician Productivity Summary</b>	Ability to generate a Daily Physician Productivity Summary report.
	<b>Exam Assignment</b>	Ability to generate an Exam Assignment report.
	<b>GRID Report</b>	Ability to generate a GRID Report from the MIPS Reporting panel.
	<b>Manual Rvu History Report</b>	Ability to generate a Manual Rvu History Report.
	<b>Out of Ordering reading</b>	Ability to generate an Out of Ordering reading report.
	<b>Peer Review Analysis</b>	Ability to generate a Peer Review Analysis report.
	<b>Peer Review Export for eRadPeer</b>	Ability to generate a Peer Review Export for eRadPeer report.
	<b>Physician Productivity Report</b>	Ability to generate a Physician Productivity Report.
	<b>Resident Productivity Report</b>	Ability to generate a Resident Productivity Report.
	<b>Tech QA Report</b>	Ability to generate a Tech QA Report.
	<b>Verbose Data Export</b>	Ability to generate a Verbose Data Export report.
	<b>Work Unit Reporting</b>	Ability to generate a Work Unit Reporting report.
	<b>Workload Balancing Report</b>	Ability to generate a Workload Balancing Report.
ANALYTICS WIDGETS (19)	<b>Auto-next Productivity</b>	Ability to generate an Auto-next Productivity widget.
	<b>Communications Turnaround</b>	Ability to generate a Communications Turnaround widget.
	<b>Current Practice</b>	Ability to generate a Current Practice widget.
	<b>ED Prelim Discrepancy</b>	Ability to generate a Prelim Discrepancy widget.
	<b>Exam Backlog</b>	Ability to generate an Exam Backlog widget.
	<b>Exam Count</b>	Ability to generate an Exam Count widget.
	<b>Exam Turnaround Time</b>	Ability to generate an Exam Turnaround Time widget.
	<b>Load Balancing Queue Count</b>	Ability to generate a Load Balancing Queue Count widget.
	<b>Non-reading task</b>	Deprecated.
	<b>Peer Review Error Rate by Time of Day</b>	Ability to generate a Peer Review Error Rate by Time of Day widget.
	<b>Peer Review Scorecard</b>	Ability to generate a Peer Review Scorecard by Time of Day widget.
	<b>Practice Productivity</b>	Ability to generate a Practice Productivity widget.
	<b>Radiologist Productivity</b>	Ability to generate a Radiologist Productivity widget.
	<b>Resource Planning</b>	Ability to generate a Resource Planning widget.
	<b>SLA Compliance</b>	Ability to generate an SLA Compliance widget.
	<b>Stacked Exam Count</b>	Ability to generate a Stacked Exam Count widget.
	<b>Stacked Turnaround Time</b>	Ability to generate a Stacked Turnaround Time widget.

	<b>Subspecialty Reading Rate</b>	Ability to generate a Subspecialty Reading Rate widget.
	<b>Unread Exam</b>	Ability to generate an Unread Exam widget.
CRM (12)	<b>Access to CRM</b>	View only access to the Contact Management, Portal User Role Management and Portal Field Editor pages in the CRM module. From Portal Exam Lists Management, the user can add/edit/delete Exam Lists.
	<b>Create and Edit Contacts</b>	Ability to add/edit/delete/merge contacts and import/export contacts in the Contact Management page.
	<b>Create and Edit Physician Groups</b>	Ability to add/edit/delete ordering physician groups. From the <b>Edit Contact</b> window, the user can click an add icon beside groups. <b>Note:</b> The user must have the <b>Edit Contacts</b> permission to access the <b>Edit Contact</b> window.
	<b>Edit Portal Forms Management</b>	
	<b>Edit Portal Settings</b>	Ability to edit portal settings for a contact in the Contact Management page.
	<b>Edit Portal User Roles</b>	Ability to add/edit portal user roles from the Portal User Role Management page.
	<b>Edit Site</b>	Ability to edit site information and import/export site details and site departments in the Contact Management page.
	<b>Portal forms field editor edit</b>	Ability to add groups and fields to forms in the Portal Field Editor page.
	<b>Portal forms field editor view</b>	
	<b>View Contacts</b>	
COMMON (5)	<b>View Portal Forms Management</b>	
	<b>View Portal User Roles</b>	
	<b>Access to MFA login</b>	Multifactor Authentication challenge at login. <b>Note:</b> Multifactor Authentication for Login must first be enabled for the system.
	<b>Bypass enterprise authentication</b>	Ability to use local/LDAP authentication when enterprise authentication has been enabled. <b>Note:</b> Permission only applicable to customers with an enterprise authentication integration enabled.
	<b>Change exam status</b>	Ability to change the status of exams. The user can select Change Status from the right-click/Actions menus to open the <b>Change Exam Status</b> window.
	<b>Concurrent user logins (multiple logins)</b>	Ability to login on different workstations at the same time. The user can login on a second workstation without being logged out on the first workstation.
	<b>View button to set away and lock system</b>	Lock System (Away) button for selection.
DASHBOARD (13)	<b>Access to Dashboard</b>	Access to the Dashboard module. The user can click on the Dashboard module but there will be no tabs available to select. Permission must be given separately to make each tab visible.
	<b>Actions on DICOM Routing Page</b>	Ability to edit DICOM routing information. It is not necessary to also give the user the <b>View DICOM Routing</b> permission.
	<b>Actions on Extender Logs Page</b>	Ability to view and delete extender logs. It is not necessary to also give the user the <b>View uploaded Extender logs</b> permission.
	<b>Actions on HL7 Sender Page</b>	Ability to resend HL7 messages. It is not necessary to also give the user the <b>View HL7 Sender</b> permission.
	<b>Actions on Linq Workflow Page</b>	Ability to activate/deactivate Linq connections. It is not necessary to also give the user the <b>View Linq Workflow</b> permission.
	<b>Actions on User Monitoring Page</b>	Ability to Log off/Pause exam assignment for users. It is not necessary to also give the user the <b>View User Monitor</b> permission.
	<b>View DICOM Routing</b>	View only access to the DICOM Routing tab.

	<b>View HL7 Sender</b> <b>View Linq Workflow</b> <b>View Services</b>	View only access to the HL7 Sender tab. View only access to the Linq Workflow tab. Access to the Service and Service Logs tabs. The user can right-click and perform actions from the PeerReviewAssignment service in the Service tab.
	<b>View uploaded Extender logs</b> <b>View User Monitor</b> <b>View Worklist Performance</b>	View only access to the Extender Logs tab. View only access to the User Monitoring tab. Access to the Worklist Performance tab.
ERROR LOG (2)	<b>Admin error logging with backtrace</b>	Access to log.php page including backtrace. It is not necessary to also give the user the <b>Show only error log info</b> permission.
MANAGEMENT (38)	<b>Show only error log info</b>	Access to log.php page.
	<b>Access Log Analyzer</b>	Access to the Log Analyzer page.
	<b>Access Peer Review Management</b>	Access to the Peer Review Management page.
	<b>Access Procedure Management</b>	Access to the Procedure Management page.
	<b>Access Registry Configuration</b>	Access to the Registry Configuration page.
	<b>Access to admin help pages</b>	Access to administrative level help pages. When the user clicks Help, they will be directed to the Clario Administration Modules Online Help in place of the Clario SmartWorklist Online Help.
	<b>Access to management</b>	Access to the Management module. The user can click on the Management module but there will be no pages available to select. Permission must be given separately to make each page visible.
	<b>Access to Navigator Management</b>	Access to the AI Navigator page.
	<b>Access to Send HL7 message</b>	Ability to send HL7 message in Assignment Rule Management. <b>Note:</b> Additional back end configuration is required.
	<b>Admin Access to Credentialing</b>	Ability to edit user credentialing in the User Credentialing page. It is not necessary to also give the user the <b>View Access to Credentialing</b> permission.
	<b>Change exam priority</b>	Deprecated.
	<b>Create new exam orders</b>	
	<b>Create or edit a user</b>	Ability to add/edit local users in the User Management page. It is not necessary to also give the user the <b>View user</b> permission. <b>Note:</b> This permission applies to adding local users and not adding users via an LDAP integration.
	<b>Create or edit user roles</b>	Ability to add/edit/delete user roles in the User Role Management page. It is not necessary to also give the user the <b>View user roles</b> permission.
	<b>Delete exam information</b>	
	<b>Delete user from User Management</b>	Ability to delete users in the User Management page. <b>Note:</b> This permission applies to deleting local users and not users added via an LDAP integration.
	<b>Edit all Exam related data</b>	Ability to edit exam data. The user can click the exam name in patient view and edit fields in the Information panel of the <b>Exam History</b> window. It is not necessary to also give the user the <b>Edit some exam related data</b> permission. <b>Note:</b> Additional back end configuration is required.
	<b>Edit all Patient related data</b>	Ability to edit patient data. The user can click the patient name in patient view and edit fields in the Information panel of the <b>Patient History</b> window. It is not necessary to also give the user the <b>Edit some patient related data</b> permission. <b>Note:</b> Additional back end configuration is required.

	<b>Edit Assign Rules</b>	Ability to add/edit assignment rules in the Assign Rule Management page. It is not necessary to also give the user the <b>View Assign Rule Management</b> permission.
	<b>Edit exam information</b>	
	<b>Edit Group Management page</b>	Ability to add/edit groups in the Group Management page. It is not necessary to also give the user the <b>View Group Management</b> permission.
	<b>Edit Notification Rules</b>	Ability to add/edit notification rules in the Notification Management page. It is not necessary to also give the user the <b>View Notification Management</b> permission.
	<b>Edit SLA Rules</b>	Ability to add/edit SLA rules in the SLA Management page. It is not necessary to also give the user the <b>View SLA Management</b> permission.
	<b>Edit some exam related data</b>	Ability to edit the following exam data fields: Status, Modality, Reason, Begin, PReport, Site Accession, Radiologist, Attending. The user can click the exam name in patient view and edit fields in the Information panel of the <b>Exam History</b> window. <b>Note:</b> Additional back end configuration is required.
	<b>Edit some patient related data</b>	Ability to edit the following patient data fields: Master Patient ID, MRN, Last Name, Gender, Has Patient Note. The user can click the patient name in patient view and edit fields in the Information panel of the <b>Patient History</b> window. <b>Note:</b> Additional back end configuration is required.
	<b>Edit Username in Password Management</b>	Ability to edit own application username in Password Management. <b>Caution:</b> Users should not edit their application usernames if using an LDAP integration.
	<b>Edit/add/import diagnostic codes</b>	Ability to import diagnostic code in Diagnostic Code tab of Procedure Management. User is able to add/edit diagnostic codes with the <b>Access Procedure Management</b> permission.
	<b>Manage System Configuration</b>	Access to the Configuration page.
	<b>Match patients</b>	Ability to confirm or reject Patient Matches and perform Manual Matches.
	<b>Merge patients</b>	Ability to perform Manual Merges. <b>Caution:</b> A merge will overwrite patient data.
	<b>Receive notification of upcoming SSL Certificate expiry</b>	A notification upon every login when the SSL certificate will expire soon (a default of 30 days).
	<b>View Access to Credentialing</b>	View only access to the User Credentialing page.
	<b>View Assign Rule Management</b>	View only access to the Assign Rule Management page.
	<b>View exam information</b>	
	<b>View Group Management</b>	View only access to the Group Management page.
	<b>View Notification Management</b>	View only access to the Notification Management page.
	<b>View SLA Management</b>	View only access to the SLA Management page.
	<b>View user roles</b>	View only access to the User Role Management page.
	<b>View users</b>	View only access to the User Management page. The user can Clone users and edit Contact Info.
MESSAGE (6)	<b>Ability to Delete Private Groups</b>	Ability to delete private groups. The user can right-click on the group name and select Delete.
	<b>Ability to Delete Public Groups</b>	Ability to delete public groups. The user can right-click on the group name and select Delete.
	<b>Access to Messaging</b>	Access to the Messaging module.
	<b>Leave Private Group</b>	Ability to leave private groups. The user can right-click on the group name and select Leave.
	<b>Leave Public Group</b>	Ability to leave public groups. The user can right-click on the group name and select Leave.
	<b>Messaging Admin</b>	Ability to create public groups and visibility of all private groups.
MULTI-SYSTEM (2)	<b>System Administrator</b>	Ability to manage users using Linq worklist sharing. The Multi System panel in User Management is not visible without this permission.

	<b>System User</b>	Access to Linq worklist sharing. <b>Note:</b> The user must have this permission in both the native and the shared Clario system.
NOTES (28)	<b>Access to PRT Notes in IntelViewer</b>	
	<b>Administrator Communication Notes</b>	Ability to delete Communication Notes.
	<b>Administrator ED Prelim Notes</b>	Ability to delete ED Prelim Notes.
	<b>Administrator Exam Notes</b>	Ability to delete Exam Notes.
	<b>Administrator Follow-up Notes</b>	Ability to delete Follow-up Notes.
	<b>Administrator Patient Notes</b>	Ability to delete Patient Notes.
	<b>Administrator Peer Review Notes</b>	
	<b>Administrator Upload file</b>	Ability to delete Upload Files.
	<b>Create Communication notes</b>	Ability to create Communication Notes. It is not necessary to also give the user the <b>View Communication notes</b> permission.
	<b>Create ED Prelim notes</b>	Ability to create ED Prelim Notes. It is not necessary to also give the user the <b>View ED Prelim notes</b> permission.
	<b>Create Exam notes</b>	Ability to create Exam Notes. It is not necessary to also give the user the <b>View Exam notes</b> permission.
	<b>Create Follow-up notes</b>	Ability to create Follow-up Notes. It is not necessary to also give the user the <b>View Follow-up notes</b> permission.
	<b>Create Patient notes</b>	Ability to create Patient Notes. It is not necessary to also give the user the <b>View Patient notes</b> permission.
	<b>Create Peer Review notes</b>	Ability to create Peer Review Notes.
	<b>Create Upload file</b>	Ability to create Upload Files. It is not necessary to also give the user the <b>View Upload file</b> permission.
	<b>ED Prelim notes Acknowledge</b>	Ability to acknowledge a discrepancy on an ED Prelim note and change the note status to complete.
	<b>ED Prelim notes Agree</b>	Ability to agree with an ED Prelim note and change the note status to complete.
	<b>ED Prelim notes Comment Discrepancy</b>	Ability to comment on an ED Prelim note in discrepancy status.
	<b>ED Prelim notes Comment Pending</b>	Ability to comment on an ED Prelim note in pending status.
	<b>ED Prelim notes Major Discrepancy</b>	Ability to flag a major discrepancy on an ED Prelim note and change the note status to discrepancy.
	<b>ED Prelim notes Minor Discrepancy</b>	Ability to flag a minor discrepancy on an ED Prelim note and change the note status to discrepancy.
	<b>View Communication notes</b>	Ability to view Communication Notes.
	<b>View ED Prelim notes</b>	Ability to view ED Prelim Notes.
	<b>View Exam notes</b>	Ability to view Exam Notes.
	<b>View Follow-up notes</b>	Ability to view Follow-up Notes.
	<b>View Patient notes</b>	Ability to view Patient Notes.
	<b>View Peer Review notes</b>	Ability to view Peer Review Notes and Peer Review worklists and search for Peer Review Notes. <b>Note:</b> The user must have the <b>Limited Peer Review Admin</b> or <b>Peer review administrator</b> permission for visibility of Peer Review Notes submitted by other users.
	<b>View Upload file</b>	Ability to view Upload Files.

PATIENT VIEW (5)	<b>Associate/disassociate exams to diagnostic codes</b>	
	<b>Final report upload admin</b>	Ability to manually upload reports and delete manually uploaded reports. The user can select the Manually Upload Report icon in the report panel. It is not necessary to also give the user the <b>Final report upload user</b> permission.
	<b>Final report upload only</b>	Ability to manually upload reports. The user can select the Manually Upload Report icon in the report panel.
	<b>Final report upload user</b>	Ability to manually upload reports and delete own manually uploaded reports. The user can select the Manually Upload Report icon in the report panel. It is not necessary to also give the user the <b>Final report upload only</b> permission.
	<b>View/Export detailed STAMP from exam history</b>	Visibility of the Detail panel of the <b>Exam History</b> window and ability to download HL7 messages from the <b>Exam History</b> window. The user can select the Show All Columns checkbox.
PEER REVIEW (4)	<b>Delete Peer Review note</b>	Ability to delete Peer Review notes.
	<b>Limited Peer Review Admin</b>	Ability to view Peer Review Notes without user information. The user will be blinded to peer reviews. <b>Note:</b> The user must have the <b>View Peer Review notes</b> permission to access Peer Review worklists and search for Peer Review notes.
	<b>Peer Review user</b>	Ability to perform peer reviews. The user will be assigned peer reviews if configured and will be blinded to peer reviews. <b>Note:</b> The user must have the <b>View Peer Review notes</b> permission to access Peer Review worklists.
	<b>Peer review administrator</b>	Ability to arbitrate peer reviews and view Peer Review notes with user information. The user will not be blinded to peer reviews. <b>Note:</b> The user must have the <b>View Peer Review notes</b> permission to access Peer Review worklists and search for Peer Review notes.
PREFERENCES (3)	<b>Prompt for shift</b>	A prompt upon login to select shift. When the Smart Worklist module is launched, the user will be presented with a <b>Select Current Shift</b> window. If this permission is selected in addition to the <b>Required Reading Location</b> permission, the user will be presented with a <b>Select Current Shift and Location</b> window.
	<b>Require Shift Selection</b>	A prompt upon login to select shift that cannot be cancelled. When the Smart Worklist module is launched, the user will be presented with a <b>Select Current Shift</b> window. The user must select a shift and click Save to proceed to the Smart Worklist. If this permission is selected in addition to the <b>Required Reading Location</b> permission, the user will be presented with a <b>Select Current Shift and Location</b> window. <b>Note:</b> The user must have the <b>Prompt for shift</b> permission.
	<b>Required Reading Location</b>	A prompt upon login to select reading location that cannot be cancelled. When the Smart Worklist module is launched, the user will be presented with a <b>Select Current Location</b> window. The user must select a reading location and click Save to proceed to the Smart Worklist. If this permission is selected in addition to the <b>Prompt for shift</b> permission, the user will be presented with a <b>Select Current Shift and Location</b> window.
PROJECT (2)	<b>Access to Project</b>	Access to the Projects module. The user can add project notes, mark tasks as complete and download task attachments.
	<b>Administrator</b>	Ability to add/edit/delete projects, add/edit/delete tasks and import/export projects.
REPORTING (8)	<b>Access to Dictation application</b>	Access to the Dictation module.



	<b>Access to Protocoling application</b>	Access to the Protocoling module. The user can add/edit/delete templates.
	<b>Access to Reporting application</b>	Access to the Reporting module. The user can add/edit/delete templates and import/export templates.
	<b>Edit template remarks</b>	Ability to add/edit/delete template remarks. The user can select Add Template Remark by right-clicking the report to open the <b>Remark</b> window. <b>Note:</b> Permission should only be enabled in demo mode, feature has not been completed.
	<b>Edit templates</b>	
	<b>Open exams for protocoling</b> <b>View exam protocols</b> <b>View templates filtered for reading radiologist</b>	Ability to launch protocoling. The user can click the Open for Protocoling icon in the worklist.  Access to reporting radiologists' templates rather than own templates when performing transcription. <b>Note:</b> Permission intended for transcriptionist user.
RESIDENT/ATTENDING (6)	<b>Ability to edit Resident/Attending queues</b>	Ability to edit the Resident Queue, Attending Queue and Attending Pool. The user can right-click the Resident Queue/Attending Queue/Attending Pool and click Edit to change the search parameters of the worklist. <b>Note:</b> Additional permission/s required to enable Edit option.
	<b>Access to Attending Pool</b>	Ability to add the Attending Pool to My Worklist. The user can click the Create Attending Pool icon.
	<b>Access to Attending Queue</b>	Ability to add the Attending Queue to My Worklist. The user can click the Create Attending Queue icon.
	<b>Access to Resident Queue</b>	Ability to add the Resident Queue to My Worklist. The user can click the Create Resident Queue icon.
	<b>Attending Radiologist</b>	Recognition as an attending radiologist and visibility in the Attending Radiologist dropdown menu in Advanced Search.
	<b>Prompt for Attending Radiologist</b>	A prompt upon login to select Attending Radiologist. When the Smart Worklist module is launched, the user will be presented with a <b>Select Attending Radiologist</b> window.
SCHEDULE (5)	<b>Resident Radiologist</b>	Recognition as a resident radiologist and visibility in the Resident Radiologist dropdown menu in Advanced Search.
	<b>Access to Calendar</b>	Access to the Calendar page.
	<b>Access to Demand Matched Scheduling</b>	Access to the Demand Matched Scheduling panel in the Shift Management page.
	<b>Access to Shift Management</b>	Access to the Shift Management page. The user can add/edit/delete shifts. It is not necessary to also give the user the <b>Read-only access to Shift Management</b> permission.
	<b>Application Section Schedule access</b>	Access to the Scheduling module. The user can click on the Scheduling module but there will be no pages available to select. Permission must be given separately to make each page available for selection.
SEARCH ACCESS (1)	<b>Read-only access to Shift Management</b>	View only access to the Shift Management page.
	<b>All radiologist access</b>	Ability to select any user from the Radiologist dropdown menu in Advanced Search and Analytics. Without this permission, the user will only have Current User available for selection. <b>Note:</b> The user is able to select any user from the Resident Radiologist, Attending Radiologist and Assigned (U) dropdown menus in advanced search without this permission.
	<b>Create a task</b>	Ability to create tasks from the status bar. It is not necessary to also give the user the <b>View tasks</b> permission.
	<b>Create alerts</b>	Ability to create alerts from the status bar. It is not necessary to also give the user the <b>View alerts</b> permission.
STATUS BAR (6)		

	<b>Create notifications</b>	Ability to create notifications from the status bar. It is not necessary to also give the user the <b>View notifications</b> permission.
	<b>View alerts</b>	Ability to view alerts from the status bar.
	<b>View notifications</b>	Ability to view notifications from the status bar.
	<b>View tasks</b>	Ability to create tasks from the status bar.
TEACHING NOTE (3)	<b>Administrator Teaching notes</b>	Ability to delete Teaching Notes.
	<b>Create Teaching notes</b>	Ability to create Teaching Notes. It is not necessary to also give the user the <b>View Tech QA notes</b> permission. <b>Note:</b> The user must also be given the <b>Use Folders</b> permission to associate Teaching Notes with Folders.
	<b>View Teaching notes</b>	Ability to view Teaching Notes.
TECH QA (3)	<b>Administrator Tech QA notes</b>	Ability to delete Tech QA Notes.
	<b>Create Tech QA notes</b>	Ability to create Tech QA Notes. It is not necessary to also give the user the <b>View Tech QA notes</b> permission.
	<b>View Tech QA notes</b>	Ability to view Tech QA Notes.
WORKLIST (40)	<b>Access to Document Viewer</b>	Access to Document Viewer tools and priors menu and ability to view documents. Without this permission the user can see the Document tab but it will read <i>No Documents Available</i> . <b>Note:</b> Permission only applicable to customers with <b>Clario Document Viewer</b> enabled.
	<b>Access to worklist</b>	Access to the Smart Worklist module.
	<b>Assign exams to one self</b>	Ability to assign exams to themselves. The user can click a Lock Exam to Me icon in the Smart Worklist.
	<b>Assign exams to users or groups</b>	Ability to assign exams to users or groups. The user can click an Assign/Lock Exam icon in the Smart Worklist or select Assign from the right-click/Actions menus to open the <b>Assign</b> window. <b>Note:</b> From the <b>Assign</b> window, the user can assign exams to themselves.
	<b>Assign to Peer Review</b>	Ability to assign exams to users to peer review. The user can select Assign to Peer Review from the right-click/Actions menus to open the <b>Assign Peer Review</b> window. <b>Note:</b> The user must also be given the <b>All radiologist access</b> permission to assign exams to peer review to users other than themselves.
	<b>Assign worklist activations</b>	Ability to configure threshold counts and active time when creating a worklist.
	<b>Cancel Exam</b>	Ability to cancel exams. The user can select Cancel Exam from the right-click/Actions menus to open the <b>Cancel Exam</b> window.
	<b>Change Priority</b>	Ability to change the priority of exams. The user can select Change Priority from the right-click/Actions menus to open the <b>Change Exam Priority</b> window.
	<b>Combine exams after dictation</b>	Ability to combine reports after dictation. <b>Note:</b> Permission only applicable to customers with Clario Reporting enabled.
	<b>Configurable HL7 send</b>	
	<b>Create My Reading Queue</b>	Ability to create/edit My Reading Queue. The user can click the Create Queue icon to open the <b>Create My Reading Queue</b> window.
	<b>Create My Support Queue</b>	Ability to create/edit My Support Queue. The user can click the Create Queue icon to open the <b>Create My Support Queue</b> window.
	<b>Create RVU/WU Credit</b>	Ability to enter RVU/WU Credit. The user can click the Enter RVU/WU Credits icon to open the <b>RVU/WU Credit</b> window.
	<b>Create Worklists</b>	Ability to create personal worklists. The user can click Add to Worklist after performing an advanced search but can only add the worklist to their My Worklist worklist group.



<b>Create exam report in patient view</b>	<p>Visibility in the Radiologist filter in the Reporting module.</p> <p><b>Note:</b> Permission only applicable to customers with Clario Reporting enabled.</p>
<b>Create/Edit group and enterprise worklists</b>	<p>Ability to create enterprise worklists. The user can click Add to Worklist after performing an advanced search and add the worklist to an enterprise worklist group.</p> <p><b>Note:</b> The user must have the <b>Create Worklists</b> permission to create worklists.</p>
<b>DICOM Route</b>	<p>Ability to trigger resending of images. The user can select DICOM Route/DICOM Route and Read from the right-click/Actions menus.</p> <p><b>Note:</b> Additional back end configuration is required.</p>
<b>Delete Exam</b>	<p>Ability to delete exams. The user can select Delete Exam from the right-click/Actions menus.</p>
<b>Disassociate combined exams</b>	<p>Ability to uncombine exams after they've been combined. The user can right-click on a combined exam and select Disassociate Exam to open the <b>Disassociate Combined Exam</b> window.</p>
<b>Edit Exam Group Tagging</b>	<p>Ability to change the subspecialty tag of exams. The user can select Edit Subspecialty from the right-click/Actions menus to open the <b>Edit Exam Group Tagging</b> window.</p>
<b>Edit exam report in transcription</b>	<p>Ability to transcribe reports and send reports for signature when Clario Reporting is enabled. The user can launch exams in Sent for Transcription status for transcription and select Complete Transcription.</p> <p><b>OR</b> Ability to launch exams in Sent for Transcription status for transcription when PowerScribe is integrated.</p> <p><b>Note:</b> Permission intended for transcriptionist user.</p> <p><b>Caution:</b> This permission conflicts with the <b>Send exam report to transcription</b> permission and the <b>Launch dictation or voice recognition</b> permission.</p>
<b>Emergency Access for Exams</b>	<p>Ability to grant temporary access of exams to portal users by site.</p>
<b>Enable subscribe for exams and worklists</b>	<p>Ability to subscribe to exams. The user can select the Subscribe Exam icon in the worklist, the Subscribe to worklist icon in Auto-Next or right-click on a worklist and select Subscribe.</p> <p><b>Note:</b> Additional back end configuration is required.</p>
<b>Launch Viewer and image viewers</b>	<p>Ability to launch viewer. The user can click the Open in Viewer icon in the worklist or select Open in Viewer from the right-click/Actions menus.</p> <p><b>Note:</b> Additional back end configuration is required.</p>
<b>Launch dictation or voice recognition</b>	<p>Ability to launch dictation. The user can click the Open in Dictation/Read Exam icons in the worklist or select Open in Dictation/Read Exam from the right-click/Actions menus.</p> <p><b>Note:</b> Additional back end configuration is required.</p>
<b>Launch patient view applications</b>	
<b>Merge exams</b>	<p>Ability to merge exams. The user can select Merge Exams from the Actions menus to open the <b>Merge Exams-Select Primary</b> window.</p>
<b>Needs Addendum</b>	<p>Ability to change exam status to Addendum. The user can select Needs Addendum from the right-click/Actions menus.</p> <p><b>Note:</b> The user can also change exam status to Addendum with the <b>Change Exam Status</b> permission.</p>
<b>Print Report</b>	<p>Ability to print reports. The user can select Print Report from the right-click/Actions menus to open the report in a printable PDF format.</p> <p><b>Note:</b> The user can select the Print PDF Report icon in the report panel without this permission.</p>

<b>Send Report</b>	Ability to fax/email reports. The user can select Send Report from the right-click/Actions menus to open the <b>Send Report</b> window. <b>Note:</b> Additional back end configuration is required.
<b>Send exam report to transcription</b>	Ability to send reports for transcription. The user can toggle between Voice Recognition (VR) and Transcription (TX) and when using Transcription (TX) can select Send to Transcription. <b>Note:</b> Permission intended for radiologist user. <b>Note:</b> Permission only applicable to customers with Clario Reporting enabled. <b>Caution:</b> This permission conflicts with the <b>Edit exam report in transcription</b> permission.
<b>Set Auto-next Columns</b>	Ability to display column headers in Auto-Next and add/remove columns. The user can click the icon to toggle column headers on/off.
<b>Set exam visibility in search (hide/unhide exams)</b>	Ability to hide exams. The user can click a Hide/Show study icon in the worklist or select Hide/Show Study from the right-click/Actions menus then select the Hide studies icon.
<b>Sign exam report</b>	Ability to sign reports when Clario Reporting is enabled. The user can select Sign Report. <b>Note:</b> The user can select Sign Prelim without this permission. <b>OR</b> Recognition as an attending radiologist when PowerScribe is integrated.
<b>Status Move</b>	Ability to change the status of exams using Status Change Icons. The user can select a Status Change Icon in the worklist or from the right-click/Actions menus.
<b>Unlimited Advanced Search</b>	Ability to perform advanced searches without applying an Exam Date filter. From the Advanced Search, the user can click Unlimited Count and Unlimited Search.
<b>Use Advanced Search</b>	Access to Advanced Search.
<b>Use Folders</b>	Access to Folder tab and ability to associate Teaching Notes with Folders. The user can add/edit/delete folders.
<b>View exam report</b>	
<b>View group and enterprise worklists</b>	Access to worklist groups.